

Diagrammatic of Components/Spares for BA Internationals BAC200E Scaler Units

COMPONENTS

- Unit OPTIMA BAC200
- Water flow rate adjustment
- Power dial
- Handpiece holder
- Handpiece connection
- External power supply connection
- Foot control connection
- Bottle
- Handpiece cord
- 9A Handpiece cord OPTIMA BAC200E
9B Handpiece cord OPTIMA BAC200EL
- Handpiece
- 10A Handpiece OPTIMA BAC200E
10B Handpiece OPTIMA BAC200EL
- Power supply
- Power cord
- Foot pedal
- Scaling tips
- Torque wrench

BA International supplies units with various accessories. The "Packing list" shows exactly what is included with your unit.

HS Part No	Desc
BA-BAC20044068	Bottle 350ml (8)
BA-BAC20044064	BAC260E Handpiece Lead (Non-optic) (9A)
BA-BAC20044066	Handpiece only (non-optic) (10A) (or BA-BAC20044067 for the LED BAC200EL unit)
BA-BAC20044072	Power supply (11)
BA-BAC20044071	Foot pedal (13)
BA-BAC251E	Scaler 3-Tips Kit (14)* (Tips also avail. individually)
BA-BAC227	Handpiece Torque Wrench (for BAC200 Tips series) (15)
BA-WO-09	Handpiece O'rings
BA-BAC261E	Supragingival Scaler Tip (type "A") Each
BA-BAC262E	Supragingival Scaler Tip (type "P") Each
BA-BAC263E	Perio/Subgingival Scaler Tip (type "PS") Each

