# Cotton buff

For high-luster polish on acrylic teeth and facings using the handpiece.

### Advantages

- Optimum take-up of polishing paste, no sticking.
- Long service life.





#### **Ordering information**

Cotton buff mounted, Ø 22 mm (0.87"), Ø shaft 2.35 mm (0.09"), max. rpm 15 000 1/min, 12 pieces

No. 205 1000

# Dia-Finish L

High-quality felt discs impregnated with fine diamond particles for polishing porcelain and facing acrylics.

#### Advantages

- Optimum combination of polishing felt and diamond paste.
- For high lustre polishing of porcelain.
- Polish without firing after occlusal grinding (Time-Saving!).





#### **Ordering information**

Dia-Finish L wheel, medium-hard, Ø 12 mm (0.47"), max. rpm 15 000 1/min, 12 pieces	No. 400001
Dia-Finish L wheel, hard, Ø 12 mm (0.47"), max. rpm 15 000 1/min, 12 pieces	No. 400101
Dia-Finish L, knife-egde, hard, Ø 12 mm (0.47"), max. rpm 15000 1/min, 12 pieces	No. 400201
In general: Ø shaft 2.35 mm (0.09")	

## Silver wire brushes

For pre-polishing cast metal frameworks and palatal surfaces.

#### Advantages

- Does not remove surface stippled effect remains intact.
- Silky matt shine of frames.
- Otimum handling due to 51 mm diameter.





### **Ordering information**

Silver wire brushes, Ø ~51 mm (2.01"), max. rpm 3000 1/min, 2 pieces

No. 1960000