

Dustex master plus

Dust Extractor Box

The non-tip and non-slip Dustex master plus extractor box protects against splinters, dust and fumes. An ergonomically well-conceived form and height-adjustable arm rests enable free, comfortable working.

Technical data

Permissible mains voltage / frequency	220-240 V, 50/60 Hz 100-120 V, 50/60 Hz
Lamp power	16 W
Ø suction fittings, interior	35 mm (1.38")
Ø suction fittings, exterior	40 mm (1.57")
Weight	5 kg (11.0 lbs)
Dimensions (W x H x D)	380 x 285 x 400 mm (15.0 x 11.2 x 15.8")
Working chamber (volume)	17 l (4.5 gal)

Ordering information

Dustex master plus, 220–240 V, Dust box, fluorescent light with lamp cover, extraction port kit, armrests	No. 26260105
Dustex master plus, 100–120 V, Dust box, fluorescent light with lamp cover, extraction port kit, armrests	No. 26261105
1 High quality magnifier with flexible swan-neck	No. 26260300
2 Gold filter	No. 26260600

The suitable complement

SILENT TS workplace extraction unit → Page 35

Further information

www.renfert.com/P934







Advantages

- Extraction directly at the object ensures a clear view.
- Maximum freedom of movement thanks to a large inner volume (17 l).
- PerfectView: Innovative LED Technology for best contour and detail recognition.

Details

- Comfortable working due to large arm inlets.
- Bright, built-in illumination with 4800 Lux.
- High safety thanks to laminated protective screen with silicone coating.
- Versatile use due to compact and sturdy, light-weight powder-coated housing.
- Practical wooden support, made of hard beech wood.