Extractor clamp

With this practical accessory your dental lab work station can be subsequently fitted with an extractor unit.

Advantages

- Robust clamp made of beech wood.
- Secure fixture to work tables of 20–40 mm thickness.
- Variable and individual connection for extraction hose on the right or on the left (fitting 40 mm / 1.58").



Ordering information

Extractor clamp (without glass shield)	No. 2925 0000
Glass shield with holder	No. 2925 1000

Grinding box

Acrylic grinding box with built-in magnifying glass.

Advantages

- Prevents ground dust from spreading.
- Very easy and particularly flexible to use.
- Optimal control of the piece of work with integrated 2x magnifying glass.



Ordering information

Grinding box No. 2625 0000

Extractor Switch

For connecting 2 dust producing devices to Renfert extraction units

Advantages

- Via selection switch individually adjustable to one or two extraction openings.
- Tool free and fast connection.



Ordering information

Extractor Switch incl. 4 mufflers (without hose) Ø suction fittings, interior/exterior 35 mm / 40 mm (1.38" / 1.57") Ø muffler interior/exterior 38 mm / 46 mm (1.50" / 1.81")

No. 29260000

binding intendifferential 30 min / 40 min (1.30 / 1.01 /	
End bushings set, 2 pieces	No. 29210003
Suction hose, 3 m (incl. 2 mufflers)	No. 900034240
Suction hose, 6 m	No. 901150823
Suction hose, 9 m	No. 902150823

H/Hepa Fine Filter

Even greater safety provided by a filtration efficiency of 99.995 %.

Advantages

- Toolless, quick change of the filter.
- Available for SILENT TS, SILENT TS2 and Vortex compact 3L.
- Each filter is delivered with an individual certificate (tested according to DIN EN/IEC 60335-2-69 A.3/AA.22.201.2 and EN 1822-5 Hepa Filters / MIL-STD 282 Type: C).



Ordering information

H+Hepa filter SILENT TS/TS2, set	No. 900034410
H+Hepa filter Vortex compact 3L, 1 piece	No. 900034420