# MT3

## Wet trimmer

The MT3 wet trimmer is impressive due to its exceptionally powerful motor. With the MT3 Renfert provides high-performance trimming at a low price.







"The stable front can be completely removed and cleaned in a few steps using the quick-release lock. The spray tube can also be removed as easily as the trimmer table and could be cleaned in plaster remover for example (we recommend Renfert's GO 2011 speed → page 189)."

Oliver Bothe, Product Management, **Master Dental Technician** 

### **Advantages**

- Quick change of the trimmer table angle using the practical tipping mechanism (90°/98°).
- Safety: when the front door is opened a safety switch stops the motor and water.
- Uniform coverage by the optimally positioned spray tube prevents clogging of the trimmer disc.

## **Details**

- A 10° inclination of the unit provides an optimum view of the model.
- Easy guidance of the plaster model with the aid of slotted angle lines on the trimmer table.

## Tip

Do not discard the old Klettfix trimmer disc. It can be used as an ideal aid to remove excess investment from a CrCo model manually.

## **Technical data**

Permissible mains voltage / frequency	220-240 V, 50/60 Hz   100-120 V, 50/60 Hz
Power consumption	1.300 W / 1.74 hp (220–240 V) 1.325 W /1.78 hp (120 V)
Motor performance P1	1.300 W / 1.74 hp (220–240 V) 1.325 W /1.78 hp (120 V)
Motor performance P2	900 W / 1.21 hp (220–240 V) 1.000 W /1.34 hp (120 V)
Rotational speed	2.880 1/min (50 Hz) / 3.400 1/min (60 Hz)
Water consumption max.	7 l/min (0.25 cfm)
Water pressure min. / max.	min. 1 bar / max. 5 bar (14.5–72.5 psi)
Ø Connection piece for drain hose	36 mm (1.42")
Dimensions (W x H x D)	305 x 330 x 410 mm (12.0 x 13.0 x 16.1")
Trimmer disc dimensions	Ø 234 mm (Ø 9.2")
Weight (without trimmer disc)	13,2 kg (29 lbs)

## **Ordering information**

MT3, 220–240 V, incl. Klettfix trimmer disc	No. 1808 0000
MT3, 220–240 V, incl. Marathon trimmer disc	No. 18080500
MT3, 100–120 V, incl. Klettfix trimmer disc	No. 1808 1000
MT3, 100–120 V, incl. Marathon trimmer disc	No. 1808 1500

## **Suitable accessories**

Trimmer discs → Page 72

## **Further information**

www.renfert.com/P806



