

**Renfert**

making work easy



# SYMPRO

Denture cleaning unit

# SYMPRO

## Denture cleaning unit

Compact, high-performance cleaning unit especially suitable for dentures, orthodontic appliances and other dental restorations.



### Hygienically Treated

Hygienic preparation of the cleaning bowl, cleaning needles and tweezers can be carried out either by means of a manual or mechanical/thermal disinfection process. The efficiency of the respective treatment process has been validated by an accredited test laboratory according to the standards EN ISO 17664, EN ISO 15883-2, ISO / TS 15883-5, AAMI TIR 30:2011.

### Advantages

- Up to 80% time savings compared to conventional methods thanks to effective and automated cleaning process.
- Easy cleaning of temporary restorations in combination with SYMPRO mini cup and temp:ex temporary cement remover.
- Validated, hygienic reprocessing of system components.

### Details

- Increased customer retention with a prophylaxis service for denture wearers.
- Optimal and fast cleaning performance due to an inclined bowl at 35°.

### Technical data

Permissible mains voltage / frequency	100–240 V, 50/60 Hz
Power consumption	90 VA
Rotational speed	1200 / 1600 / 2000 1/min
Dimensions (W x H x D)	150 x 240 x 280 mm (5.91 x 9.45 x 11.02")
Weight w/o bowl	approx. 3 kg (6.6 lbs.)

### Ordering information

**SYMPRO**, 100–240 V | 50/60 Hz

No. 6500 0000

### Scope of delivery SYMPRO

SYMPRO, tweezer, cleaning vessel with cover for use in cleaning device, cleaning bowl with cover for post-treatment, retaining magnet, cleaning needles, manual for hygienic preparation

### Accessories

1	Cleaning pins, 75 g (2.63 oz.)	No. 6500 0550
2	SYMPROfluid Universal, 2 x 2 l (2 x 0.53 gal)	No. 6500 0600
3	SYMPROfluid Nicoclean, 2 x 2 l (2 x 0.53 gal)	No. 6500 0610
4	SYMPRO mini-cup set	No. 6500 0410

1



2



3



4



### Further information

[www.renfert.com/sympro](http://www.renfert.com/sympro)

