

NOW **eASIER** THAN EVER

4 SHADES =

Simple, brilliant restorations, every time!



Your Smile. Our Vision.
www.sdi.com.au



aura eASY

ultra universal restorative material

aura eASY it's the easy choice for your practice

2.21%

Low volumetric shrinkage

aura eASY delivers reduced stress to provide long-lasting restorations.

115 MPa

Excellent flexural strength

aura eASY resists major deformations without fracturing.

375 MPa

High compressive strength

aura eASY resists the high forces of compression in the mouth.

- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- No slumping and holds its shape
- Quick and easy placement in seconds



It will truly allow you to deliver more while stocking less.

A1/B1 ≈ ae1	A3/B3 ≈ ae3
A2/B2 ≈ ae2	A3.5/B4 ≈ ae4

Approximate equivalence to VITA classical shade guide

PART NO.	QTY	SYRINGES
8565106	4x1	aura eASY Syringe Kit 4g syringe (ae1, ae2, ae3, ae4)
	1	aura eASY shade guide
aura eASY Syringe Refills		
8560014	1	ae1 – 4g
8560015	1	ae2 – 4g
8560016	1	ae3 – 4g
8560017	1	ae4 – 4g



PART NO.	QTY	COMPLETS
8565107	20x4	aura eASY Complet Kit 0.25g complet (ae1, ae2, ae3, ae4)
	1	aura eASY shade guide
aura eASY Complet Refills		
8565014	20	ae1 – 0.25g
8565015	20	ae2 – 0.25g
8565016	20	ae3 – 0.25g
8565017	20	ae4 – 0.25g



Your Smile. Our Vision.
www.sdi.com.au

Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Ireland 00800 0225 5734

Italy 00800 0225 5734
New Zealand 0800 734 034
Spain 00800 0225 5734
United Kingdom 00800 0225 5734
USA & Canada 1 800 228 5166

aura eASY
ultra universal restorative material

