


SAND

#4220012 (1kg bottle)
#4240012 (4kg bottle)



LIGHT GREY

#4220010 (1kg bottle)
#4240010 (4kg bottle)



IVORY

#4220009 (1kg bottle)
#4240009 (4kg bottle)



KeyModel Ultra®


KeyModel Ultra is a material designed for 3D printing of dental and orthodontic models.



ULTRA-FAST
TECHNOLOGY



FLAWLESS
SURFACE DETAIL



THERMOFORMING
RELEASE AGENT



TPO FREE

FEATURES

- HIGH-SPEED PRINT FORMULATION
- INTEGRATED RELEASE AGENT FOR EASY THERMOFORM SEPARATION
- SMOOTH SURFACE FINISH
- STABILITY IN PRODUCTION WORKFLOW
- INTEGRATED WITH A WIDE RANGE OF 3D PRINTERS

PROPERTY	METHOD	*KEYMODEL ULTRA
Maximum Tensile Strength	ASTM D638	>50 MPa
Tensile Modulus	ASTM D638	1700 MPa
Elongation at Break	ASTM D638	5%
Maximum Flexural Strength	ASTM D790	>70 MPa
Flexural Modulus	ASTM D790	>1800 MPa
Viscosity @ 25°C	ASTM D7867	500-600 cP

*Typical Values



keystone[®]
industries

keyprint.com | sales@keystoneind.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH, Stockholzstr. 11, 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01

REV 1 09/2025