

**AH,
THAT'S
BETTER**




CATTANI – FOR PEACE AND QUIET, AND PEACE OF MIND

Cattani's world-leading engineers created SMART Cube, the suction system that's HALF AS LOUD*

Have Confidence in Cattani

For more than 50 years, Cattani's global business has been dental air technology – no distractions. We stand for quality, performance and reliability with the aim of keeping the dental practice running day after day, year after year.

*A conservative approximation based on a reduction in decibel level between the standard suction unit and the Cube of about 10 decibels, roughly in the 70 down to 60 range. In that range, a reduction of 10 decibels halves perceived loudness.

A woman with dark hair in a braid, wearing a red t-shirt and black pants, is sitting in a meditative pose on a grassy field. She is looking to the right with her eyes closed. The background is a blurred landscape of a field with some trees in the distance under a bright sky.

**It's time for Cattani
SMART Cube, the quieter,
smarter, more economical
suction system**

Quieter – acoustic housing cuts noise by half. Easy opening for service and maintenance. Unit has its own cooling system.

Smarter – Cattani-pioneered inverter control ensures power output matches chair use. Parallel units also provide built-in redundancy to protect and secure surgery workflow.

More economical – consistent air flow = optimal energy use = real savings. Unit is also environmentally friendly, being more compact and using fewer raw materials.

3-year warranty – we stand by our product, and nothing says reliability like a long warranty.



Cattani SMART Cube

- Acoustic housing cuts noise by 50%
- Inverter control slashes energy consumption
- Consistent vacuum performance
- Easy to add parallel units
- Amalgam retention optional
- Compact, all-in-one design
- Generous 3-year parts and labour warranty
- Unrivalled technical service back-up



WE LOVE WHAT WE DO.

**It's going to be a lot quieter
from now on.**

SMART Cube suction systems from Cattani

Quieter • Smarter • More efficient



Find tranquillity at Cattani

Australia

+61 3 9484 1120

cattani@cattani.com.au

www.cattani.com.au

280 Dundas Street
Thornbury VIC 3071

New Zealand

0800 68 22 88

cattani@cattani.co.nz

www.cattani.co.nz