()

TISSUE MANAGEMENT

ViscoStat™

20% FERRIC SULFATE











- Provides profound hemostasis
- Stops moderate bleeding
- Saves chair time
- Does not impede hard or soft tissue healing
- Eliminates sulcular fluid contamination for optimal bonding
- Decreases costly impression remakes

ViscoStat hemostatic is a 20% ferric sulfate equivalent hemostatic agent with inert binding agents in a viscous, aqueous solution. It contains fumed silica to limit the acidic activity, making it kind to hard and soft tissue.

ViscoStat hemostatic solution is suited for a variety of dental and oral surgery procedures to arrest surface capillary bleeding. Such procedures include fixed prosthodontics, restorative-operative, periodontal treatment, etc. ViscoStat hemostatic is also recommended for retrofillings, canine impactions, gingivectomies, and as a "fixative" for pulpotomies.

Tip: Prevent leakage caused by sulcular fluid contamination during direct bonding procedures. Soak an Ultrapak™ knitted cord in a hemostatic and isolate the tissues. Follow with a firm air/water spray.

"ViscoStat hemostatic has allowed me to take clean, dry, accurate impressions even in some of the most hemorrhagic situations. Thank you! I appreciate it, my lab man appreciates it, and my patients appreciate it."

—DR. MONTE PERSON – FRESNO, CA

"ViscoStat hemostatic can stop bleeding in the cervical area of a crown prep like nothing else I have ever seen! It is a miracle solution! I wish I knew how many thousands of hours it has saved me over the years. Thank you, Ultradent!"

—DR. SCOTT J. HADLEY – HAXTUN, CO

"Without a doubt I would be lost without ViscoStat hemostatic! It quickly stops sulcular bleeding, allowing me to place the Ultrapak #000 cord and get an excellent impression on the first try. It is also great to use when bonding composites near gingival tissues. ViscoStat hemostatic is a vital part of my crown prep and composite armamentarium."

—DR. JULIE ANN ROUTHIER – SAVANNAH, GA



3093 - ViscoStat Dento-Infusor Syringe Kit

4 x 1.2 ml syringes 20 x Metal Dento-Infusor tips



647 - ViscoStat Dento-Infusor IndiSpense™ Syringe Kit

1 x 30 ml IndiSpense syringe 20 x Metal Dento-Infusor tips 20 x 1.2 ml empty syringes



3070 - ViscoStat Syringe 20pk
1.2 ml syringes



645 - ViscoStat IndiSpense Syringe 1pk
30 ml syringe



1278 - ViscoStat Empty Syringe 20pk
1.2 ml empty syringes

Note: Do not use epinephrine preparations with ferric sulfate products (ViscoStat, Astringedent), as blue/black precipitate will occur.

1. realityesthetics.com.