

> FACTS & FEATURES

What?

• lifelike, elaborate premium anterior teeth made of MRP composite

Benefits?

- simple fabrication of highly esthetic dentures due to lifelike tooth moulds with "golden" proportions
- brilliant play of shade and light thanks to excellent light dynamics and harmonious shade progression
- good shade reproduction with natural remaining dentition due to ideal shade accuracy to the shade guide
- excellent durability
 due to high load capacity and abrasion stability
- excellent gingiva design thanks to broadly dimensioned tooth neck

What for? • ideal solution for highly esthetic full, partial and implant-retained dentures.

VITAPAN EXCELL®





> BASIC INFORMATION

Characteristics

Mould

- lifelike three-dimensional tooth geometry for brilliant results
- distinctive incisal characteristics for special vitality
- wide palatal bars for ideal interdental spacing/papilla creation

Special features

- pronounced enamel portion for good light transmission
- elaborate surface structure for natural diffusion and reflection of light
- integrated mamelon/opalescence for natural shade effect
- 3D layered structure for natural shade gradient

Recommended indications

Ideal prosthetic solution for:

- full and partial dentures
- implant-retained dentures

Occlusal/setup concept

Suitable for all established concepts of lingual and balanced occlusion:

e.g., according to Gysi, Prof. Dr. A. Gerber, APF/TiF and ABC


Available shades

- 11 x VITA SYSTEM 3D-MASTER® shades 0M3, 1M1, 1M2, 2L1.5, 2M1, 2M2, 2R1.5, 3L2.5, 3M2, 3R2.5, 4M2
- 15 x VITA classical A1–D4® shades A1–D4 (without B1)

Available moulds/sets

- 17 x upper anterior tooth moulds (divided into four shape groups)
- 5 x lower anterior tooth moulds

> **OVERVIEW OF SHAPES**





> RECOMMENDED COMBINATIONS

Combinations according to indications:

For full/partial dentures and implant-retained dentures:

VITAPAN EXCELL Anterior in combination with VITAPAN LINGOFORM Posterior

| VITAPAN EXCELL Anterior | | VITAPAN LINGOFORM Posterior | |
|--------------------------------------|---------------------|-----------------------------------|--------------------|
| UJ / No. | LJ / No. | UJ / No. | LJ / No. |
| 041 41.3 | L33 33.2 | 21L 30.1 | 21L 30.7 |
| 043 43.1 | L33 33.2 | 21L 30.1 | 21L 30.7 |
| 045 45.1 | L35 35.3 | 22L 31.6 | 22L 32.1 |
| 047 47.3 | L37 37.3 | 22L 31.6 | 22L 32.1 |
| 049 49.1 | L39 39.0 | 23L 34.2 | 23L 34.8 |
| 051 51.4 | L 41 41.2 | 24L 35.8 | 24L 36.4 |

upper jaw (UJ) / lower jaw (LJ)

| VITAPAN EXCELL Anterior | | VITAPAN LINGOFORM Posterior | |
|--------------------------------------|--------------------|-----------------------------------|--------------------|
| UJ / No. | LJ / No. | UJ / No. | LJ / No. |
| T42 42.0 | L33 33.2 | 21L 30.1 | 21L 30.7 |
| T44 44.1 | L33 33.2 | 21L 30.1 | 21L 30.7 |
| T46 46.2 | L35 35.3 | 22L 31.6 | 22L 32.1 |
| T47 46.9 | L37 37.3 | 22L 31.6 | 22L 32.1 |
| T48 48.1 | L37 37.3 | 22L 31.6 | 22L 32.1 |
| T50 50.0 | L39 39.0 | 24L 35.8 | 24L 36.4 |
| T54 54.4 | L41 41.2 | 24L 35.8 | 24L 36.4 |
| R45 45.4 | L35 35.3 | 22L 31.6 | 22L 32.1 |
| R47 47.0 | L37 37.3 | 22L 31.6 | 22L 32.1 |
| R49 49.1 | L39 39.0 | 23L 34.2 | 23L 34.8 |
| \$46 46.3 | L35 35.3 | 22L 31.6 | 22L 32.1 |

> ASSORTMENTS AND STORAGE





VITAPAN EXCELL mould charts

Metal case with all moulds for selecting teeth:

VITAPAN EXCELL Anterior


VITAPAN EXCELL tooth storage

- X-small*: 67 Sets / 402 Anterior
- Small*: 134 Sets / 804 Anterior
- Medium*: 268 Sets / 1,608 Anterior
- Large*: 400 Sets / 2,400 Anterior


Available in VITA classical A1-D4® and VITA SYSTEM 3D-MASTER® shades **

- *) Note: You will receive the most frequently used moulds and shades with the packaged tooth storage. We are also happy to put together storage systems according to your individual specifications.
- **) See available shades

> SYSTEM COMPONENTS







youthful, lifelike moulds