

iSense Extra Nitrile® Gloves

Intended Use

To protect the user from harmful bio hazards and substance contamination within medical and pharmaceutical work environments.

Trade Name		UOM	Materials	Storage Condition		
iSense Extra Nitrile Gloves		rile	100 gloves	100% Nitrile. No powder lubricant added.	Store in a dry, ventilated area t maintain the properties. Avoid direct sunlight, fluorescent lighting, head and moisture.	
Type Nitrile Examination Glove, 290mm extended cuff, Powder Free, Online Single Chlorinated, Non-sterile, Ambidextrous. 5.2g		Colour White	Surface Finish Textured fingertips.	Shelf Life The gloves shall have shelf life of 5 years from the date of manufacture with the above storage conditions. Do not store above 38°C as this will lead to accelerated aging.		
	Dimensions Millimetres	Size X Small Small Medium Large X Large	Palm width 73 – 82 83 – 92 93 – 102 103 – 112 113 – 122	Length (mm) 290 (± 10mm) 290 (± 10mm) 290 (± 10mm) 290 (± 10mm) 290 (± 10mm)	Physical PropertiesTensile Strength (Mpa)Elongation (%)Force at Break: (ASTM D 6319Country of Origin Malay	,
	Thickness	Palm	Finger	Cuff		

Min 0.12 – 0.14

Min 0.05





Marking of gloves protecting against virus, bacteria and fungi. Protection against bacteria & fungi- Pass Protection against viruses - Pass

iSense Extra meets ASTM D6319 Current Revision and Customer Requirement

Pack Part No.

Small	230011
Medium	230012
Large	230013
Extra Large	230014

Australian Government Department of Health

TGA approved Class 1 medical device

Before Aging	After Aging
Min 14	Min 14
Min 550	Min 400
	Min 14

Intry of Origin Malaysia



Designing and manufacturing the world's best disposable and reusable protective gloves

Packaging Configuration & Labelling

Min 0.07

Millimetres

Packaging style	100 gloves x 10 boxes x 1 carton	
Packaging	100 gloves per inner box (by weight) 10 boxes per outer carton.	
Labeling	The inner boxes and outer cartons have the following information: Product description, brand, size, and net quantity of gloves. Country of origin, name and address of distributor, and single use only symbol.	
Barcodes	GS1 format	

100



www.theglovecompany.com.au