## **Water Test Protocol / Record**

Protocol for documenting and recording dental unit waterline compliance. To meet NZ drinking water standards water quality must be below 200 CFU/m for when treating immune compromised patients.

Sampling location (room)							
Sample No.		Da	te	//	Time .	am/pm	1
Treatment Unit:							
Name							
Input Water Quality:							
Please take sample from in	nput supply e.g. tap water	r bottle	is fille	ed from			
Sampling location							
Name of decontamination (e.g. Dentapure cartride	-	l water	<b>:</b>				
Brand Name							
Place of sampling water	r: Please Tick (🗸)						
Tumbler Filler							
Multifunction syring	e (Dentist side)		Multif	unction syrin	ige (Ass	istant side)	
Turbine (Highspeed	l line)						
Motor 1	Motor 2						
Scaler Device							
Mixed sample (Den	tist side)		Mixed	d sample (As	sistant s	side)	
Sampling by: Please Tic	k ( <b>√</b> )						
Staff member		U Te	echni	cian/Sales R	epresen	tative	
Name	Się	gnature	)				
Result: Please Tick (  )							
This result (Colony-forming Unit	s - CFU) is based on the exact	procedu	ire acc	cording to the p	roduct inf	formation	
No CFU Less th	nan 200 CFU/mL 🔲 Les	ss than	500	CFU/mL	More	than 500 CFU/r	nl
Test result determined Da	ate/ Time		. am/	pm No.of da	ays after	sampling	
<b>Evaluation:</b> Please Tick	( <b>/</b> ) Negative	Ţ	<b>]</b> 5	Shock treatm	ent requ	uired	