

EXCLUSIVE TO HENRY SCHEIN

PLANMECA

IMAGING SALE NOW ON



ProX Intraoral X-ray

With over **800 installations** to NZ practices over the last decade, we're not embellishing when we say the Planmeca ProX is one of the most popular intraoral x-ray units on the market.

\$5,700 ~~\$6,750~~



ProSensor HD Intraoral Digital Sensor

Combining the ProX intraoral X-ray unit with the Planmeca ProSensor HD sensor elevates user-friendliness to a new level – allowing images to be displayed on-screen just seconds after exposure. With a true resolution of over 20 lp/mm, the intraoral sensor delivers exceptional image quality.

FROM
\$8,230



ProMax CBCT Promax 3D Classic

Package comes with all you need to start scanning patients right away.

Includes: SmartPan System, Ultra Low Dose Technology, ProMax Extended FOV, Planmeca CALM and a Romexis Software Package.

\$81,500 ~~\$85,550~~

Need a Ceph Arm?

Promax 3D Classic with Ceph Now Only \$99,900



ProOne OPG 2D Imaging

Easy Patient Positioning - monitor patients from the front and the side, with positioning further assisted by a triple laser beam system that indicates the correct anatomical positioning points.

Guided User Interface - full-colour graphical user interface that helps guide you through procedures with clear texts and symbols.

Autofocus - this feature automatically positions the focal layer using a low dose scout image of your patient's central incisors as well as landmarks in the patient's anatomy to calculate placement **resulting in perfect panoramics every time.**

\$27,950 ~~\$33,980~~



Terms and Conditions

PRICES ARE EFFECTIVE FROM 1/9/2024 TO 01/11/2024 AND EXCLUDE GST, FREIGHT, INSTALLATION AND SOFTWARE UNLESS OTHERWISE SPECIFIED. PRICES OR OFFERS IN THIS PUBLICATION ARE NOT APPLICABLE TO GOVERNMENT OR CORPORATE ACCOUNTS. EQUIPMENT CONFIGURATIONS MAY DIFFER FROM THE IMAGES SHOWN IN THIS FLYER.