

 **HENRY SCHEIN**<sup>®</sup>  
DENTAL | INSTRUMENTS

**Rely** on Us<sup>™</sup>

# RELY ON US

every smile counts



2018  
**INSTRUMENTS**  
specialty catalog

# Welcome to the Henry Schein Dental Instrument Catalogue

## Page Guide

3 - 8

Diagnostics

---

9 - 15

Hygiene

---

16 - 22

Operative Instruments

---

23 - 39

Surgical

---

40 - 41

Surgical Accessories

---

42 - 46

Specialty Instruments

---

47

Laboratory Instruments

---

48 - 55

Orthodontic Instruments

---

56 - 59

Instrument Cassettes & Accessories

---

60

Sharpening Accessories

---



# Diagnostics

## Cotton & Dressing Pliers



How the best perform

HF-DP2 College - Serrated



HF-DP18L Self-locking



HF-DP1 Perry



HF-DPU17 Utility Pick-Up



AE-DA240 College - Smooth 150mm



AE-DA241 College 150mm



AE-DA251 College - Smooth 150mm



AE-DA271 Self-locking - Serrated 150mm



AE-DA204 Perry - Serrated 130mm



AE-DA225 Meriam - Serrated 160mm



## Anatomic Pliers



### Dressing Pliers

- HS-9002011 13cm
- HS-9003147 14.5cm
- HS-9002012 18cm



### Adson Tissue Pliers

- HS-1005659 12cm Serrated



## Tweezers



Premium quality. Satin finish with precise nonslip tips and serrated handles. Made of surgical-grade stainless steel.

### Air Tweezer

- HS-9008013 #317 Serrated



### Meriam Air Tweezer

- HS-9008014 Serrated



### Meriam Air Tweezer

- HS-9008016 Self Locking Serrated



### College Tweezer

- HS-1003313 #317 Serrated 15cm



### Meriam Cotton Dressing Tweezer

- HS-1005706 Serrated 16cm



### Tweezer

- HS-1011514 Self Locking Serrated



### College Tweezer

- HS-9002059 Smooth 15cm



## Explorers - Single Ended



How the best perform

HF-EXS6



HF-EXS6L



HF-EXS23



HF-EXS54



HF-EXS17



AE-DA406



AE-DA416



Beach Explorer 23B - light weight (10g) and slim (3.5mm diameter) enhances tactile evaluation of enamel terrain and reduces finger fatigue.



M-101121



## Explorers - Double Ended



How the best perform

HF-EXD5



HF-EXD6



HF-EXD11



HF-EXD2



HF-EXS23-6



HF-EXTU17-23



## Probes

### Qulix Colour Coded Probes

Durable Qulix colour coded surface accurately marked and smoothly finished. Designed for repeated ultrasonic cleaning and autoclaving. Will not wear down, chip, fade or catch in periodontal pockets.



How the best perform

HF-PCP115B WHO



HF-PCP11 3-6-8-11



HF-PCP12 3-6-9-12



HF-PQOW Michigan O with Williams Markings 1-2-3-5-7-8-9-10



HF-PCPUNC15 University of North Carolina 1-2-3-4-5-6-8-9-10-11-12-13-14-15



HF-PCP2 2-4-6-8-10-12  
HF-PCP8 3-6-8-11  
HF-PCP12 3-6-9-12



### MAINTENANCE TIP


Hu-Friedy QULIX Colour-Coded Probes may be ultrasonically cleaned, autoclaved, chemiclaved, or dry-heat sterilised and will withstand temperatures up to 350°F/177°C.

# Diagnostics

## Colorvue Perioscreen Probes



Easy to read, ergonomically designed, eco-friendly. Hu-Friedy's Colorvue Perioscreen probe provides a fast and easy-to-use method of screening periodontal pocket depth prior to definitive probing. The universal colour-coding of green and red indicates absence or presence of periodontal disease. The Colorvue Perioscreen features a metal handle that is not only more ergonomic, but can be reused or recycled.




- HF-PPS10KIT12 Kit #10 2 Handles & 12 Tips
- HF-PPS10KIT6 Kit #10 1 Handle & 7 Tips
- HF-PPS10PT Refill Tips #10 Pack of 12
- HF-PPS12KIT12 Kit #12 2 Handles & 12 Tips
- HF-PPS12KIT6 Kit #12 1 Handle & 7 Tips
- HF-PPS12PT Refill Tips #12 Pack of 12
- HF-PH6 Colorvue Satin Steel Handle

## Colorvue Probes



Feature a vivid yellow tip with black markings for enhanced contrast to intraoral structures. Safe for use with implants and feature twist on design for easy replacement of worn tips.



- HF-PCV11PT 3-6-8-11 Pkt 12 Tips
- HF-PCV12PT 3-6-9-12 Pkt 12 Tips



- HF-PH6 Colorvue Satin Steel Handle
- Kits of Handle + Tips Also Available.

## Air Collection Single Ended Probes



#23 SE

- HS-9008005 Air Handle
- HS-1003620 Solid Round Handle



#8 SE

- HS-9008007 Air Handle
- HS-9003163 Solid Round Handle



#9 SE

- HS-9008008 Air Handle
- HS-9003164 Solid Round Handle



#17 SE

- HS-9001662 Solid Octagonal Handle



## Air Collection Double Ended Probes



#8/9 DE

- HS-9008589 Air Handle



#5 DE

- HS-9004403 Air Handle
- HS-1008008 Solid Round Handle



#16 DG

- HS-9004404 DE Air Handle
- HS-1009279 DE Solid Round Handle



## Air Collection Single Ended Probes



CPITN SE

- HS-9008017 Air Handle
- HS-9001659 Solid Round Handle



PCP/14 SE

- HS-9008020 Air Handle



Williams SE Air

- HS-9008022 Solid Round Handle



PCP/15 SE

- HS-9008021 Air Handle
- HS-9002071 Solid Round Handle



# Diagnostic

## Periodontal Pocket Probes



HF-PW Williams  
Single Ended  
1-2-3-5-7-8-9-10

HF-PO Michigan "O"  
3-6-8

HF-PGF  
1-2-3-5-7-8-9-10

HF-P2N Double Ended Nabers Furcation



## Novatech Right Angled Probes



Periodonal pocket probes with a unique right angled parallel tip design facilitating placement of the tip into the interproximal and lingual aspects of posterior teeth.

HF-PCPNT2 3-6-9-12

HF-PCPNT115B WHO

HF-PCPNT11 3-6-8-11



## Beach W.H.O. Periodontal Probe



World Health Organisation endorsed periodontal probe allows for tactile analysis of calculus, while facilitating visual confirmation of pocket depth. Light weight (4.5g); rounded tip will not cut healthy tissue; angular relationship of probe to handle allows operator to work in a balanced posture. Ideal for rapid survey.

M-101970




## Periodontal Probes



AE-DB767 Who 160mm Single Ended

AE-DB770 Fox 155mm Single Ended



## Explorer-Probe Combinations



HF-XP23-0

HF-XP23-W  
HF-XP23-11.5B



## Mirrors

### HD Mirrors



Feature a proprietary coating that provides superior brilliance and colour for increased visibility in all dental procedures. Also scratch-resistant, our HD mirrors are the clear choice for optimal patient treatment.

- 113% reflection factor for exceptional image clarity
- 38.5% brighter than rhodium coated mirrors
- 50% brighter than other front surface mirrors
- Scratch-resistant for product longevity
- Large ergonomic handle options

### Single-Sided Cone Sockets

- #3  
HF-MIR3HD Single
- #4  
HF-MIR4HD Single  
HF-MIR4HD3 Box 3  
HF-MIR4HD12 Box 12
- #5  
HF-MIR5HD Single  
HF-MIR5HD3 Box 3  
HF-MIR5HD12 Box 12



## Double Sided Cone Socket

- #4  
HF-MIR4DSHD Single  
HF-MIR4DSHD6 Box 6
- #5  
HF-MIR5DSHD Single  
HF-MIR5DSHD6 Box 6



## Single-Sided Simple Stem Socket

- #3  
HF-MIR3HDE Single
- #4  
HF-MIR4HDE Single  
HF-MIR4HD3E Box 3  
HF-MIR4HD12E Box 12
- #5  
HF-MIR5HDE Single  
HF-MIR5HD3E Box 3  
HF-MIR5HD12E Box 12



## Double-Sided Simple Stem Socket

- |                     |                     |
|---------------------|---------------------|
| #4                  | #5                  |
| HF-MIR4DSHDE Single | HF-MIR5DSHDE Single |
| HF-MIR4DSHD6E Box 6 | HF-MIR5DSHD6E Box 6 |

## Front Surface Mirrors (Stainless Steel) Cone-Socket Single-Sided

- Pack of 3
- HF-MIR43 Size 4
- HF-MIR53 Size 5

## Front Surface Mirrors (Stainless Steel) Cone-Socket Single-Sided

- Pack of 12
- HF-MIR412 Size 4
- HF-MIR512 Size 5



## Front Surface Mirrors (Stainless Steel) Cone-Socket Double-Sided

- Pack of 3
- HF-MIR4DS3 Size 4
- HF-MIR5DS3 Size 5

## Front Surface Mirrors (Stainless Steel) Cone-Socket Double-Sided

- Pack of 6
- HF-MIR4DS6 Size 4
- HF-MIR5DS6 Size 5

## Maxima Mouth Mirrors



Premium quality mouth mirrors made from high grade stainless steel. Simple Stem, Stainless Steel, Pkt 12.

### Plane Surface Mirror

Plane surface mirrors with a longer lifetime and less scratches on the surface.

- HS-9000753 Size #3
- HS-9000754 Size #4
- HS-9000755 Size #5

### Front Surface Mirror

Front surface mirrors for a true and clear mirror reflection.

- HS-9000747 Size #3
- HS-9000748 Size #4
- HS-9000749 Size #5

### Magnifying Mirror

Plane surface mirrors that are magnifying and have a rhodium coating.

- HS-9000798 Size #4

## Cone Socket Mouth Mirrors



Stainless Steel, Pkt 12

- |                             |                             |
|-----------------------------|-----------------------------|
| <b>Plane Surface Mirror</b> | <b>Front Surface Mirror</b> |
| HS100-5243 Size #3          | HS104-9213 Size #3          |
| HS100-6714 Size #4          | HS104-2350 Size #4          |

## Beach Mirror



Small (5g), ergonomically designed autoclavable mirror. Proper angulation allows for balanced operation. Small face maneuvers better in tight spots.

- M-107130 Set - 5 Handles & 5 Mirrors #2
- M-107150 Beach Mirror #2 Plain Box 10
- M-107000aFS2 Beach Mirror Front Surface #2 Pack

# Diagnostic

## Surface Mirror

The Standard Mirror

- Fine image without refractions under mid luminous conditions
- Ideal for standard treatment like diagnostics, scaling
- Pack of 12

AE-DA026 Surface 22mm  
AE-DA027 Surface 24mm



## Rhodium Mirror

The Premium Mirror

- Outstanding image without any refractions, even under higher luminous conditions
- Excellently for surgery, treatment of deep cavities, endodontic, implantology etc.
- Pack of 12

AE-DA036 Rhodium 22mm  
AE-DA037 Rhodium 24mm



## Mirror Handles



How the best perform



HF-MHE6 Satin Steel #6 Simple Stem  
HF-MHE7 Satin Steel #7 Simple Stem



HS-1006633 Round Simple Stem



M-107145 Beach

## Diagnostic Instruments

### Air Diagnostic Instrument Sets



The Air Diagnostic Instrument set #23 consists of a Probe SE #23, mouth mirror with handle and a College Tweezer #317, serrated. All light as Air in environmental friendly packaging.

HS-9008184 Set #23 Air



The Air Diagnostic Instrument set #9 consists of an Explorer #9, mouth mirror with handle and a College Tweezer. All light as Air in environmental friendly packaging.

HS-9008401 Set #9 Air



AE-DA80 Octagonal/Knurled



AE-DA90 Ergonomic Shape



HA-413 Hahnenkratt - Octagonal Smooth



HA-426 Hahnenkratt - Round Serrated Handle



## Gracey Curettes

### Double-End Handle Options



How the best perform

Most popular instruments available in most handle options.



#9 EverEdge 2.0™ Technology



#8 Resin Eight Colours™



#7 Satin Steel® Colours™



#6 Satin Steel®



#4 Round



#2 Octagonal

### Satin Steel Colours Signature Series Grips



How the best perform

- Pack of 16
- HF-SSG1 Grey
- HF-SSG3 Orange
- HF-SSG5 Yellow
- HF-SSG6 Red
- HF-SSG7 Green
- HF-SSG8 Blue
- HF-SSG9 Black
- HF-SSG10 Pink
- HF-SSG0 Assorted



### BioGent



How the best perform

Bio-type based periodontal therapy is based on maintaining the integrity and conservation of the patient's tissue. This therapy provides greater patient comfort and promotes optimal oral health. The proprietary design of Hu-Friedy BioGent Curettes provides a more tissue-friendly approach. These instruments are slim enough to provide enhanced pocket access. All BioGent Curettes are designed with EverEdge Technology for super-sharp edges and efficient scaling.



- HF-SG129LC 1/2 All Surfaces Anterior Grey
- HF-SG789LC 7/8 Buccal & Lingual Surfaces Posterior Green
- HF-SG11129LC 11/12 Mesial Surfaces Posterior Orange
- HF-SG13149LC 13/14 Distal Surfaces Posterior Blue

### Double Ended - Area Specific



1/2 Slight contra-angle for incisors and canines.

HF-SG1-2



3/4 Short contra-angle for incisors and canines.

HF-SG3-4



5/6 Medium contra-angle for anteriors.

HF-SG5-6



7/8 Medium contra-angle for lingual surfaces of premolars and molars.

HF-SG7-8



9/10 Long contra-angle for molars.

HF-SG9-10

# Hygiene



11/12 For mesioposterior surfaces. Designed to reach mesial surfaces of all posterior teeth with minimum cheek distention.

HF-SG11-12



13/14 For distoposterior surfaces.

HF-SG13-14



15/16 Same contra-angle bends as 13/14, but blade positioned to reach mesial posterior surfaces.

HF-SG15-16



17/18 Accentuated angles improve access to distal posterior surfaces.

HF-SG17-18

Also available in Rigid and Extra Rigid.

## Resin 8 Colors



Hu-Friedy's Resin 8 Colors feature lightweight colour handles for easy identification combined with exclusive EverEdge Technology for the quality and performance.

Eight distinct colours allow clinicians to quickly and easily select instrument types that follow a colour-coding system based on the specific area of the mouth where it is used.



HF-SG1-2C8E2 1/2 Grey



HF-SG3-4C8E2 3/4 Grey



HF-SG5-6C8E2 5/6 Yellow



HF-SG7-8C8E2 7/8 Green



HF-SG9-10C8E2 9/10 Green



HF-SG11-12C8E2 11/12 Orange



HF-SG13-14C8E2 13/14 Blue



HF-SG15-16C8E2 15/16 Orange



HF-SG17-18C8E2 17/18 Blue



HF-SG11-14C8E2 11/14 Red




HF-SG12-13C8E2 12/13 Red

## After Five



Patterned on basic Gracey designs but with terminal shank 3mm longer - for deeper periodontal pockets - 5mm or more. Finer blade for smooth subgingival insertion.



HF-SRPG1-2



HF-SRPG7-8

# Hygiene



HF-SRPG1-2C8E2 1/2 Grey



HF-SRPG3-4C8E2 3/4 Grey



HF-SRPG7-8C8E2 7/8 Green



HF-SRP11-12C8E2 11/12 Orange



HF-SRP13-14C8E2 13/14 Blue



HF-SRPG11-12



HF-SRPG13-14



HF-SRPG15-16

Also available in Rigid & Extra Rigid Shank Options

## Micro Mini - Five



How the best perform

Ultra slend blades, extended shank length, reduced blade length, precise bends.



HF-SMS1-29



HF-SMS7-89



HF-SMS11-129



HF-SMS13-149

## Mini Five



How the best perform

Same shank length as AFTER FIVE series but with reduced blade length for deep, narrow pockets.



HF-SAS1-2C8E2 1/2 Grey



HF-SAS3-4C8E2 3/4 Grey



HF-SAS5-6C8E2 5/6 Yellow



HF-SAS7-8C8E2 7/8 Green



HF-SAS11-12C8E2 11/12 Orange



HF-SAS13-14C8E2 13/14 Blue



HF-SAS15-16C8E2 15/16 Orange



HF-SAS1-29



HF-SAS3-49



HF-SAS5-69



HF-SAS7-89



HF-SAS11-129



HF-SAS13-149




HF-SAS15-169

Also available in Gracey and Rigid.

## Gracey Currettes



### Air Hollow



HS-9008047 1/2

HS-9008050 3/4

HS-9008051 5/6

HS-9008052 7/8

HS-9008108 9/10

HS-9008048 11/12

HS-9008049 13/14

HS-9008109 15/16

HS-9008110 17/18

### MaxiGrip Hollow

HS-1048375 1/2	HS-1048381 11/12
HS-1048377 3/4	HS-1048399 13/14
HS-1049091 5/6	HS-1048383 15/16
HS-1048378 7/8	HS-1048384 17/18
HS-1048380 9/10	

### Anatomic Colours




Coloured handles for clear identification. Special stainless steel working tips with satin finish to avoid reflections.



AE-DB380 1/2 Grey	AE-DB384 9/10 Turquoise
AE-DB381 3/4 Black	AE-DB385 11/12 Purple
AE-DB382 5/6 Yellow	AE-DB386 13/14 Blue
AE-DB383 7/8 Green	AE-DB387 15/16 Black

## Universal Currettes

### Barnhart



HF-SBH5-6

HF-SBH1-2RC8E2 1/2 Pink Rigid

HF-SBH5-6C8E2 5/6 Pink


### Columbia



AE-DB362 2R/2L 170mm, Ergonomic Grey Handle

AE-DB363 4R/4L 170mm, Ergonomic Black Handle

### Columbia



HF-SC2R-2LC8E2 2R/2L Pink

HF-SC4R-4LC8E2 4R/4L Pink

HF-SC13-14C8E2 513/14 Pink

HF-SC2R-2L Anteriors

HF-SC4R-4L Posteriors

HF-SC13-14 Universal Application

## McCall

**Hu-Friedy**

How the best perform



HF-SM13-14C8E2 13/14 Pink



HF-SM13-14SC8E2 13/14 Pink



HF-SM17-18C8E2 17/18 Pink

## Langer

**Hu-Friedy**

How the best perform



HF-SL1-2C8E2 1/2 Pink



HF-SL3-4C8E2 3/4 Pink



HF-SL5-6C8E2 5/6 Pink



HF-SL17-18C8E2 17/18 Pink

## Younger-Good

**Hu-Friedy**

How the best perform



HF-SYG7-8C8E2 7/8 Pink

## Vision Curettes

**Hu-Friedy**

How the best perform

An enhancement to the original Gracey design. Increased blade curvature and reduced blade length.



HF-SGC1-2 Anterior Premolar



HF-SGC11-12 Mesial Surfaces of Molars



HF-SGC13-14 Distal Surfaces of Molars



HF-SGCSUB-0 Sub Zero - For Anterior Premolar Surfaces

## Sickle Scalers

### Sickle

**Hu-Friedy**

How the best perform

Particularly for interproximal areas. Slight blade angulation provides accessibility to all mesial and distal surfaces.



HF-SH6-7H2

Unique angulation for less accessible posterior interproximal areas.



HF-SM13-14 McCall



HF-SCNEVI1 Nevi 1 Anterior



HF-SCNEVI2 Nevi 2 Posterior



HF-SCNEVI39E2 Nevi 3 Posterior



HF-SCNEVI49E2 Nevi 4 Posterior



HF-SCNEVI1C8E2 Nevi 1 Anterior Purple



HF-SNEV1H5C8E2 Nevi 1 Anterior Purple



HF-SCNEVI2C8E2 Nevi 2 Posterior Red



HF-SCNEVI3C8E2 Nevi 3 Posterior Red



HF-SCNEVI4C8E2 Nevi 4 Posterior Red



HF-SH5-33C8E2 SH5/33 Anterior Purple



HF-SH6-7C8E2 SH6/7 Anterior Purple



HF-SJ34-35C8E2 SJ34/35 Posterior Red



HF-SM0-00C8E2 SM0/00 Anterior Purple



HF-SN135C8E2 SN135 Anterior Purple



HF-SN137C8E2 SN137 Anterior Purple



HF-S204SC8E2 S204S Posterior Red



HF-S204SDC8E2 S204SD Posterior Red



HF-SUI1530C8E2 SUI15/13 Anterior Purple



HF-SUI15-33C8E2 SUI15/33 Anterior Purple

## Sickle



Anatomic colours / yellow handle. H6/7.



AE-DB370



AE-DB512 Cumine



AE-DB515 Mitchell

## Beach Scaler #3



Sharp and durable blades are all that are generally needed to treat all surfaces. Short light-weight handle for enhanced tactile sensation and reduced finger fatigue.

M-101133-3

## Implant Scalers

### Titanium Implant Scalers



How the best perform



PERFORM AT YOUR BEST with Hu-Friedy's new, state of the art line of Titanium Implant Scalers. Made from the same titanium alloy as implants and abutments, these scalers are expertly designed and manufactured for implant maintenance, debridement, biofilm removal and can be used both Supra and Subgingivally.

HF-IMP204SDT	204SD
HF-IMPLG1-2T	Langer 1/2
HF-IMP1-2T	Mini Five 1/2
HF-IMP11-12T	Mini Five 11/12
HF-IMP13-14T	Mini Five 13/14

## Implacare II



How the best perform

Implacare is made of PLASTEEL, a high grade, unfilled resin that causes the least amount of alteration to implant abutment surfaces. 20% thinner tip design which allows for improved access to implant abutment surfaces and allows for more effective maintenance care. Two new universal curette designs are now available. Barnhart 5/6 and Langer 1/2 provide additional curette design options so you can select the instrument best suited for your patient's needs.



### Implacare II Tips (12 Pairs) Including Handle Kit

HF-IMPL4R-4L6 4R/4L6  
 HF-IMPL 204S6  
 HF-IMPLH6-76 204S6 H6/76  
 HF-IMPLSORT6 Assorted  
 HF-IMPLINTRO6 Implacare II Introductory Kit  
 HF-IMPCKITII Implant II Maintenance Starter Kit

### Implacare II Refill Tips (12 Pairs)

HF-IMPL4R-4L 4R/4L6  
 HF-IMPL204S 204S6  
 HF-IMPLH6-76 H6/76  
 HF-IMPLSORT Assorted

## Double Ended Implant Scalers




### Implant Scaler Barnhart 5/6

HS-9007561 MaxiGrip Handle




### Implant Scaler 204

HS-9007564 MaxiGrip Handle



### Implant Scaler 4R/4L

HS-9007562 MaxiGrip Handle



### Implant Scaler Langer 1/2

HS-9007565 MaxiGrip Handle



### Implant Scaler H6/H7

HS-9007563 MaxiGrip Handle



## Periodontal Scalers



### Single Ended U15

HS-9008059 Air Handle



### Double Ended 204

HS-9008053 Air Handle  
 HS-1048370 MaxiGrip Handle



### Double Ended H6 H7

HS-9008055 Air Handle  
 HS-9792182 Solid Round Handle



# Operative Instruments

## Amalgam Carriers

### Del-ject Plastic Amalgam Carriers



#### Features:

- Grey
- Autoclavable to 130 Degs C
- Disassemble/unscrew before autoclaving
- Light weight: only 13 g




DJ-355036 Curved	DJ-355029 Tip Straight
DJ-355037 Straight	DJ-355027 Plunger Nylon
DJ-355028 Tip Curved	DJ-355026 Spring

### Amalgam Carriers



How the best perform

Syncote coated barrel ensures years of clog-free use. Improved, durable spring. Cold sterilisation not recommended, but fully autoclavable. Guaranteed for 12 months against clogging and defects in material and workmanship.



#### Single Ended

HF-AC5101 Mini (1.5mm)
HF-AC5102 Regular (2.0mm)
HF-AC5103 Large (2.8mm)



#### Double Ended

HF-AC5201 Mini/Regular
HF-AC5202 Regular/Large
HF-AC5203 Regular/Jumbo

#### MAINTENANCE TIP


Promptly remove excess material before sterilising with autoclave, chemiclave, or dry heat. Cold sterilisation solutions are NOT recommended. They contain chemicals that may adversely affect the performance of the CF II Carrier.

## Amalgam Packers

### Amalgam Packers



How the best perform



## Burnishers

### Burnishers



How the best perform

Designed to condense, smooth and polish amalgam.



### Double Ended Burnishers



#### Amalgam Burnisher, Air

HS-9008062 Large
HS-9008063 Medium
HS-9008064 Small



#### Ball Burnisher XL (2.7/3.2)

HS-9008069 Air Handle
HS-9001654 Octagonal Solid





# Operative Instruments

## Ball Burnisher L (2.0/2.5)

HS-9008070 Air Handle  
HS-9001673 Octagonal Solid



## Ball Burnisher M (1.2/1.6)

HS-9008071 Air Handle  
HS-9001672 Octagonal Solid



## Ball Burnisher S (0.8/1.0)

HS-9008072 Air Handle  
HS-9001671 Octagonal Solid



## Amalgam Burnisher Hollenback #1

HS-9008076 Air Handle  
HS-1012851 Round Hollow



## Amalgam Burnisher Hollenback #2

HS-9008077 Air Handle  
HS-1010579 Round Hollow



## Carvers

### Carvers



How the best perform



HF-CVIPC Interproximal



HF-CVW18 8 Wiland, Minnesota, IPC1



HF-CVHL1-2 1-2 Hollenback



HF-CVHL3 Hollenback 3



HF-CVWL1 Walls #1



HF-CVWL3 Walls #3



HF-CVWR2 Ward #2



HF-CD89-92 Cleoid/Discoïd

## Carvers



### Plastic Filling Instrument 49 DE

HS-9008096 Air Handle



### Plastic Filling Instrument 154 DE

HS-9008082 Air Handle



### Plastic Filling Instrument 155 HOPSON

HS-9008089 Air Handle



### Plastic Filling Instrument 156 HOPSON

HS-9008090 Air Handle



### Plastic Filling Instrument 178 DE

HS-9008083 Air Handle



### Plastic Filling Instrument 179 HOPSON

HS-9008091 Air Handle



### Plastic Filling Instrument CaOH, SE

HS-9008023 Air Handle



### Plastic Filling Instrument Heidermann M#1

HS-9008079 Air Handle  
HS-9002047 Octagonal Solid




## Aesthetic Carving Knives

### Aesthetic Carving Knives



How the best perform

Thin blade design for carving all restorative materials.  
Maintains sharpness even after autoclaving.



HF-CR21 For Posterior Teeth  
HF-CR20 For Anterior Teeth

# Operative Instruments

## Calcium Hydroxide Placement

### Calcium Hydroxide Placement



How the best perform

HF-PICH



## Modified Cleoids

### 1/ Hollenback /2 Modified Cleoid



How the best perform

HF-CVH1-MC2 Small  
HF-CVH1-MC3 Large



### Beach Amalgam Carver



Double ended "fan" end ideal for hoeing occlusal anatomy of lingual aspect, while other end is designed for occlusal anatomy of buccal aspect and interproximal carving. Light touch, slim diameter, knurled for positive finger control.

M-101200

#### DID YOU KNOW?

CLEOID-spade shape with point for carving occlusal anatomy.  
DISCOID-round shape for reducing excess material.

## Cement Spatulas



How the best perform



Combined mini spatula/placement instrument.

HF-SP6061



Combined long cement spatula/scoop loading instrument.

HF-CSNT5

### Cement Spatula L 7.5mm DE

HS-9008029 Air Handle  
HS-9001666 Octagonal Solid




### Cement Spatula M 6mm DE

HS-9008030 Air Handle  
HS-9001665 Octagonal Solid



### Cement Spatula S 5mm DE

HS-9008031 Air Handle  
HS-9001664 Octagonal Solid



AE-DF163


## Chisels

### Chisels



How the best perform

For use in planing and cleaving enamel.



HF-CP1-2 Wedelstaedt 1/2  
HF-CP11-12 Binangle 15-8-8

## Composite Brush

### Brush Handle



How the best perform



Autoclavable. Double ended - Ideal for fissure sealant and acid etch procedures.

HF-HCB1 (Brushes not included)

### Composite Manipulation Brush



GC-0104F4400000 Handle Pack 2  
GC-0104F4380000 Tips Pack 10

# Operative Instruments

## Composite Instruments

### Composculp by Dr. Dietschi



How the best perform



**Accurate Restorations** - Highly polished stainless steel working ends resist sticking to composite materials.

**Durable and Lasting** - Working ends made with Immunity Steel®, an optimal blend of carbon and chromium providing strength and resisting corrosion.

**Easy Cleaning** - The smooth handle's sleek design aids in preventing the build up of composite material.

- HF-PFIDD1/2 1/2 Dietschi Smooth Handle
- HF-PFIDD7/8 7/8 Dietschi Smooth Handle
- HF-PFIDD9/10 9/10 Dietschi Smooth Handle
- HF-PFIDDASST Dietschi Kit 5 Instruments & IMS Cassette

### Flexi Thins



How the best perform

**High Polish, Flexibility and Thinness.**

Suitable for composite resin, glass ionomer or acrylic temporary restorations.



HF-CIGFT1 # 1



HF-CIGFT3 # 3 (Extra Flex) for Class 3 & 4 Restorations



HF-CIGFTMINI1 Mini 1



HF-CIGFTMINI3 Mini 3

### Felt Composite Instruments



How the best perform

Lightweight, non stick anodised aluminium black colour instruments. Ultra-sonic cleaning not recommended.



HF-CI0130 # 4 Medium Blades For Class 3 & 4 Restorations



HF-CI0140 # 6 Smaller Blades For Class 1 & 3 Restorations with Limited Access

### Satin Steel XTS



How the best perform

Aluminium Titanium Nitride Coating resists scratches and sticking. Unique black finish offers enhanced contrast between instrument, tooth structure and composite. Smooth, large, lightweight handle design.



HF-TNCIGFT1 #1 Goldstein Flexi Thin



HF-TNCIGFT2 #2 Goldstein Flexi Thin  
HF-TNCIGFT3 #3 Goldstein Flexi Thin  
HF-TNCIGFT4 #4 Goldstein Flexi Thin



- HF-TNPFI8A # 8A PFI
- HF-TNPFI6 # 6 PFI
- HF-TNPFIW3 # W3 PFI
- HF-TNCVIPC Interproximal Carver
- HF-TNPCCI Duckhead

### Composite Instrument



#### Composite Instrument Air, Anodised

- HS-9008065 F4
- HS-9008066 F5
- HS-9008067 F6



# Operative Instruments

## CompoRoller

KAVO  Kerr

CompoRoller is an innovative composite modelling instrument that is designed to give you complete control in layering and contouring a direct composite restoration to its final form.



**HAW-5300 Assorted Kit**  
Contains: 1 CompoRoller Handle, 49 Tips (7 each of Cylindrical, Conical, Cylindrical 3mm, Disc, Point, Oval, Spatula)



### Refills Pack 2

HAW-5301 Handle

### Pack of 100

HAW-5302 Cylindrical Tip 4mm

HAW-5303 Conical Tip

HAW-5304 Cylindrical Tip 3mm

HAW-5305 Disc Shape Tip

HAW-5306 Point Shape Tip

HAW-5307 Oval Shape Tip

HAW-5308 Spatula Tip

## Composite Instruments

AESCULAP 

### Features:

- Ergonomic plastic handles
- Gold plated working tips for easy identification



For exact moulding of the cusp

AE-DF256



For adapting composite to edges.

AE-DF259



For removing excess material.

AE-DF258



For application of composite in the cavity.

AE-DF270



For application of composite in the cavity.


AE-DF086

## Cord Packers

### Cord Packers

**Hu-Friedy**

How the best perform



HF-GCP113 Serrated  
HF-GCPG7 Guyer Serrated  
HF-GCPYD1 Yardley



Thin blade and rounded contour for atraumatic and accurate cord placement. Bilateral notch allows placement in limited access areas.

HF-GCPBN1

### Cord Packer

**HENRY SCHEIN**  
CORPORATE BRAND

### Cord Packer Serrated

HS-9008033 Air Handle



### Fischer's Ultrapak Packer

**ULTRADENT**  
PRODUCTS, INC.



Thin serrated packing instrument.

ULT-0171 Small  
ULT-0172 Regular

## Crown Removers

### Christensen Crown Remover

**Hu-Friedy**

How the best perform



Christensen Crown Remover allows for easy crown removal after sectioning of the crown with bur (e.g. Komet H34). Two notches at the tip provide a secure grip and as the instrument is turned the notches engage the crown and break the seal of the cement.

HF-CRCH1 Straight Shank For Anteriors  
HF-CRCH2 Right Angled Shank For Posteriors

# Operative Instruments

## Crown Remover Pliers

Designed for easy and safe removal of direct restorations during trial fitting, temporary crowns and metal cores. Features a rasped tip for positive grip on temporary crowns or rubber pads for careful removal of final restorations.




GC-KYPLIERSET Crown Removal Pliers  
GC-KYRUBBERPACK Replacement Rubber Pads Pack 50

## Crown Removers



AE-DP770 Crown and Bridge Remover With 3 Points  
AE-DP781 Handle  
AE-DP783 Tip  
AE-DP784 Tip  
AE-DP787 Tip

## Excavators



HF-EXC246 #246

HF-EXC127-8 #127-8

HF-EXC220-1 Darby Perry #220-1

HF-EXC17 #17

HF-EXC31L Endodontic #31L

HF-EXC33L #33L

HF-EXCGL2 Glick #2

## Double Ended Excavators



### S (1.0mm)

HS-9003423 CT Handle  
HS-9002044 Octagonal Solid

### M (1.4mm)

HS-9003422 CT Handle  
HS-9008045 Air Handle  
HS-9002043 Octagonal Solid

### L (1.7mm)

HS-9003421 CT Handle  
HS-9008044 Air Handle  
HS-9002042 Octagonal Solid

### #125/126

HS-9008038 Air Handle

### #127/128

HS-9008039 Air Handle

### #129/130

HS-9008040 Air Handle

### #131/132

HS-9008041 Air Handle



## Gingival Retractors

### Meinershagen



For mandibular incisors.


HF-GRM1

For maxillary and mandibular pre-molars and canines. Also maxillary lateral incisors.

HF-GRM2

For maxillary central incisors and wide canines.

HF-GRM3



# Operative Instruments

## Hatchets

For preparation of retentive areas, sharpening internal line angles and removing hard caries.



HF-CP8-9H 10-7-14

## Margin Trimmers

To produce correct bevel on gingival enamel.



HF-MT28 #28



HF-MT77-78 #77-78



HF-MT79-80 #79-80



HF-MT29 #29

## Plastic Filling Instruments



HF-PFI1 #1



HF-PFIHF2 #2




HF-PFIA6 #A6



HF-PFI179 #179



HF-PFIG4-5 Gregg



AE-DE117 170mm Double Ended



AE-DE120 170mm Double Ended

# Surgical

## Bone Files

### Bone Files



HF-BF12 Howard



HF-BF67 Howard



HF-BF45 Miller



HF-BF2X Miller/Coburn (Cross Cut)



HF-BF1X Miller-Colburn 1X



HF-BF3X Miller-Colburn 3X

### Miller



AE-D0771



AE-D0772



AE-D0773

## Chisels



HF-CC2 Chandler



HF-CG2 Gardener 2



HF-CW1 Wakefield 1W



HF-CW2 Wakefield 2W



HF-CC7 Chandler 7



HF-C36-379 Back-Action EverEdge



HF-S13K-13KL9 #13K/13K EverEdge



HF-S13K/TG9 #13K/TG DE EverEdge



HF-C01



HF-C02



HF-C04



AE-D0711 Lucas 160mm



# Surgical

## Curettes

**Hu-Friedy**  
How the best perform



HF-SKRA39 Kramer EverEdge  
HF-SPR1-29 Prichard EverEdge



Larger and heavier curettes for the removal of granulation tissue and tenacious subgingival deposits.

HF-SK8-9

HF-CL84 Lucas 84

HF-CL85 Lucas 85

HF-CL86 Lucas 86

HF-CM10 Miller 10

HF-CM11 Miller 11

HF-CM12 Miller 12

HF-CM9 Miller 9

### Black Line


**Hu-Friedy**  
How the best perform



HF-CM2-4X 2/4 Molt  
HF-CM2-4SX 2/4 Molt Short Handle



HF-CL86X 86 Lucas  
HF-CL86SX 86 Lucas Short Handle



HF-CL87X 87 Lucas  
HF-CL87SX 87 Lucas Short Handle

**AESCLAP**  
Aesculap is a B. Braun company



AE-D0608 Lucas



AE-D0609 Lucas



AE-D0651 Jansen



AE-D0655 Hemingway



AE-DX220 Lucas Ergonomic Handle

## Elevators

**Hu-Friedy**  
How the best perform




HF-EWJS Warwick James - Straight



HF-EWJL Warwick James - Left



HF-EWJR Warwick James - Right



HF-EA1C Coupland A1C



HF-EA2C Coupland A2C



HF-EA3C Coupland A3C



HF-E44 Cryer (Small) 44



HF-E45 Cryer (Small) 45



# Surgical

HF-E1C Cryer 1C



HF-E2C Coupland 2C



HF-E3C Coupland 3C



HF-E21 Mini Cryer 21



HF-E22 Mini Cryer 22



AE-DL236 Heidbrink



AE-DL237 Heidbrink



AE-DL221 Cryer, Wide



AE-DL222 Cryer, Wide



AE-DL15 Flohr



AE-DL52 Warwick James - Straight



AE-DL263 Warwick James - Left



AE-DL264 Warwick James - Right



AE-DL46 Coupland



AE-DL47 Coupland



AE-DL48 Coupland



AE-DL211 Cryer, Narrow



AE-DL212 Cryer, Narrow



AE-DL2385 Flohr



## Elevators



### Luxator Perioste



Preserves bone integrity and performs less traumatic extractions. Highly acclaimed, renowned and proven functions.


- |                    |                    |
|--------------------|--------------------|
| DR-1S 1mm Straight | DR-3C 3mm Curved   |
| DR-2S 2mm Straight | DR-5S 5mm Straight |
| DR-3S 3mm Straight | DR-5C 5mm Curved   |

## Luxating Elevators

**Hu-Friedy**  
How the best perform



HF-ELKIT Kit Contains: 4 x Luxating Elevators 3C, 3S, 5C, 5S and Sharpening Stone




HF-EL3C Luxating Elevator Curved 3mm  
HF-EL3S Luxating Elevator Straight 3mm  
HF-EL5C Luxating Elevator Curved 5mm  
HF-EL5S Luxating Elevator Straight 5mm

HF-EL3CSM 3mm Curved Small Handle  
HF-EL3SSM 3mm Straight Small Handle  
HF-EL4C 4mm Curved  
HF-EL4S 4mm Straight  
HF-EL5CSM 5mm Curved Small Handle  
HF-EL5SSM 5mm Straight Small Handle



HF-SSLUX Sharpening Stone

## Black Line 2014



HF-EL3SX Straight 3mm    HF-EL5CX Curved 5mm  
HF-EL3CX Curved 3mm    HF-EL5SX Straight 5mm

## Luxator Hybrids

**Hu-Friedy**  
How the best perform

Luxator Hybrids were designed with thin, precisely angled tips to aide clinicians in preserving the buccal plate during atraumatic extractions. Enhanced oral adaptation with 7 manufactured tip options. Easily severs the periodontal ligament with finely honed, sharp cutting edges.



HF-EPTSSP Straight Spade



HF-EPTSMS Small Straight



HF-EPTSMC Small Curved



HF-EPTSMD Small Distal



HF-EPTSMM Small Mesial



HF-EPTLGD Large Distal



HF-EPTLGM Large Mesial

HF-EPTKIT Luxating Hybrid Kit (7) Includes:  
SSKC Sharpening Stone & IMS Cassette IM5080

## Periosteal Elevators



How the best perform

### 9 Molt

For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 4.3mm/8.4mm

HF-P9



### 9A Allen

Hole allows tissue to be held during suturing.

HF-P9A



### Grisdale

HF-PGR1



### 9H Howard

HF-P9H



### 3 Mead

HF-P3



### 3 Pritchard

Featuring two blades, one large flat blade and one narrow sharp blade on round tip. 4.7mm/11.0mm

HF-PPR3



### 149 Periosteal

HF-P149



### Goldman-Fox #14

HF-P14



### #23 Seldin

HF-P23

## Black Line Periosteal Elevators



How the best perform

### Buser

HF-PPBUSERX Buser  
HF-PPBUSERSX Buser Short Handle



### Allen

HF-PPAELX Allen Elevator  
HF-PPAELSX Allen Short Handle  
HF-PPAELAX Allen Elevator Anterior  
HF-PPAELASX Allen Elevator Anterior



### Prichard

HF-PPR3X PR-3 Prichard  
HF-PPR3SX PR-3 Prichard Short Handle



### Molt

HF-P9X 9 Molt  
HF-P9SX 9 Molt Short Handle



### 152K-N

HF-PKN152X 152 K-N  
HF-PKN152SX 152 K-N Short Handle

## Root Picks



Flohr - For upper roots, spoon shaped.

AE-DL15

Heidbrink - For lower roots

AE-DL236

Heidbrink - For lower roots

AE-DL237

## Heidbrink

For removal of roots in restricted areas.



How the best perform

HF-EHB1 #1



HF-EHB2 #2



HF-EHB3 #3



HF-EHB13-14 Heidbrink 13/14



## Root Tip Elevators



How the best perform

HF-EH02 Howard #2



HF-EH03 Howard #3



West Apical - for removal of extremely small pieces.

HF-EW2 West #2



HF-EW3 West #3



HF-EW4 West #4



HF-EW5 West #5



HF-EW6 West #6



## Extraction Forceps



### Upper Incisors & Canines

HS-9001047 EF7



HS-9001046 EF2



### Upper Bicuspids

HS-9001047 EF7



### Upper Molars

HS-9001042 EF17 Right Side



HS-9001043 EF18 Left Side



HS-9001052 E67A 3rd Molar



### Upper Roots

HS-9001044 EF51A



HS-9002016 EF51A Fine



### Lower Bicuspid/Cuspid

HS-9001050 EF13



### Lower Molars

HS-9001048 EF22 Left/Right



HS-9001051 EF79 3rd Molar Left



### Lower Incisal Roots

HS-9001049 EF33A



# Surgical



Anatomically shaped for a safe and reliable grip on teeth and roots. Slender shape and curved handles for optimum handling. Guaranteed precision instruments.

## Upper Incisors & Canines

AE-DG5 Fig 2



AE-DG6 Fig 1



## Upper Canines with Short Beaks

AE-DG32 Fig 107



## Upper Bicuspid

AE-DG101 Fig 136



AE-DG107 Fig 7



## Upper Molars - Right Side

AE-DG117 Fig 17



## Upper Molars - Left Side

AE-DG118 Fig 18



## Upper Wisdom Teeth - Narrow

AE-DG250 Fig 67



## Upper Wisdom Teeth - Wide

AE-DG252 Fig 67A



## Upper Anterior Roots

AE-DG305 Fig 29



AE-DG320 Fig 30



AE-DG330 Fig 76



## Upper Roots, Short, Narrow Beaks

AE-DG351 Fig 51



## Upper Interior Roots

AE-DG300 Kaiser



## Lower Incisors & Bicuspid

AE-DG413 Fig 13



## Lower Canines & Bicuspid

AE-DG440 Fig 75



## Lower Molars

AE-DG555 Fig 73



AE-DG522 Fig 22



AE-DG570 Fig 86



# Surgical

## Lower Roots & Incisors

AE-DG700 Fig 74N



## Lower Roots - 4.5mm Wide Beaks

AE-DG704 Fig 74



## Lower Roots - Narrow Beaks

AE-DG732 Fig 33A



## Lower Roots - Wide Beaks

AE-DG733 Fig 33



## Children's Patterns

AESCLAP<sup>®</sup>  
Knowledge - A & B team company

### Upper Incisors

AE-DK1 Fig 37



### Upper Molars Right

AE-DK17 Fig 39R



### Upper Molars Left

AE-DK18 Fig 39L



### Lower Molars

AE-DK60 Fig 22S



## Lower Roots & Incisors

AE-DK70 Fig 162



## Atraumair Forceps

**Hu-Friedy**

How the best perform

Atraumair Forceps feature enhanced beak geometries that allow for optimal access and grip when extracting a tooth. The matte finished handle reduces light reflection, while the unique hole design reduces the weight of the handle.

HF-F35XS #35



HF-F36XS #36



HF-F13XS #13



HF-F1XS 1 Standard



HF-F23XS 23 Cowhorn



HF-F67XS #67



HF-FX79XS #79



HF-FX73XS #73



HF-FX51XS #51 Root Forceps



HF-F45XS #45 Serrated



HF-FAF35XS #35



# Surgical

HF-FAF36XS #36



HF-FAF1XS 1 Apical Forceps

## Forceps



### Upper Molars

HF-F18LS 18LS Harris



HF-F18RS 18R Harris



HF-F90 90 Cook



HF-FX67A 67A European Serrated



HF-FX17 17 European Serrated



HF-FX18 18 European Serrated



### Lower Molars

HF-F16 16 Cowhorn



HF-F23 23 Cowhorn



HF-F217 217



### Lower Molars

HF-F222 222



HF-FX22 22 European Serrated



HF-FX73 73 European Serrated



HF-FX79 79 European Serrated



HF-FMD4 MD4 Mead Serrated



### Upper Anteriors & Molars

HF-F1 1 Serrated



HF-FX1 1 European Serrated



HF-FX7 7 European Serrated



HF-F99C 99C Kells



### Upper Molars & Pre-Molars

HF-F150 150 Cryer Universal



HF-FMD1 MD1 Mead Serrated



### Lower Anteriors & Pre-Molars

HF-FX13 13 European Serrated



HF-FMD3 MD3 Mead Serrated



# Surgical

## Upper Roots

HF-F300 300 Serrated



HF-FX51 51 European Serrated



HF-FAFX74N AF74N European Serrated



## Lower Roots

HF-FX74N 74N European Narrow Beaks Serrated



HF-FX46L 46L European Serrated



HF-FX233 233 European Serrated



HF-FX74 74 European Serrated



## Paedodontic

HF-F151S 151S Pedo



HF-F39 39 Pedo



HF-F40 40 Pedo



HF-F150S 150S Pedo



## Rainbow Paedodontic Forceps

HF-F150SR F150SR




HF-F151SR F151SR



## Surgical Forceps

**HENRY SCHEIN®**  
CORPORATE BRAND



### Children's Forceps with Spring

HS-9002017 #1 For Upper Incisors  
HS-9002018 #3 For Upper Molars  
HS-9002019 #5 For Lower Incisors  
HS-9002020 #6 For Lower Molars

### Forceps Dunhill Curved

HS-9008670 12.5cm



### Stieglitz Root Tip Forceps

HS-1014772 45°  
HS-1015107 90°



## Haemostats

**Hu-Friedy**  
How the best perform



HF-H1 Kelly - Straight  
HF-H2 Kelly - Curved  
HF-H3 Halsted - Mosquito Curved  
HF-H4 Straight



# Surgical

**HENRY SCHEIN®**  
CORPORATE BRAND



HS-1004862 Halsted Mosquito Curved 12.5cm  
HS-1004662 Halsted Mosquito Straight 12.5cm




HS-9003589 Rongeur Mini-Friedman 12.5cm  
HS-9009616 Rongeur Friedman 14cm




HS-1000107 Backhaus Towel Clamp 9cm  
HS-1000341 Backhaus Towel Clamp 14cm


**AESCULAP®**  
HEINZ + HEIMANN



AE-BH108 Micro-Halsted 125mm  
AE-BH109 Micro-Halsted 125mm



Halsted-Mosquito 125mm  
AE-BH110  
AE-BH111



Crile 140mm  
AE-BH144  
AE-BH145

## Needle Holders

**Hu-Friedy**

How the best perform



HF-NH5020 Castro Straight 14cm/5.5"  
HF-NH5020R Castroviejo Straight Round 14cm/5.5



HF-NH5021 Castroviejo Curved 14cm/5.5"  
HF-NH5021M Castro Micro Curved 14cm/5.5"  
HF-NH5024 Castroviejo Straight 18cm/7"  
HF-NH5024R Castroviejo Straight Round 18cm/7"

# Surgical



HF-NH5037 Halsey 13cm/5"  
 HF-NH5038 Crile-Wood 15cm/6"  
 HF-NH5042 Mayo-Hegar 16cm/6.25"



HF-NH5044 Mayo-Hegar Delicate 18cm/7"  
 HF-NHM5026R Castro Micro Straight  
 Diamond Dusted Round 18cm/7"

**HENRY SCHEIN**  
 CORPORATE BRAND



HS-9006214 Mathieu TC Inserts 14cm  
 HS-9006215 Mathieu TC Inserts 17cm



HS-9008637 Castroviejo TC Inserts 14cm




HS-9002079 Mathieu 14cm  
 HS-9002080 Mathieu 17cm



HS-1008670 Mayo Hegar 14cm

**AESCULAP**  
 Aesculap - a B. Braun company



AE-BM235 Mayo-Hegar 150mm  
 AE-BM361 Mathieu 170mm  
 AE-BM13 Baby Crile-Wood with Carbide Inserts 150mm  
 AE-BM003 Durogrip 180mm with Tungsten Carbide Inserts

# Surgical

## Rongeurs

**AESCULAP**  
Aesculap - A B. Braun company

AE-F0415 Bleumenthal 155mm



AE-F0409 Luer-Friedmann 145mm



AE-F0411 Luer 150mm



AE-F0412 Luer 155mm



**Hu-Friedy**  
How the best perform

HF-R1 #1 16.5cm



HF-R1A Mead 1A



HF-R5S Cleveland 55



## Scissors


**Hu-Friedy**

How the best perform

HF-S5 Wagner Straight 11.5cm/4.5"



HF-S5001 Kelly Curved 16cm/6.25"



HF-S5054 Metzenbaum Straight Blunt 14.5cm/5.5"



HF-S5056 Metzenbaum Straight Pointed 14.5cm/5.5"



HF-S5057 Metzenbaum Curved Pointed 14.5cm/5.5"



HF-S5080 Goldman-Fox Straight 13cm/5"



HF-S5083 Iris Curved 11.5cm/4.5"





HF-S6 Wagner Curved 11.5cm/4.5"  
 HF-S8 Quinby 12.5cm/5"  
 HF-S9 Dean 16.5cm/6.5"

## Black Line Scissors



The Black Line Scissors feature a performance engineered coating for optimal edge retention, enriched contrast and visual acuity at the surgical site with reduced light reflection. The scissors are in a Super-Cut design featuring a razor sharp blade and serrated edge for reduced tissue slippage.



HF-S14SCX LaGrange #14 Double Curved Super Cut  
 HF-S16SCX Goldman-Fox #16 Curved Super Cut  
 HF-S17SCX Iris #17 Straight Delicate Super Cut  
 HF-S18SCX Iris Super Cut Curved



HF-S19SCX Metzenbaum Super Cut  
 HF-S20SCX Goldman-Fox Straight Super Cut  
 HF-SPVX Castroviejo Microsurgical



**Spencer Suture Scissor**

HS-9003089 Straight 90mm



**La Grange Gum Scissor**

HS-9002084 Curved 11.5cm



**Iris Scissors**

HS-1008820 Curved 11.5cm  
 HS-9002062 Curved Serrated 11.5cm



**Iris Scissors**

HS-1002767 Straight 11.5cm  
 HS-9002061 Straight Serrated 11.5cm



**Kelly Scissor**

HS-9002025 Curved 16cm



**Mayo Scissor Blunt/Blunt**

HS-1003525 Curved 14cm HS  
 HS-1003943 Curved 17cm HS



**Mayo Scissor Blunt/Blunt**

HS-1001106 Straight 14cm HS  
 HS-1001237 Straight 17cm HS



**Crown Scissor**

HS-9001998 Curved Serrated

# Surgical

## Crown Scissor

HS-9001997 Straight Serrated



## Goldman Fox Scissor

HS-9002083 Curved 13cm



**AESCLAP**  
Advanced Endoscopic Surgery

AE-BC110 Straight 110mm



AE-BC111 Curved 110mm



AE-BC140 Straight 130mm



AE-BC141 Curved 130mm



AE-BC606 Metzenbaum Disecting 180mm Curved




AE-BC601 Metzenbaum Disecting 145mm Straight



AE-BC802 Spencer 115mm



AE-BC154 La Grange 110mm  
 One Blade Serrated



AE-BC881 Bandage/Cloth/  
 Rubber Dam Scissors



## Scalpel Handles

**Hu-Friedy**

How the best perform



HF-10-130-03 #3



HF-10-130-05 #5 Straight Round



HF-10-130-05E European Style Round #5



HF-10-130-07 #7



HF-SHDPV Microsurgical Stainless Steel

Swann-Morton®



MX-8-136 Handle #3 2122 Grad. Stainless Steel Each



MX-8-137 Handle #4 9081 Grad. Stainless Steel Each

AESCLAP®  
Aesculap - A.B. Braun company



AE-BB73

## Scalpel - Complete

Swann-Morton®



MX-8-154 Disposable Scalpel Sterile Steel #11 0503 Pkt/10



MX-8-155 Disposable Scalpel Sterile Steel #12 0504 Pkt/10



MX-8-156 Disposable Scalpel Sterile Steel #15 0505 Pkt/10

AESCLAP®  
Aesculap - A.B. Braun company



Sterile scalpel with plastic handle and blade protector.

- Pack 10  
 AE-BA210 # 10    AE-BA212 # 12  
 AE-BA211 # 11    AE-BA215 # 15

## Scalpel Blades

Swann-Morton®



### Sterile Carbon Steel Box 100

- MX-8-114 Blade #10 0201  
 MX-8-116 Blade #11 0203  
 MX-8-117 Blade #12 0204  
 MX-8-118 Blade #15 0205

### Sterile Stainless Steel Box 100

- MX-8-128 Blade #11 0303  
 MX-8-129 Blade #12 0304  
 MX-8-130 Blade #15 0305



AESCLAP®  
Aesculap - A.B. Braun company

Sterile Packaged Blades Box 100

## Tissue Forceps

Hu-Friedy®

How the best perform



HF-TP31 Semkin-Taylor Straight 12.5cm/5"



HF-TP33 Semkin-Taylor Straight 1 x 2 12.5cm/5"



HF-TP34 Semkin-Taylor Curved 1 x 2 12.5cm/5"



HF-TP41 Adson Plain 12cm/4.75"



HF-TP42 Adson 1 x 2 12cm/4.75"



HF-TP43 Adson-Brown 7 x 7 12cm/4.75"

# Surgical



HF-TP5042 Adson 1 x 2 12cm/4.75"



HF-TP5070 General 15cm/6"



HF-TP5071 General 1 x 2 15cm/6"



HF-TP5090 DeBakey 15cm/6"



HF-TPG1 Gerald 1 x 2 Serrated/Straight 18cm/7"



AE-BD45 Standard 130mm



AE-BD700 Adson-Brown 120mm



AE-BD660 Gillies 155mm



AE-BD512 Adson 120mm

## Surgical Pliers



### Tissue Pliers 1 x 2

HS-9002009 13cm HS-1009871 15cm HS-9002010 18cm

### Adson 1 x 2

HS-1005891 12cm



## Pliers

### Crown Pliers



MD-0417 Crown Crimping Pliers



MD-0421 Small Crown Crimping Pliers



MD-0416 Crown Contouring Pliers

## Scissors



How the best perform

### Crown & Gold Scissors



HF-SCGC Curved 12cm



HF-SCGCP Curved 9cm Paediatric



HF-SCGS Straight 12cm

### Crown Scissors




Deluxe Curved Festooning.

MD-1202

# Surgical Accessories

## Mouth Props



**Hu-Friedy**

How the best perform

Markel Pack 2

HF-MP54 Pair Adult Latex-Free

HF-MP52 Pair Child Latex-Free

HF-MP48 Pair Child Small Latex-Free



**University of Minnesota**

HF-CRM 5.5" Long x 20.5mm Wide



**Cawood-Minnesota**

HF-CRM2 15cm Long



**Columbia**

HF-CRC 5.4" Long

Silicone mouth wedge with metal core. Particularly advantageous for longer implantations and also provide support for the patient.



**AESCLAP**  
Aesculap - A.B. Braun company

AE-DX113 Pair, Small

AE-DX114 Pair, Medium

AE-DX115 Pair, Large

**AESCLAP**  
Aesculap - A.B. Braun company

AE-BT20 125mm



**Roux**

AE-BT33 28 x 29mm



## Retractors



**Hu-Friedy**

How the best perform

**Adult Plastic**

Suitable for intra-oral photography. Cold sterilisation only.

HF-CRPA



**Child Plastic**

HF-CRPC 4.4" Long

**Langenbeck**

AE-BT320 210mm



**Rake**

AE-BT118 165mm



**Bowdler-Henry Rake**

AE-DO470 190mm



**Laster**

AE-DO472





# Surgical Accessories



## Free Access

### Features:

- Provides cheek and lip retraction without a second pair of hands.
- Ideal for crown and bridge impressions, bite registrations, implant or edentulous impressions and denture try-in.
- Provides unobstructed view during light-curing, bleaching, porcelain veneer try-in or tooth preparation.
- A convenient assistant during oral prophylaxis, periodontal charting, scaling and root planing.
- An essential tool for sealant applications, composite resin bonding procedures or fluoride applications.
- Maintains accessibility during air abrasion, micro-etching, periodontal surgery or laser techniques.
- Safe for autoclave, dry heat or chemical sterilisation.



**Helps guide the patient into maximum intercuspation (MI) when taking impressions.**

Ideal for impression taking, bite registration and denture try-in.



**Transparent material for aesthetic clinical photography.**

Provides cheek and lip retraction without a second pair of hands.

M-780455-L Large  
M-780455-S Small

## Optiview



KAVO ■ Kerr

OptiView is an innovative lip and cheek retractor that is designed to provide even retraction of lips and cheeks for great buccal and gingival access at the same time. Thanks to its anatomic shape, OptiView offers patient comfort even during longer dental treatments.

HAW-5500 Assorted Kit Contains: 2 Lip and Cheek Retractors with Mount Cushions  
HAW-5501 Cushion Refills Pack 10

## Sterilising Forceps



### Cheattle Sterilising Forceps

Angled to side, 280mm.

AE-JG820



### Sterilising Forceps

Straight, 200mm.

AE-JG822

# Specialty Instruments

## Endodontics

### Explorers & Excavators



How the best perform

HF-EXDG16 Explorer G16

HF-EXC31L Excavators #31L

HF-EXC33L #33L

HF-EXCGL2 Glicks Spoon #2

### Spreaders

ISO sizing with 5mm laser markings.

HF-RCS30  
HF-RCS40

HF-RCS3

### Nickel Titanium Spreaders

HF-RCSD11TNT

HF-RCSD11TSNT

HF-RCSD11T

HF-RCSD11TS  
HF-RCSD11

### Pluggers - Single Ended

ISO sizing with 5mm laser markings.

HF-RCP30  
HF-RCP40

### Pluggers - Double Ended

HF-RCP1-3

HF-RCP5-7

HF-RCP9-11  
HF-RCPGL1 Glick Blade #1

#### DID YOU KNOW?

ISO Pluggers have an ISO colour-coded band on the handle for ease of identification during usage: #30-Blue, #40-Black, #50-Yellow, #60-Blue

### Vertical Condensers

Sizes 8-12 for anterior & posterior, 1/2 sizes also available.

HF-RCP8

### Buchanan Hand Pluggers



Featuring both a tapered nickel-titanium end and a stainless steel end, Buchanan Hand Pluggers can be used in any common obturation technique, including warm and cold condensation. The pluggers are also ideal for backfilling canals. Use the NiTi end for finishing off condensation procedures in even the most radically curved canals. Either end condenses and controls the movement of gutta percha under pressure. In addition, the Buchanan Hand Pluggers are fabricated of the highest grade materials for excellent durability and sterilisation.

KE-9740058 Size #1  
KE-9740059 Size #2

# Specialty Instruments

## S Kondensers



### Features:


- Designed specifically for condensing warm gutta percha
- Colour-coded, light-weight anodised aluminium handles
- Tapered NiTi tips at one end for greater compressive strength and flexibility
- Stainless steel tips at the other end, ideal for condensing backfill material
- Three instrument sizes available individually or as a complete set: Black 40/80, Yellow 50/100 and Blue 60/120
- Ergonomic handle design with finger rests
- NiTi condenser end has laser depth markings at 5mm intervals to gauge depth in canal

M-BCONSET Set of 3 Instruments

## Nickel Titanium Spreader



Very flexible, eliminates pre-bending.



K-NTD11T

## Implantology

### Ergoplant Range



AE-DX505 Applicator for Bone Splitting Retractor



AE-DX220 Bone Curette (Lucas)



AE-DX555 Bone Material Applicator



AE-DX252 Bone File, Miller 195mm  
AE-DX253 Bone File, Miller 200mm



AE-DX600 Membrane Placement Instrument



AE-DX201 Raspatory (Buser)




AE-DX560 Sinus Elevator 200mm  
AE-DX561 Sinus Elevator 220mm



Titanium Implant Curette for Diamonds  
AE-DX360 4.0mm 170mm Long  
AE-DX361 5.5mm 170mm Long

AE-DX601R Bone Fixation Forcep



## Periosteal Elevator



How the best perform

PP Buser spear and paddle shaped ends are ideal for delicate tissue retraction.



HF-PPBUSER

## Periotomes



How the best perform

Thin sharp blades used to facilitate removal of teeth with minimal damage to surrounding alveolar bone by severing periodontal membrane.



HF-PT1 Posterior



HF-PT2 Anterior

HF-PEREUR6 Wise Double Ended  
HF-PT3 #3 Posterior Right Single-Ended  
HF-PT4 #4 Posterior Right Single-Ended  
HF-PT5 #5 Anterior Angled Single-Ended  
HF-PT6 #6 Anterior Straight Single-Ended  
HF-PJ01 Jovanovic 13K/TG 1mm/4mm

# Specialty Instruments

## Interchangeable Periostomes



### Features:

- Thin sharp blades available in four designs to easily sever the periodontal ligament
- Large diameter handle for enhanced ergonomics, tactile feel and optimal clinical performance
- Longevity and corrosion resistance afforded by handcrafting from the finest grade surgical stainless steel



HF-PTI7 Straight Wide



HF-PTI8 Straight Narrow



HF-PTI9 Angulated



HF-PTI10 Mesial/Distal Angulated



HF-PTIH Handle  
 HF-PTIKIT Interchangeable Periostome Kit  
 Contains: One Handle, Four Periostome Tip Designs  
 HF-IMS-1372TH IMS Bur Cushion

## Implant Chisels



HF-PPBUSER1-2 Buser 1/2



HF-PPBUSER2-3 Buser 2/3



HF-PPBUSER1-3 Buser 1/3



HF-PJ01 Jovanovic 13K/TG 1mm/4mm



HF-PJ02 Jovanovic 13K/TG 2mm/5mm

## Bone Pluggers



HF-PLGLABANCA LaBanca 3 x 1.5mm



HF-BPACK Bone Packer DE 5-7mm



HF-RASP Bone Rasp

## Bone Scoops



HF-CSEUWES6 Euwe Bone Scoop #6 Short Satin



HF-BSJ05 Jovanovic Bone Scoop 5mm



HF-IMP6522 Bone Graft Scoop



HF-CSLABANCA6 LaBanca Bone Curette 7mm

## Bone Preparation



HS-9001745 Periosteal Elevator Williger



HS-1003086 Periosteal Elevator Goldman Fox #14

## Implant Periosteals



HF-PGR1 Grisdale



HF-PPLABANCA 6#1 LaBanca

# Specialty Instruments



HF-PPLABANCA2 #2 LaBanca Modified



HF-PPBUSER6 Buser



HF-P9PA6 P9/PA



HF-P20SIC #20 Sicilia

## Periodontal Knives



How the best perform

### #15/16 Kirkland, EverEdge

For initial bevel incision for gingivectomy or gingivoplasty procedures. Especially effective in retro-molar region.

HF-KK15-169



### #1/2 Orban DE, EverEdge

Angled blade and contra-angled shank for posterior use to extract upper roots.

HF-K01/29



### 1/2 Allen Orban, EverEdge

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue.

HF-K012KP03A9



### #15/16 Kirkland

For initial incision for gingivectomy, especially effective in retro-molar region.

HF-HF-KK15-16



### #5/6 Buck

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue.

HF-KB5-6



### Orban

Angled blade and contra angled shank particularly suitable for posterior application. Useful for excising interproximal tissue.

HF-K01-2



### #1/2 Sanders

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue.

HF-KS1-2



## Black Line Periodontal Knives



How the best perform

### Tunneling Knives

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue. 1.8mm/2.4mm

HF-TKN1X 1 HF-TKN1SX 1 Short Handle



HF-TKN2X 2 HF-TKN2SX 2 Short Handle



### 1/2 Allen Orban

30% thinner than the standard Orban, features an angled blade and contra-angled shank for posterior use. Also useful for excising interproximal tissue.

HF-K012KP03AX HF-K012KP03RX



### Allen End-Cutting Intrасulcular

For incisions that remove or recontour soft tissue. Also useful to excise interproximal tissue.


HF-KPAX



### 15/16 Kirkland

For initial bevel incision for gingivectomy or gingivoplasty procedures. Especially effective in retro-molar region.

HF-KK15-16X 15/16 Kirkland  
HF-KK15-16SX 15/16 Kirkland Short Handle



# Specialty Instruments

## 1/2 Orban

Angled blade and contra-angled shank for posterior use. Also useful for excising interproximal tissue.

HF-KO1-2X

## GF-11 Goldman-Fox

Spear Point, similar to KGF8, except blades are offset for use in posterior region.

HF-KGF11X

## Implant Kits

### Dr Allen Oral Plastic Kit

The Dr. Allen Oral Plastic Kit provides precision micros-instruments specifically designed for today's minimally invasive soft tissue grafting procedures

HF-ALLENKIT

### Allen-Johnston Suture Removal Kit

Designed for removing sutures placed during periodontal plastic surgery procedures, when microsurgical techniques are utilised.

HF-AJSRKIT

## Osteotomes

HF-OSTPAL1 #1 2.7mm

HF-OSTPAL2 #2 3.2mm

HF-OSTPAL3 #3 3.7mm

HF-OSTPAL4 #4 4.2mm

**Hu-Friedy**  
How the best perform

**Hu-Friedy**  
How the best perform

HF-OSS6518S Straight 4mm

HF-OSS6520S Straight 7.5mm

## Sinus Retractors

**Hu-Friedy**  
How the best perform

HF-SGRT Sinus Retractor

HF-SRM10 Sinus Retractor Misch 10

## Sinus Lift Instruments

**Hu-Friedy**  
How the best perform

HF-IMP6345-71 #710 Sinus Lift Instrument 3mm

HF-SINC2 Sinus Curette #2

HF-IMP6577SC Kramer Nevins Acute

HF-IMP6578 Kramer Nevins Obtuse

HF-IMP300 #300 Sinus Lift Instrument

HF-PESIMI Sinus Lift Instrument 3mm S/E

# Laboratory Instruments

## Plaster Knives

### Plaster Knives



AE-DS60 Plaster Knife 160mm  
AE-DS62 Plaster Knife 180mm

## Wax Knives

### Wax Knives



AE-DS103 180mm

### Wax Knives



HF-LS7 7 Tapered  
HF-LS8R Standard Rigid

## Laboratory Spatulas



HF-WS7 Wax Spatula #7  
HF-LWSS Waxing Spoon/Spatula Large Black Plastic

## Mallets



AE-DF114 190mm / 165g



HF-MAL1 Mead 19cm / 7 1/4oz

## Waxing Instruments



**Pack Thomas**  
Used to flow on molten wax.

AE-D0472



HF-PKT2



HF-PKT3



HF-PKT4

# Orthodontic Instruments

## Orthodontic Pliers

 HENRY SCHEIN®  
CORPORATE BRAND

**Adams Pliers**

HS-9007625



**Crown & Band Contouring Pliers**

HS-9007633



**Hard Wire Cutter**

HS-9007602 15°



**Adhesive Removing Pliers**

HS-9007610



**Debonding Pliers**

HS-9007612



**Hollow Chop Contouring Pliers**

HS-9007624



**Angulated Debonding Pliers**

HS-9007613



**Detailing Step Pliers**

HS-9007629 1mm



**Hook Crimping Pliers**

HS-9007632 Angled



**Arch Bending Pliers**

HS-9007623



**Detailing Step Pliers**

HS-9007628 3/4mm



**How Pliers**

HS-9007609



**Bird Beak Pliers**

HS-9007615



**Distal End Cutter (Safety Hold)**

HS-9007600 Long Handle



**Jarabak Pliers**

HS-9007614



**Bird Beak with Cutter**

HS-9007616



**Distal End Cutter Safety Hold & Flush Cut**

HS-9007601



**Light Wire Bird Beak Pliers**

HS-9007617





# Orthodontic Instruments



**Light Wire Bird Beak with Groove**

HS-9007618



**Lingual Arch Forming Pliers**

HS-9007630



**Loop Tie Back Pliers**

HS-9007626



**Micro Cutter**

HS-9007604



**Micro Mini Pin & Ligature Cutter**

HS-9007605



**Micro Mini Pin & Ligature Cutter**

HS-9007606 Long Handle



**Micro Mini Pin & Ligature Cutter**

HS-9007607 15°



**Mini Three Jaw Pliers**

HS-9007620



**Multi-Use Cutter**

HS-9007603



**Nance Loop Forming Pliers**

HS-9007627



**Optical Pliers**

HS-9007622



**Posterior Band Removing Pliers**

HS-9007611 Long



**Slim Weingart Pliers**

HS-9007608



**Stop (V-Bend) Pliers**

HS-9007631



**Three Jaw Pliers**

HS-9007619



**Torquing Pliers Set**

HS-9007635



**Tweed Loop Forming Pliers**

HS-9007621

# Orthodontic Instruments

## Orthodontic Pliers



With Tungsten Carbide inserts. Max wire size .022 x .028"

AE-DP512 Hard Wire Cutter 125mm



With Tungsten Carbide inserts. Max wire size .014"

AE-DP511 Ligature Cutter



With Tungsten Carbide inserts. Universal pliers.

AE-DP036 Weingart Pliers



Max wire size .019 x .025"

AE-DP301 Loop Forming Pliers 135mm



Suitable for all orthodontic wires.

AE-DP335 De La ROSA, Arch Bending Plier 125mm



Steps of 3,4,5 and 6mm. For utility arches.

AE-DP333 Nance, Step Loupe Bending Plier 125mm



With Tungsten Carbide insert.

AE-DP744 Band Remover 140mm



Maximum wire size .036"

AE-DP316 Aderer Three-Finger Pliers 125mm



AE-DP1 Flat Nose 140mm



Bends wire up to 0.7mm.

AE-DP311 Angle Bird Beak



Bends wire up to 0.7mm

AE-DP312 Adams Bird Beak

# Orthodontic Instruments

## Wire Cutters



How the best perform

With our diamond-honed cutting edges and orbit-formed joint, you will experience sharp and accurately-aligned cutting edges for every procedure, every time. For each individual area of focus, we provide a specifically engineered cutter. To meet these stringent standards, all of our blades are hand finished to maximise cutting performance for the appropriate wire size.

### Universal Cut & Hold Distal End Cutter

When used intraorally it safely holds the cut portion to be discarded. Patient comfort is enhanced by the smooth, tapered tip. Wire Size  $\leq$  .012" to .021" x .025" (.30mm - .53mm x .64mm)

HF-678-101

### Universal Cut & Hold Distal End Cutter With Long Handle

Same as 678-101, with long handle. Wire Size  $\leq$  .012" to .021" x .025" (.30mm - .53mm x .64mm)

HF-678-101L

### Flush Cut Distal End Cutter

Wire Size  $\leq$  .021" x .025" (.53mm x .64mm)

HF-678-102

### Hard Wire Cutter, 15°

15° offset tip enhances access in difficult areas. Wire Size  $\leq$  .022" x .028" (.56mm x .71mm)

HF-678-103

### Hard Wire Cutter

Wire Size  $\leq$  .022" x .028" (.56mm x .71mm)

HF-678-104

### Multi-Use Cutter

Perma Sharp tips ensure long life. Wire Size  $\leq$  .015" (.38mm)

HF-678-105

### Mini Pin & Ligature Cutter

Slim profile is a plus in difficult to access areas. Perma Sharp tips ensure long life. Wire Size  $\leq$  .012" (.30mm)

HF-678-106

### Micro-Mini Pin & Ligature Cutter

Ultra slim profile is a plus in the most difficult to access areas. Perma Sharp tips ensure long life. Wire Size  $\leq$  .012" (.30mm)

HF-678-107

### Micro-Mini Pin & Ligature Cutter with Long Handle

Same as 678-107, with long handle. Wire Size  $\leq$  .012" (.30mm)

HF-678-107L

### Micro-Mini Pin & Ligature Cutter, 15°

15° offset tip with ultra slim profile is a plus in the most difficult to access areas. Perma Sharp tips ensure long life. Wire Size  $\leq$  .012" (.30mm)

HF-678-109

### Micro Cutter

Reduced tip size. Perma Sharp tips ensure long life. Wire Size  $\leq$  .012 (.30mm)

HF-678-110

### Flush Cut & Hold Distal End Cutter

Cuts and holds wire distally flush to the buccal tube while reducing debondings. Wire Size  $\leq$  .020" (.51mm)

HF-678-111

### Slim Flush Cut & Hold Cutter

Cuts and holds wires flush to the buccal tube. 50% thinner end than 678-111 Ability to cut up to .021" X 0.25". Wire Size up to .021 x .025. (.53mm x .64mm)

HF-678-113

### Slim Line Micro

Designed to cut soft wire pins and ligature. Wires Size up to .012" (.30mm).

HF-678-500

## Wire Cutters

Angle cutter for soft and hard wire.

AE-DP540



# Orthodontic Instruments

## Utility Pliers

### Clear Collection

Clear Collection consists of innovative instruments designed to accent, individualise and optimise the biomechanics of the "invisible" aligner experience.

#### The Tear Drop

An instrument created to add a notch at the gingival margin of clear aligners anywhere elastic hooks are required within the arch.

HF-678-800

#### The Vertical

Emphasises overcorrection of rotations by producing an indentation at the mesial or distal line angles on the facial and/or lingual of a specific tooth.

HF-678-801

#### The Horizontal

Produces an indentation to accent individual root torque.

HF-678-802

#### The Hole Punch

Creates half-moon cutouts for relief of clear aligners around bonded buttons used for elastics. No need to attempt to attach buttons to plastic trays.

HF-678-803

## Utility Pliers

### Ultra Slim Weingart

Versatile utility pliers with serrations, accommodating wires up to .030" (.76mm).

HF-678-501

### Slim Hammerhead NiTi

Multi-use pliers for bending NiTi up to .025" (.64 mm).

HF-678-502

**Hu-Friedy**

How the best perform



### Slim Band Remover

Slim band removing pliers with tapered beak providing improved visibility and access.

HF-678-503

### Weingart Pliers

Versatile utility pliers with serrations. Tapered and angled at 32° for placement and removal of archwires.

HF-678-201

### Slim Weingart Pliers

Versatile utility pliers with serrations.

HF-678-202

### How Pliers

How utility pliers have serrated tips. Useful for placement and removal of archwires, pins and other auxiliaries.

HF-678-203

### How Pliers, Offset

How utility pliers have serrated tips. Useful for placement and removal of archwires, pins and other auxiliaries.

HF-678-204

### Adhesive Removing Pliers

Reversible, replaceable blade is used to remove the bulk of composite material after debonding attachments.

HF-678-206

### Posterior Band Removing Pliers, Long

The alignment of the tip and occlusal pad is well suited.

HF-678-207

### Posterior Band Removing Pliers, Short

The offset alignment of the tip and occlusal pad is well suited for band removal.

HF-678-208



# Orthodontic Instruments

## Force Module Separating Pliers

Used for placement of elastic separators.

HF-678-210

## Coon Style Ligature Tying Pliers

Most popular pliers for tying stainless ligature wires. Wire Size  $\leq$  .015" (.38mm)

HF-678-211

## Bracket Removing Pliers

Works well for steel, ceramic and plastic brackets of all types.

HF-678-219

## Angulated Bracket Removing Pliers

High angle design for removal of all bracket styles from anterior and posterior teeth. With long handle.

HF-678-220L

## Crown & Band Contouring Pliers

Designed for fitting of stainless steel crowns.

HF-678-221

## Slim Crown & Band Contouring Pliers

Same as 678-221, with slimmer beak.

HF-678-221M

## Band Crimping Pliers

Thin tips allow for easy access into molar and bicuspid bands.

HF-678-225

## Utility Instrument

**Hu-Friedy**

How the best perform

## Wire Forming Distal Cinch Back Hand

Bends stainless steel archwire distal to the 1st or 2nd.



HF-678-506

## Wire Forming Pliers

**Hu-Friedy**

How the best perform

Wire forming is more than a skill, it is an art. These instruments perform as natural extensions of your professional artistry. Our impeccably smooth finished tips allow you to bend the wire without risk of damage. Whether you are shaping stainless steel or nickel titanium, Hu-Friedy has a hand finished, precision instrument to help bring your brilliant workmanship to life.

### Jarabak Pliers

Versatile light wire bending pliers. Grooves assist with precise bending and closing of loops. Wire Size  $\leq$  .020" (.51mm)

HF-678-301

### Three Jaw Pliers

Useful for contouring and bending wire. Wire Size  $\leq$  .030" (.76mm)

HF-678-302

### Omega Loop Forming Pliers

Graduated cone is designed to form precise omega loops: .045", .060" and .075". Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-303

### Bird Beak Pliers

Versatile loop forming pliers for wires. Wire Size  $\leq$  .030" (.76mm)

HF-678-304

### Light Wire Bird Beak With Groove

A longer cone than a standard bird beak is useful for forming precise loops in light wire. The groove allows for locating and duplicating loops. Wire Size  $\leq$  .020" (.51mm)

HF-678-305

### Tweed Loop Forming Pliers

Loop forming pliers with .045", .060" and .075" diameters for creating bull and helical loops in wires. Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-306

### Arch Bending Pliers

Designed for placing first, second and third order bends, .050" blade width. Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-307

# Orthodontic Instruments

## Tweed Pliers, Slim

Designed for placing first, second and third order bends, .025" blade width.  
Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-307S



## Rectangular Arch Bending Pliers

Designed for placing first, second and third order bends, .070" blade width.  
Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)


HF-678-308



## Step Pliers

1/2mm Surgery Set Up

HF-678-340



3/4mm Surgery Set Up

HF-678-341



1mm Surgery Set Up

HF-678-342



## Hollow Chop Contouring Pliers

Perfect for subtle contouring of arch wires.  
Wire Size  $\leq$  .030" (.76mm)

HF-678-310



## Mini Three Jaw Pliers

Tips are rounded for ease of contouring light wire. Wire Size  $\leq$  .030" (.76mm)

HF-678-312



## Light Wire Bird Beak

A longer cone than a standard bird beak is useful for forming precise loops in light wire. Wire Size  $\leq$  .020" (.51mm)

HF-678-318



## Adams Pliers

Excellent precise right angle bends and appliance application.  
Wire Size  $\leq$  .040" (1.01mm)

HF-678-320



## V-Stop Pliers

Places V bend in archwires, including nickel titanium, to prevent archwire from traveling.  
Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-321



## NiTi Three Jaw Pliers

Designed for forming and contouring of all archwire, especially for NiTi. Wire Size  $\leq$  .020" (.51mm)

HF-678-302



## Optical Pliers

Excellent for forming helical closing loops and arch adjustments. Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-323



## Bird Beak With Cutter

Smooth bending surfaces for a range of loop sizes while providing a sharp blade for cutting.  
Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-325



## Long Tapered Bird Beak Pliers

Bends a variety of smooth loops.  
Wire Size  $\leq$  .022" x .025" (.53mm x .64mm)

HF-678-326



## Hammerhead NiTi Tie-Back Pliers

Multi-use NiTi bender, no heat.  
Wire Size  $\leq$  .025" (.64mm)

HF-678-327




# Orthodontic Instruments

## Orthodontic Pliers

**Hu-Friedy**

How the best perform



Slimmer tip makes this ideal for gripping and placing elastic auxiliaries.

**Ortho Mathieu, Narrow Tip Pliers**

HF-678-330



With wider, stronger tips.

**Ortho Mathieu Wide Tip Pliers**

HF-678-331



Thin tips made with carbide inserts.

**Ortho Mathieu, Perma Sharp**

HF-678-332

## Orthodontic Tweezers

**Slim Direct Bond Bracket**

**Hu-Friedy**

How the best perform



Slim tipped bracket tweezers with tapered tail to aid in bracket manipulation.

HF-678-504

**Slim Buccal Tube Bonder**

**Hu-Friedy**

How the best perform



Ultra slim tipped buccal tube tweezers to aid in placement.

HF-678-505

## Orthodontic Scissors

**Wire Cutting Scissors Straight**

**Hu-Friedy**

How the best perform



Serrated blade, notch for soft wire-cutting.

HF-WCS

## Orthodontic Accessories

**Replacement Parts**

**Hu-Friedy**

How the best perform

HF-678-001 Replacement Blade for 678-206 Qty 1  
HF-678-014 1/4" Replacement Pads for 678-207, 678-208 Qty 6  
HF-678-036 3/16" Replacement Pad for 678-207, 678-208 Qty 6  
HF-678-014A 1/4" Replacement Pad for 678-206 Qty 6  
HF-678-011 Silicone Insert for 678-111 Qty 3  
HF-678-002 Replacement Blade for 678-303 Qty 1  
HF-678-113 Replacement Silicone Inserts  
HF-678-013 Pack of 3 each

# Instrument Cassettes & Accessories

## IMS INFINITY Series™ Cassettes



IMS Infinity Series™ Cassettes are made of long-lasting, low maintenance stainless steel. Soft, coloured silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and colour options, cassettes can be specialised to meet your needs.

### Exam Cassette - 5 Instruments

Exterior Dimensions: 3" x 8" x 1.25" (76mm x 203mm x 32mm). Holds 5 instruments.




### Small Cassette - 8 Instruments

Exterior Dimensions: 5.5" x 8" x 1.25" (140mm x 203mm x 32mm). Holds 8 instruments with an accessory area.




### Small Cassette - 10 Instruments

Exterior Dimensions: 5.5" x 8" x 1.25" (140mm x 203mm x 32mm). Holds 10 instruments.



### Signa-Stat Cassette - 12 Instruments

Exterior Dimensions: 6.5" x 10.5" x 1.25" (165mm x 267mm x 32mm). Holds 12 instruments with an accessory area.




### Large Cassette with Accessory Area - 12 or 16 Instruments

Exterior Dimensions: 8" x 11" x 1.25" (203mm x 280mm x 32mm). Also available in 8" x 10.5" x 1.25" size (IMN416xQ) to fit unique ultrasonic tanks.




### Large Cassette - 20 Instruments

Exterior Dimensions: 8" x 11" x 1.25" (203mm x 280mm x 32 mm). Holds 20 instruments.




### Double-Decker® Cassette - 8 Instruments

Exterior Dimensions: 3" x 8" x 1.5" (76mm x 203mm x 38mm). Holds 8 instruments.




### Double-Decker® Cassette - 14 Instruments

Exterior Dimensions: 4.5" x 8" x 1.5" (114mm x 203mm x 38mm). Holds 14 instruments.




### Double-Decker® Cassette - 20 Instruments

Exterior Dimensions: 6" x 9.7" x 1.5" (152mm x 246mm x 38mm). Holds 20 instruments with an accessory area. Has a built-in needle cap holder and slot for a winged syringe.



### Oral Surgery Cassette

Exterior Dimensions: 8" x 14.5" x 1.25" (203mm x 368mm x 32mm). Accommodates 12 instruments, 4 elevators, 3 hinged instruments in the lid and a separate accessory area.



### Ultrasonic Insert/Handpiece Cassette

Small: Exterior Dimensions: 1.5" x 8" x 1.25" (38mm x 203mm x 32mm)  
Large: Exterior Dimensions: 3" x 8" x 1.25" (76mm x 203mm x 32mm)



### Radiography Cassette

Exterior Dimensions: 5.5" x 8" x 1.25" (140 mm x 203 mm x 32 mm)



### Double-Decker® Small Orthodontic Cassette

Exterior Dimensions: 4.5" x 8" x 1.5" (114mm x 203mm x 38mm). Holds up to 4 pliers, a Mathieu needle holder and 2 hand instruments






# Instrument Cassettes & Accessories



**Large Thin Orthodontic Cassette**  
Exterior Dimensions: 8" x 11" x 0.75"  
(203mm x 280mm x 19mm). Holds up  
to 7 pliers and 4 hand instruments.



**Small Thin Orthodontic Cassette**  
Exterior Dimensions: 7" x 8" x 0.75"  
(178mm x 203mm x 19mm) Holds up  
to 4 pliers and 3 hand instruments.

## Instrument Cassettes

- Stainless steel locking instrument cassette tray with silicone holding rails.
- Fully vented design for optimal sterilisation.
- Upper and lower rails are interchangeable.
- Suitable for thermo-disinfector and autoclave.




- 5 Instrument Cassette Blue**  
HS-9008543 175mm x 69mm x 34mm
- 10 Instrument Cassette Blue**  
HS-9008544 175mm x 144mm x 34mm
- 20 Instrument Cassette Blue**  
HS-9008545 279mm x 182.5mm x 32mm



## Aluminium Instrument Trays

- Aluminium anodised
- Solid or perforated trays
- Size: 18cm x 28cm
- Not suitable for thermo disinfector



- |                            |                       |
|----------------------------|-----------------------|
| <b>Perforated Trays</b>    | <b>Solid Trays</b>    |
| HS-9001276 Silver          | HS-1005433 Silver     |
| HS-9001270 Green           | HS-9001274 Green      |
| HS-9001271 Red             | HS-9001273 Yellow     |
| HS-9001272 Blue            | HS-1008936 Blue       |
| <b>Perforated Tray Lid</b> | <b>Solid Tray Lid</b> |
| HS-9001433 Silver          | HS-1025760 Silver     |

## Accessories

### IMS Parts Box



Designed to keep small items such as rubber dam clamps and matrix band retainers organised during the sterilisation process.

- HF-IMS-1273 IMS Large Parts Box  
Exterior Dimensions: 3.5" x 1.12" x 1" (89mm x 29mm x 25mm)
- HF-IMS-1271 IMS Small Parts Box  
Exterior Dimensions: 1.56" x 1.25" x 1" (40mm x 32 mm x 25mm)

### IMS Cassette Barcode Clip



Snap-on barcode for outside of cassettes.  
Dimensions: .75" x 2.26" x .067" (19mm x 57.4mm x 1.7mm)

HF-IMS-CABL

### IMS Endodontic Stand



Hinged, double-lid stand holds endodontic files and reamers. Lid is etched with a ruler to measure stop positions easily from either direction.



- HF-IMS-1275 Holds 24  
Exterior Dimensions: 2.6" x 2.25" x 1.2" (67mm x 57mm x 30mm)
- HF-IMS-1275D Holds 48  
Exterior Dimensions: 2.6" x 2.25" x 1.2" (67mm x 57mm x 30mm)

# Instrument Cassettes & Accessories

## Air and Water Clips



Designed for IMS Cassettes. Air and water clips can be permanently installed in your cassettes.



- HF-IM1005N Infinity Series Air/Water Clip
- HF-IM10005 Signature Series Air/Water Clip
- HF-IM1006 Infinity and Signature Series Double-Decker Air/Water Clip

## IMS Bur Cushions



A hollow, rectangular, silicone base with a stainless steel lid that is designed to hold burs securely in place during cleaning and sterilisation. Hollow design allows drainage and proper cleaning.

### IMS Bur Cushion, Short Lid

HF-IMS-1372S Holds 12 burs.  
Exterior Dimensions:  
3.3" x 1.3" x 0.9"  
(84mm x 33mm x 22mm)

### IMS 1/2 Bur Cushion, Short Lid

HF-IMS-1372SH Holds 6 burs.  
Exterior Dimensions:  
1.75" x 1.3" x 0.9"  
(44mm x 33mm x 22mm)

### IMS Bur Cushion, Tall Lid

HF-IMS-1372T Holds 12 burs.  
Exterior Dimensions:  
3.3" x 2.44" x 0.9"  
(84mm x 62mm x 22mm)

### IMS 1/2 Bur Cushion, Tall Lid

HF-IMS-1372TH Holds 6 burs.  
Exterior Dimensions:  
1.75" x 2.44" x 0.9"  
(44mm x 62mm x 22mm)

## IMS Rubber Dam Clamp Organising Board



Dimensions: 4.6" x 2.4" x .375" (117mm x 61mm x 10mm). Holds up to 8 clamps during sterilisation. Fits into most cassette accessory areas.



HS-RDCOB

## Hinged Instrument Clips



Secure scissors, hemostats and other hinged instruments in lid of cassette.



- HF-IM1000 Hinged Instrument Clips, 2 Pack
- HF-IM1002 Hinged Instrument Clips, 6 Pack

## Needle Cap Holder



Maintains cap in upright position for one-handed needle recapping.



- HF-IM1001 Needle Cap Holders, 2 Pack
- HF-IM1003 Needle Cap Holders, 6 Pack
- HF-IM1004 Combination Pack\*\*

## IMS Cassette Autoclave Wrap



IMS Cassette Wrap is engineered to withstand the harsh conditions of sterilisation. Our unique compilation of textiles is designed to be strong and soft for easy wrapping and durability. It offers the best combination of softness and added strength. This low-lint, memory free wrap is available in two grades to meet your needs. Manufactured for single use only.



- HF-IMS-1210 24" x 24" (610 mm x 610 mm) For Large Cassettes, 500 Sheets
- HF-IMS-1212 20" x 20" (508 mm x 508 mm) For Large Cassettes, 500 Sheets
- HF-IMS-1211 15" x 15" (380 mm x 380 mm) For Small Cassettes, 1,000 Sheets
- HF-IMS-1213 12" x 12" (305 mm x 305 mm) For Small Cassettes, 1,000 Sheets

# Instrument Cassettes & Accessories

## IMS Monitor Tape



Strong 3/4" (19 mm) tape can be written on and visually differentiates processed from non-processed cassettes. Pre-printed tape is available with specific procedural set-ups, colour-coded or blank. 60 yards (55 m) per roll

HF-IMS-1243 Ortho	HF-IMS-1256 Exam
HF-IMS-1244 Archwire	HF-IMS-1257 Extraction
HF-IMS-1245 Banding	HF-IMS-1258 Perio
HF-IMS-1246 Bonding	HF-IMS-1259 Propphy
HF-IMS-1247 Restorative	HF-IMS-1260 Pedo
HF-IMS-1248 Surgery	HF-IMS-1262 Gray
HF-IMS-1249 X-Ray	HF-IMS-1263 Orange
HF-IMS-1250 Blank	HF-IMS-1264 Purple
HF-IMS-1251 Amalgam	HF-IMS-1265 Yellow
HF-IMS-1252 C&B	HF-IMS-1266 Red
HF-IMS-1253 Cement	HF-IMS-1267 Green
HF-IMS-1254 Composite	HF-IMS-1268 Blue
HF-IMS-1255 Endo	HF-IMS-1269 Black

## IMS Tape Dispensers



Rubber suction cups on the base secure dispenser to table or counter. Dispenser can also be mounted under cabinet.

HF-IMS-1270 IMS Tape Dispenser 10.75" x 6.25" x 7.6" (273 mm x 159 mm x 194 mm) Holds up to 9 rolls of IMS Monitor Tape.
HF-IMS-1270H IMS Mini Tape Dispenser 4.4" x 6.25" x 7.6" (111 mm x 159 mm x 194 mm) Holds up to 4 rolls of IMS Monitor Tape.

## Bagette Self-Sealing Sterilisation Pouches



Bagette® Self-Sealing Sterilisation Pouches are constructed with durable materials and thick seals. Built-in external and internal quality control indicators relay accurate time and temperature verification, allowing you peace of mind. The full line of pouches include those custom fit especially for our cassettes and smaller sizes for all other sterilisation needs.

HF-IMS-1346 2-1/4" x 4" 200 Rubber Dam Clamps or Burs
HF-IMS-1347 2-3/4" x 9" 200 1 Handpiece or 1-3 Hand Instruments
HF-IMS-1236 3-1/2" x 9" 200 2-5 Hand Instruments
HF-IMS-1348 5-1/4" x 10" 200 Handpiece or Multiple Instruments
HF-IMS-1237 7-1/2" x 13" 100 Exam Cassette†
HF-IMS-1238 10" x 14" 100 Small Cassette†
HF-IMS-1239 12" x 17" 100 Large Cassette†
HF-IMS-1345 13" x 20" 100 Oral Surgery Cassette†

† Please see Hu-Friedy.com for complete list of suitable cassettes for this pouch.

## Instrument Lubricants

### Instrument Lubricant

Prevents instrument discolouration and build up of mineral deposits and keeps hinged instruments moving freely.

HF-ILS Non-Aerosol  
Pump Bottle 8oz



### Shine Renew

Shine Renew instrument wipes effectively remove stains and rust, restoring the lustre to stainless steel, brass, copper and anodised aluminium. Keeps hinged instruments moving freely by dissolving built-up residue.

HF-IMS1455 Tub 20



# Sharpening Accessories

## Sharpening Stones & Accessories

### Sharpening Stones



A uniform coating of industrial diamonds bonded to high quality stainless steel provides an extremely hard, flat surface for instrument sharpening.

#### Extra Fine Grit Diamond Sharpening Card

Used to create very fine edges on the blades of instruments.

HF-DSCXFINE



#### Fine Grit Diamond Sharpening Card

For routine or regular maintenance sharpening.

HF-DSCFINE



#### Medium Grit Diamond Sharpening Card

For instruments badly worn or needing reconditioning.

HF-DSCMED



#### Extra Fine, Fine and Medium Grit Diamond Sharpening Cards

One card of each grit to handle a variety of instrument sharpening needs.

HF-DSCSET



#### Arkansas Wedge Shape

combines flat surfaces with rounded edges for sharpening straight and curved edges. 11.4 x 4.8 x 0.9cm

HF-SS6A



#### Arkansas Flat #4

Restores the cutting edge on dull instruments. 10 x 2.5 x 0.6cm

HF-SS4



#### Arkansas Bates #5

Used to sharpen curettes, bone curettes, elevators and explorers.

HF-SS5



#### Dual Grit Ceramic Stone, Blue Coarse Grit with White Fine Grit

Features a blue coarse side for recontouring and a white side for finishing and honing. 12 x 2.5 x 1.6cm

HF-SSDGRIT



### Sidekick Sharpener



For routine maintenance sharpening of scalers and curettes. Cordless, battery operated, compact and requires minimal space

HF-SDKIT Kit Contains: Sidekick Unit, User Manual & CD ROM, Ceramic Sharpening Stone, Plastic Test Stick, Magnifying Lens and 2 AA Batteries



### Plastic Test Sticks



Hard plastic test stick used to determine the sharpness of an instrument. Pack 6.

HF-PTS



### Sharpening Oil



Pre-sharpening application to stone, floats away metal particles and prevents damaging heat build up. 1oz bottle.

HF-SSO

