





# RELY ON US

every smile counts










INFECTION
PREVENTION
BEST PRACTICE

specialty catalogue




# The perfect decontamination room

It's best practice – starting from the dirty area and finishing in the clean area, to have a distinctive work flow with a set down, wash and rinse sink, washer/disinfector, inspection/packaging, sterilisation and cooling areas.



Reference: AS/NZS4815:2006 Standards New Zealand; Office based health care facilities — Reprocessing of reusable medical and surgical instruments and equipment, and maintenance of the associated environment. Reproduced for training and reference purposes only. Figure 1.1 Suggested layout for a reprocessing area.





<sup>\*\*</sup>Optional instrument washer/disinfector (\*or provision for future upgrade).

### Welcome to the Henry Schein Dental Infection Prevention Catalogue

### Page Guide



2 - 3 The Perfect Decontamination Room



**5** - **10** Working Environment



11 - 12 Waterline Maintenance



13 Hand Hygiene



14 - 18 Personal Protection



19 - 20 Surface Management



21 - 35 Instrument Reprocessing



36 - 39 Sterilisation



40 Storage



### **Barrier Sleeves and Wraps**

### **Disposable Headrest Covers**





#### **Features and Benefits**

- Disposable poly barrier protects patient and dental chair headrest
- These economical covers come in a convenient space-saver dispenser box of 250

HS112- 4871 Regular 10"x 11.5", Clear, Box of 250 HS112- 4872 Jumbo 10" x 14.5" Clear, Box of 250

## Air Water (Triplex) Syringe Covers





### **Features and Benefits**

- Syringe covers with opening
- For 3-way air/water syringes

HS100-9860 Clear, 7cm x 26cm, Box of 500

### **Disposable Tray Inserts**





#### **Features and Benefits**

- Premolded disposable tray inserts that offer an organised, efficient and effective method for eliminating cross contamination
- Made of man-made fibres that are completely degradable through the absorption of light
- Standard size

HS900-0937 111/4" x 71/4" x 9/16. White. Pkt of 50

### **Tray Barriers**





### **Features and Benefits**

- · Made of high quality clear plastic
- Covers the entire instrument tray
- Minimises time spent cleaning trays
- · Overlapping folds keep the instrument trays from falling out

HS101- 4275 Clear, 10.5 x 14", Box of 500

# T-Bar MENRY SCHEIN® CORPORATE BRAN Light Handle Cover



HS100-9864 Clear, Box of 500

### X-Ray Covers





HS101-6507 Clear, 38 x 67cm, Roll of 250

### Barrier Cover-All Film



### **Features and Benefits**

 Cover-all Film disposable sheets with adhesive back

HS-100-9869 Clear, 6"x 4" 1200 per roll HS-1011838 Blue, 6"x 4" 1200 per roll



### **Cover-All<sup>™</sup> Dispenser**

### **Features and Benefits**

- Table top dispenser makes removal of films quick and effortless
- Dispenses adhesive films up to 4" wide
- Dispenses Cover-All infection control film and Tube Sleeve barrier film

KE-P1600 W 6" x H 7" x D 9"

### Slip-Nots™ Bib Chain

### **Features and Benefits**

- Colourful, autoclavable, bib clips made from a special soft flexible plastic
- Extra wide neck band prevents patient's bib from sliding out of place

KE-P2000A Assorted Colours Pkt of 4

### Disposable Head Rest Cover K△VO ■ Kerr

### **Features and Benefits**

• Disposable platic cover for headrests

KE-P3511 11" x 9.5" Pkt of 250 KE-P3514 14" x 9.5" Pkt of 250



K4VO Kerr

### K4VO ■ Kerr : Cover-All™

### **Features and Benefits**

- This thin sheet of adhesive plastic effectively covers anything you might touch during a procedure
- It is disposable, easy to apply and remove, and will not leave a residue



K4VO Kerr

KAVO Kerr

KE-P3700C Clear Roll 4"x 6" Pkt of 1200 Perforated Sheets KE-P3710B Blue Roll 4"x 6" Pkt of 1200 Perforated Sheets

\*Dispenser not included.

### **HP Sleeve**®

### Features and Benefits

- Disposable cover for handpieces, ultrasonic scalers, jet polishers
- Unique super soft plastic allows users to hold handpieces easily and maintain tactile feel
- Unit dose package makes it simple to insert and apply to handpiece



KE-P4050 2" x 10" Fits Contra Angle Slow Speed Handpieces & Midwest Shorty Pkt of 500

### **Disposable Tube Sleeve**

### **Features and Benefits**

 Disposable barrier for dental unit tubing

KE-P3720 2" 1200 feet per roll KE-P3740 4" 1200 feet per roll



K4VO Kerr

### **Disposable Light Sleeve**

### **Features and Benefits**

- Allows the user to quickly slip the plastic sleeve over the handle
- Patented design keeps the sleeve from falling off while in use

KE-P3600 'T' Style Pkt of 500 KE-P3650 'L' Style Pkt of 500







### **Disposable Tray Sleeve**

#### **Features and Benefits**

- Instruments are placed directly on top of the disposable plastic sleeve, which is made to replace the present paper tray cover used on instrument trays
- Unique fold lock end keeps the tray safely locked inside the sleeve
- Minimises staff time spent on cleaning the tray

KE-P3000B Small 10.5"x 14" Fits Ritter B, A-dec Pkt of 500 KE-P3100-A1 11 5/8" x 14.5" Fits Chayes, Divided Ritter B Pkt of 500 KE-P3200E Large 11.5" x 16" Fits Midwest, Weber Pkt of 500 KE-P3300F 7.5" x 10.5" Fits Mini Tray Pkt of 500 KE-P3500M 14" x 19.5" Fits Mayo Tray Pkt of 500

### **Disposable Chair Sleeve**

### K4VO∎Kerr

### **Features and Benefits**

- Tapered design conforms neatly to your chair
- It is disposable, no chair clean up is required
- Eliminates the need for a headrest cover and with the reduced usage of disinfectant, it extends the life of your chair

KE-P3800 27.5" x 24" Pkt of 225 KE-P3825 32" x 32" Pkt of 200

### Disposable Curing Light Tip Cover K♦VO■\*Kerr

### **Features and Benefits**

- Disposable cover for curing light tips
- Decreases possibility of cross contamination
- Prevents composite from bonding to end of tip, which decreases curing depth





### **Disposable Sensor Cover**

### **Features and Benefits**

 Sensor Covers are disposable covers for the sensor on digital X-ray imaging systems

KE-P3735 3.5cm x 20.5cm Pkt of 500 KE-P3736 4cm x 21cm Pkt of 500



### **Disposable Syringe Sleeve**

#### **Features and Benefits**

- Standard disposable sleeve made to cover the 3-way syringe, HVE and saliva ejector valves
- Sleeve with adhesive tape strip fastens bag to tubing to prevent sliding
- Sleeve with exclusive pre-cut opening allows the user to simply slide the syringe tip through the end of the bag



K4VO Kerr

KE-P3760 Standard Sleeve 2.5" x 10" Pkt of 500 KE-P3763WA With Adhesive Strip 2.5"x10" Pkt of 500 KE-P3767WO With Opening 2.5"x 10" Pkt of 500

### Disposable Impression Gun Cover

### K4VO Kerr

### **Features and Benefits**

 Disposable plastic sleeve made to cover the impression mixing gun used in cartridge type mixing systems

 Offers improved infection prevention, while keeping the gun clean of impression material



KE-P4600 Pkt of 400



### **Disposable X-Ray Sleeve**

### **Features and Benefits**

 The X-Ray sleeve's tapered design protects your x-ray from the harmful drying effects of surface disinfectant

KE-P3900 15" x 26" Pkt of 250 KE-P3950XL 23" x 31" Pkt of 250



### Safe-Tips® EZ

### **Features and Benefits**

- Directly interchangeable with A-dec<sup>®</sup> style metal tips
- Identical performance
- Constructed with plastic and stainless steel
- No syringe conversion required before use
- White colour

KE-P44150 Dispenser Pouch Pkt of 150

KE-P44600 Bulk Kit Pkt of 1600



K4VO Kerr

### Dispos-a-Trap®

### **Features and Benefits**

- The Dispos-a-Trap filter is a totally disposable filter designed to replace the present trap filter in your evacuation unit
- Reduces the possibility of infection and contamination, it improves the performance of your suction system and reduces infectious aerosol
- It does all this while eliminating one of the most hazardous and unpleasant tasks in the dental office, cleaning the waste that is present in your filter



K4VO\_Kerr

KE-P5500 Fits A-dec, Cascade and Performer Pkt of 144

KE-P5501 Fits A-Dec Models with Black Plastic Canisters, Plus Newer Models by Ampco, Belmont, DCI, Forest, Marus, MDT, Schein Pkt of 144

### **Cure Sleeve®**

### **Features and Benefits**

- Disposable cover for curing light handles
- Protects from cross contamination and the harmful effects of surface disinfectant
- Patented design prevents overheating of the curing light gun

KE-P4550 Fits 5" x 10" Handle Pkt of 400

### Seal-Tight® Adapter

### K4VO Kerr



### **Features and Benefits**

- All DentalEZ syringes manufactured after January 1997 will use an A-dec adapter
- After the adapter is installed, there is no need to loosen and tighten the thumbnut for each tip placement, simply push on the yellow auto lock ring and remove the tip and then push a new tip into the adapter and it's ready to use
- The inner yellow seal acts like a new 0-ring with each tip replacement, insuring dry air on demand
- The tip can be bent to a 90° angle for certain procedures

KE-P77021 Adaptor A-dec Styles KE-P77031 Adaptor Snap-On KaVo KE-P77200 Disposable Syringe Tips Pkt of 200



### **Disposables**

### Disposable Syringe Tip (Suretip)

### **Features and Benefits**

- Durable all plastic construction
- Allows visible water flow, no leaking of water and dry air
- Fits most models of air/water syringes that have not already been converted (A1 and P1 type syringes)

HS102-8981 Bag of 250



✓ Henry Schein®

### **Saliva Ejectors**

### **Features and Benefits**

- Non-removable tip with smooth, rounded edges and extra-soft tips
- Perfectly positioned suction ports to maximise fluid removal and protect delicate mucosa
- Latex-free with flexible bonded tips
- Will not aspirate tissue
- Standard size

HS100-4092 Clear with Blue Tip, 15cm, Bag of 100 HS100-5205 White, 15cm, Bag of 100 HS900-4272 Pink, 15cm, Bag of 100

### **Gauze Sponge**

### **Features and Benefits**

- Soft, pliable and absorbent
- Non-sterile
- 8-ply

HS100-6458 Unfilled, 3"x 3", Box of 200



### **Henry Schein Sterile Drape**

### **Features and Benefits**

- Fenestrated
- Sterile

HS900-4685 Box of 50



### 3 Way Syringe

### **Features and Benefits**

- Metal
- Triplex
- Autoclavable

HS104-8066 Pkt of 6



HENRY SCHEIN

### **Evacuator Tip**

#### **Features and Benefits**

- Autoclavable, plastic, high-volume evacuator tips
- Thick-walled design provides excellent cheek retraction

HS101-4856 Vented, Pkt of 50 HS100-4568 Non-Vented, Pkt of 50



HENRY SCHEIN

### **Disposable Air/Water Flash Tips**

### **Features and Benefits**

### **Fast**

- Quick and easy changeover after each patient
- No converters required (on most common syringes)

### Safe

- Disposable... to help prevent cross contamination
- Locking indicator and wide locking groove ensure tip is in place
- Smooth rounded edges to prevent harm to patients and equipment

### **Solid construction**

- Inner tube does not shift under pressure (unlike leading competitors)
- Separate air and water chambers to minimise crossover
- Built for strong cheek retraction

#### Sleek

• 5 great colours in each pack!

SU- 23071 Pkt of 250 SU- 23075 Pkt of 1200





Comfort Plus® Premium



Saliva Ejectors



The kinder, gentler saliva ejector... takes in fluid, not tissue!

### **Features and Benefits:**

- Patented tip reduces patient discomfort
- Optimises fluid removal
- Minimises trauma to mucosal tissue
- Facilitates uninterrupted fluid and light debris removal
- Prevents tissue from blocking suction

CX-ZCWCP Clear White 100 Pkt CX-ZWGCPM White/Green Minted Scented 100 Pkt CX-ZWPCPBG White/Pink Bubble Gum 100 Pkt CX-ZWWCP White/White 100 Pkt

### **Premium Plus® Saliva Ejectors**



### **Features and Benefits:**

- Wire-reinforced tube forms and maintains any desired shape
- Provides optimal suction without aspirating tissue
- Soft, non-removable tip bonded to the tube for maximum patient safety.

CX-ZCB Clear / Blue 100 Pkt CX-ZCMC Multi-Colored Tips 100 Pkt CX-ZCW Clear / White 100 Pkt CX-ZWW White / White 100 Pkt CX-ZPP Pink with a Purpose 100 Pkt



### **Advantage Saliva Ejectors**

### **Features and Benefits:**

- Wire-reinforced tube forms and maintains any desired shape
- Provides optimal suction without aspirating tissue
- Soft, non-removable tip bonded to the tube for maximum patient safety.

CX-ZBBI Saliva Ejectors Blue 100 Bag CX-ZCBI Saliva Ejectors Clear / Blue 100 Pkt CX-ZCCI Saliva Ejectors Clear / Clear 100 Pkt CX-ZWWI Saliva Ejectors White / White 100 Pkt



### MaxVac Tips

### Features and Benefits:

- Available in 5 colours blue, green, pink, purple, yellow
- Two tips in one vented end on one side and nonvented end on the other
- Unique "S" tip design with NO sharp edges for maximum patient comfort

CX-ZETMVBL Blue CX-ZFTMVGR Green

**CX-ZETMVP** Assorted Colours

CX-7FTMVPK Pink CX-ZETMVLV Lavender CX-ZETMVYE Yellow

### **Crosstex HVE Tip**

### **Features and Benefits:**

- Two tips in one vented end on one side and nonvented end on the other.
- Unique "S" tip design with NO sharp edges for maximum patient comfort.

**CX-ZET** White CX-7FS50 Slotted White CX-ZETNV Non Vented White S Curve

### **Evacuation**

### **Purevac® Evacuation System Cleaner**



CROSSTEX

### **Features and Benefits:**

- Trusted Brand the Number One Evacuation System Cleaner in the US
- Powerful cleaning dissolves and destroys a wide range of dental debris
- Avoid expensive downtime
- Removes calcium and deposits that can hinder your system's performance
- Fresh lemon fragrance leaves your system smelling great!
- Non-foaming and biodegradable
- Economical concentrated formula 1oz (30ml) makes 1 quart (33.8oz) of solution
- Dispenses easily with 1 oz (30ml) pump (included)

SU- 21113 Liquid 2 Ltrs SU- 21115 Liquid 5 Ltrs





### Aspirmatic® Cleaner, **Disinfectant and Deodouriser** for Suction and Spittoon

### **Features and Benefits:**

- Aldehyde free
- Simple and safe handling
- Proven by extensive practical tests with leading dental units
- For all kinds of amalgam traps
- No disruptive foam
- Economical use

SK-150106 Liquid 2 Ltrs

### schülke -> Orotol® Plus Disinfection



### **Concentrate Solution Features and Benefits:**

- Super concentrate for the simultaneous disinfection, deodorisation, cleaning and care of all dental suction systems and amalgam separators
- Aldehyde free and foam free
- Range of action: bactericide, tuberculocide, fungicide, anti-viral (covers viruses, adeno-viruses HBV, HCV, HIV)

DU-CDS110P61 Liquid 2.5 Ltrs





### Majac CliniJet®-WS (Weekly)

### **Features and Benefits:**

- CliniJet-WS should be used once or twice per week as part of the routine maintenance and sanitation of aspiration lines and to assit the removal of bioburden
- Highly effective sanitiser and descaler helps remove soils that are resistant to other detergents
- Incorporates unique de-foaming surfactants that ensure penetration of biofilm and contaminating soils
- Mildly acidic descaling agents based on natural Lactic acid and Citric acid are incorporated for their excellent compatibility with metal, plastic and rubber components

MJ-CJWS1 Liquid 1 Ltr





### **Orotol Ultra Powder Concetrate**



### **Features and Benefits:**

- Powdery concentrate for the automatic and simultaneous disinfection, deodorisation, cleaning and care of all suction systems, amalgam separators and spittoon bowls.
- Broad effective spectrum: bactericidal, tuberculocidal, fungicidal, virucidal inclusive of HBV, HCV and HIV.
- 1% application concentration only;
- Kind to the environment biodegradable according to OECD guidelines



### Majac CliniJet®-DS (Daily)

### **Features and Benefits:**

- CliniJet-DS should be used each day as part of the routine maintenance and sanitation of aspiration lines and to assit the removal of bioburden
- Passes the TGA Disinfectant Test Option B
- Highly effective sanitiser incorporating an advanced de-foaming formulation for rapid sanitation of the aspiration lines
- Incorporates unique de-foaming surfactants to ensure detergency of contaminating soils and penetration of biofilm
- Mild alkaline auxiliary agents enhance the activity of the sanitising agents and aid in the solubilisation of proteins and soils

MJ- CJDS5 Liquid 5 Ltrs

### **PULIJET Gentle Suction Cleaner**

Daily disinfection of your suction system is vital to keep your suction lines clean and maintain your equipment and suction unit. Puli-jet GENTLE is suitable for the following dental suction systems:

### **Features and Benefits:**

- Semi-wet / semi-dry systems
- Dry suction systems
- Wet suction systems
- 1L bottle makes 30 x 2.5L solutions

Puli-jet GENTLE has a simple and accurate dispensing cup, making it easier to measure the solution and to pour in the Pulse Cleaner.

CA-M060970 Liquid 1 Ltr









# 365 days of of compliant water with 1 simple cartridge installation

### Water Treatment Technology Developed for Space Now Available for Your Dental Practice

Independent Water Bottle Cartridge



Municipal Water Cartridge



### simple

- Installation in minutes
- Extremely cost-effective
- Compatible with all dental units and related equipment

### effective

- Fast-acting
- One install for 365 days of pure water
- · Broad-spectrum anti-microbial action

### reliable

- Eliminates human error
- Works for 365 days
- Testing is easy with available iodine test strips
- · Will not interfere with dental materials and bonding
- We'll remind you when it's time to change cartridges

### safe

- Non-toxic, no allergenic iodine
- No hazardous silver
- · No harsh chemicals
- Dispose of cartridges in regular trash



### **Dentapure Waterline Cleaner**

CR1-DP365B DentaPure Cartridge 365 Day Bottle CR1-DP365M DentaPure Cartridge 365 Day Municipal





Discover the problem you didn't know you had!

DON'T MISS IT!

DentaPure.com

Learn the water purification secret NASA has known for years and revolutionise the way your office handles water treatment.







### **Hand Washing & Moisturising**

### Sensiva® Hand Wash Lotion

### schülke -}

### **Features and Benefits:**

- Syndet based on non-irritant tensides
- Gentle on the skin
- Effective cleaning and free of colourings, perfume and wool wax
- Allantoin to protect the skin against irritations
- pH skin neutral

SK- 114006 500ml SK-114002 1Ltr

SK- 180303 3ml Hand Pump for 500ml and 1 Ltr



### Sensiva® Protective Emulsion

### schülke -}

### **Features and Benefits:**

- Moisturises and protects the skin while absorbing quickly
- Contains allantoin and bisabolol, thus heals wounds and is anti-inflammatory
- Does not reduce the disinfection effect, e.g. with hand disinfection
- Free of colourings and perfume, easy-spreading




SK- 116003 500ml

### Avagard™ General Hand & Body Wash



### **Features and Benefits:**

- For use where a fully compatible, general hand and body wash is required
- Cleans visibly soiled hands as a routine handwash
- Contains high levels of moisturisers and emollients
- Soap-free and gentle on the skin
- Inter-compatible with full range of Avagard handwashing and skin hygiene products
- Avagard D antiseptic



MD-9230P 500ml

### **Hand Disinfecting**

### **Desderman® Pure Hand** Disinfectant Gel Rub

### Features and Benefits:

- Bactericidal (incl. TB), fungicidal, virucidal against lipophilic and hydrophilic viruses: polio, rotaviruses, adenoviruses, vaccinia and herpes simplex viruses, HBV, HCV, HIV, noroviruses
- Norovirus efficacy tested (15 sec.) in accordance with latest test methods (murine norovirus) with bonus skin protection and care through re-fatting agents
- Hypoallergenic and quick drying due to ethanol base

SK- 126902 500ml SK-126903 1 Ltr

### Sensiva® Touchless **Dispenser Refills**

schülke -1-

schülke -1-

SK-114012 Sensiva Wash Lotion Touchless Refill 1L

SK-116812 Desderman Touchless Pure Refill 1 Ltr

SK-126911 Desderman Touchless Pure Gel Refill 1L

SK-668300 Touchless Dispenser For 1L Disinfectant and Wash Lotion Refills



### schülke -1-

### Microshield **Angel Blue**



#### **Features and Benefits:**

- Hygienic Hand Rub
- Kills bacteria without water
- pH balanced
- Enriched with emollient and moisturiser

SK-70000379 500ml Bottle

### schülke -1-

### Microshield **Angel Clear**



#### **Features and Benefits:**

- Colour-free
- Perfume-free
- pH balanced
- Kills bacteria without water
- Enriched with emollient and moisturiser

SK-70000723 500ml bottle



### **Safety Glasses**

### **Rush Plus**

### **Features and Benefits:**

- Ultra light weight
- Built in brow protection
- Adjustable PVC nose bridge plus co-moulded temple arms for superior comfort and grip
- Highlight Feature: New
   Platinum Coating A high
   performance anti-fog + anti-scratch
   coating system applied on both sides
   of the Clear / Smoke lenses

BL-1662301 Rush Plus Clear BL-1662302 Rush Plus Smoke

# ens ip

bollé

### **Bandido**

### **Features and Benefits:**

- Brilliant coverage for a wide variety of face shapes
- Single piece optically correct (10 base curve) distortion free lense

BL-1667201 Bandido Clear BL-1667202 Bandido Smoke





bollé

### **B Line BL11**

### **Features and Benefits:**

- Single lens wrap around
- Polycarbonate over spetacles
- Moulded brow gaurd
- · Vented side shields
- Anti-scratch and anti-fog coated lenses



BL-1681101 Clear Lens BL-1681102 Smoke Lens

### **Rush Plus Small**

### **Features and Benefits:**

- Ultra lightweight Built in brow protection
- Adjustable PVC nose bridge plus co-moulded temple arms for superior comfort and grip
- Smaller frame
- Highlight Feature: New
   Platinum Coating A high
   performance anti-fog + anti-scratch
   coating system applied on both sides
   of the Clear / Smoke lenses

BL-1672301 Rush Plus Clear Small BL-1672302 Rush Plus Smoke Small

### **B Line BL10**

### **Features and Benefits:**

- Lightweight spec
- Wrap around
- Anti-scratch and anti-fog coated lenses


BL-1653901 Clear Lens BL-1653902 Smoke Lens



### **B Line BL13**

### Features and Benefits:

- Comfortable spec
- Wrap around
- Tiltable
- Adjustable pins
- Anti-scratch and anti-fog coated lenses





bollé

### **Aprons and Gowns**

### **Disposable Gowns**

### **Features and Benefits:**

- Full-length gown made of fluidresistant, non-woven polypropylene material
- Full cut, open neck, open back, and waist ties
- Sewn neck binder and knit cuffs

HS100-5870 Blue Med/Large Pkt of 10 HS100-6784 Pink Medium Pkt of 10 HS100-9392 Pink Large Pkt of 10



### **Med-X Disposable Apron**

#### **Features and Benefits:**

- Applications: Heavy duty, wash-up, packing, healthcare
- · Low cost, limited life
- Size: 66cm x 110cm





FB-APRON White Apron Box of 100

### Masks and Face Shields

### **Earloop Mask - Level 1**

### **Features and Benefits:**





- 3-ply glass-free filter media surrounded by soft, non-irritating inner and outer facing
- Full-width flexible nosepiece provides fit and security
- Manufactured to exceed all current ASTMF 2100-03a standards

HS-1043809 Blue Box of 50 HS-1046611 White Box of 50 HS-1043730 Pink Box of 50 HS-104-8600 Lavender Box of 50

### Tie-On Masks - Level 1

### **Features and Benefits:**

- · Light, comfortable and latex-free
- A glass-free filter media surrounded by soft, non irritating inner and outer facing
- Full-width flexible nosepiece provides fit and security
- Manufactured to exceed all current ASTMF 2100-03a standards
- Superior fluid protection, passed standard Fluid Penetration Resistance test

HS-1046809 Blue Triple Layer Box of 50

### **Face Shield**

#### Features and Benefits:

- Latex-free, optically clear, full-length visor provides increased protection against splash, spatter and aerosols
- Features a 11/2" thick, 9" long foam headband that keeps the face shield off the face
- Ultrasonically welded, 13" elastic headband for added strength and increased wearer comfort
- Special fog-free coating on inside and outside of visor

HS101-2254 W 13" x L 7 1/2" Box of 24

### Comfit Earloop Mask - Level 1

### **Features and Benefits:**

- Two-way DELNET® film barrier provides fluid protection with excellent breathability
- 99% at 0.1 micron, fluid resistant at 80mm/Hg
- Foam strip keeps glasses fog free
- Delta P of 2.7
- Hypoallergenic and latex-free
- Available in earloop or tie-on design with anti-fog face shield

HS100-835 Blue Super High Filtration Box of 40







HENRY SCHEIN

### **Maxi-Gard Earloop Mask**

- Level 2

#### **Features and Benefits:**

- Level 2 barrier resistance 120mmHa
- Anti-fog and earloop
- Available in pink, blue and white

HS-HRSUMK2279 Pink Box 50





- Level 2

**Features and Benefits:** Level 2 fluid resistance

Health Fluidshield Tie-on Mask

- Surgical mask with ties
- Anti-fog



KC-62113 Surgical Mask with Ties

HS-HRSUMK2282 Blue Box 50 HS-HRSUMK2283 White Box 50

### Procedural Earloop Mask - Level 2





### **Features and Benefits:**

- Level 2 fluid resistant outer layer
- Patent pending design with aluminum strip under chin for custom fit
- Fluid resistant white inner layer
- Extra-long aluminum nose piece
- PFE = 99.7% @ 0.1 micron

CX-GCPBLSF Earloop SecureFit Blue CX-GCPLVSF Earloop SecureFit Lavender CX-GCS Surgical Tie-On Laces

### Isofluid Earloop Mask - Level 1



M HALYARD



### **Features and Benefits:**

- Level 1 fluid resistant outer layer
- Patent pending design with aluminum strip under chin for custom fit
- White tissue inner layer
- Extra-long aluminum nose
- PFE = 98.9% @ 0.1 micron

CX-GCIBLSF Earloop SecureFit Blue 50 Box CX-GCILVSF Earloop SecureFit Lavendar 50 Box CX-GPLUSWHSF PlusPure SecureFit White 50 Box CX-GCICXBSF FogFree Earloop with SecureFit Blue 50 Box

### **Ultra Sensitive Earloop Mask - Level 3**

### **Features and Benefits:**

- Level 3 fluid resistant outer laver
- PFE = 99.8% @ 0.1 micron
- Patent pending design with aluminum strip under chin for custom fit
- White hypoallergenic inner cellulose layer
- Extra soft will not lint, tear or shred
- Free of chemicals, inks and dyes
- Extra-long aluminum nose piece

CX-GCFCXSSF Earloop 50 Pack White CX-GCFCXSFSF FogFree Earloop 40 Pack White CX-GCPWSSF FogFree Earloop W/Shield W/SecureFit





### **Latex Gloves**

### **Criterion Latex Powder-Free Gloves**

#### **Features and Benefits:**

- Higher elasticity and higher tensile strength than most comparable gloves on the market
- New improved powder-free manufacturing process ensures minimum chlorine smell
- Upper limit of 50mcg/dm2 of extractable latex protein which minimises the chances of users developing a latex allergy reaction
- The textured fingers ensure better grip when handling instruments

HS-102-5354 Extra Small Box 100 HS-102-8296 Small Box 100 HS-102-5355 Medium Box 100

HS-102-7763 Large Box 100 HS-102-5357 X-Large Box 100

### **Aloe Latex Textured Powder-Free Gloves**

#### **Features and Benefits:**

- Contain Aloe Vera Gel to moisturise and hydrate the skin
- Manufactured to EN455123 and ISO standards
- Green with excellent tactile sensitivity and textured for better grip
- Powder-Free latex
- Aloe content and purity in this product is certified by the international Aloe Science Council

HS-9002871XS Extra Small Box 100 HS-9002872S Small Box 100 HS-9002873M Medium Box 100

HS-9002874L Large Box 100 HS-9002875XL X-Large Box 100

11111 (36)

### ✓ HENRY SCHEIN®



### Latex Lightly Powdered Gloves Henry Schein



### **Features and Benefits:**

- Economy lightly powdered latex gloves with beaded cuff
- Non-sterile
- 9.5" rolled beaded cuff
- Smooth exterior finish
- Thiuram-free



HS-HSC8010XS Extra Small Box 100 HS-HSC8010S Small Box 100

HS-HSC8010M Medium Box 100 HS-HSC8010L Large Box 100

### DentaGlove Latex Powder-Free Gloves Ansell

#### **Features and Benefits:**

- Powder free, natural rubber latex gloves offering softness, comfort and excellent grip
- Ambidextrous
- Low allergenicity glove due to both low chemical residues and latex protein (Total Protein Levels <50 ua/a)

MICROTOUCH

- Thiurams and mercaptobenzothiazoles free reduces the risk of contact dermatitis
- Layer pack for easier dispensing and reduced wastage
- Full size range: XS-XL (5 sizes)

AN-465 Sizes: XS - XL Box 100

### Synthetic Gloves

### Nitrile Powder-Free Gloves



### **Features and Benefits:**

- Excellent fit
- Powder-free, nitrile gloves
- Super elastic for excellent comfort and fit
- Textured fingers for improved grip
- · Pink with Bubblegum scent



HS-9003964 Extra Small Box 100 HS-9003964 Small Box 100

HS-9003965 Medium Box 100 HS-9003966 Large Box 100

### **Nitrile Powder-Free Gloves**



### Features and Benefits:

- Latex free
- Powder free
- Soft nitrile formulation helps give a superior range of motion around thumb area



- Textured palm and fingers for excellent tactile sensitivity
- Odourless
- Box of 200

HS-90086\_\_ XS-XL White Powder Free Box 200

### **DentaGlove Blue Nitrile**

#### **Features and Benefits:**

- ERGOFORM ergonomic design technology reduces hand fatigue during long procedures
- Innovative film technology provides comfort and grip
- Robust strength helps avoid rips and tears while donning

AN-31350XS X-Small 100 Box AN-31350S Small 100 Box AN-31350M Medium 100 Box AN-31350L Large 100 Box AN-31350XL X-Large 90 Box



### **DentaGlove Neoprene**

### **Features and Benefits:**

- Unique formulation for extra sensitivity with flexibility and comfort
- Strong and durable to avoid nagging and tearing through long procedures and donning
- ERGOFORM™ ergonomic design technology reduces hand fatigue during long procedures
- · Low to no odour

AN-31360XS Extra Small 100 Box AN-313600S Small 100 Box AN-313600M Medium 100 Box AN-313600L Large 100 Box AN-31360XL X-Large 90 Box






### Ansell

### **Sterling Nitrile Exam Gloves**

### M HALYARD

### **Features and Benefits:**

- Halvard's manufacturing process boosts the force required to break through the glove and increases its already superior barrier protection qualities
- The process also allows a reduction in glove thickness, which means that less force




HALYARD

KC-139\_ Silver XS - XL Box of 200

### **Criterion Pure Nitrile Powder-Free Gloves**

#### **Features and Benefits**

- Accelerator-free, powder-free, latex-free
- Prevent type IV allergic reactions as the gloves do not contain thiurams, carbamates or thiazole



- Conforms with FDA's 1000ml watertight test and En455 (1 & 2)
- Meets ASTMD 6319-00a standards
- Beaded cuff and extended palm cut reduces hand and thumb fatigue

HS-9004483 Extra Small Box 100 HS-9004484 Small Box 100 HS-9004485 Medium Box 100

HS-9004486 Large Box 100 HS-9004487 X-Large Box 100

CRITERION PURE NITRILL

### **Lilac Utility Gloves**



### **Features and Benefits:**

- Flocklined nitrile utility gloves for use when cleaning and processing instruments
- Puncture and chemical resistant with a non-slip grip
- Steam autoclavable to 250°F (121°C) up to five times

HF-4006\_ S, M and L Pkt of 3

### Gammex® Latex



### **Features and Benefits:**

- Gammex® PF provides excellent donning and comfort through the use of a hydrophilic/hydrophobic polymer lining allowing for both wet and dry donnina
- The wide thumb ball area combined with Ansell's soft latex formulation ensures extended comfort and less hand fatigue
- Gammex® PF has a textured finger and palm yet a smooth back providing excellent grip and easy double donning
- Gammex® PF, you can feel the difference



### Alcohol-Free Cleaning & Disinfection Solutions & Wipes

# Mikrozid® Alcohol-Free schülke --Cleaning & Disinfecting Wipes & Liquid



### **Features and Benefits:**

- Broad spectrum of efficacy including bacteria, fungi and viruses (inc HBV, HCV and HIV)
- Ideal for all surfaces especially dental chairs and alcohol sensitive surfaces
- Alcohol free, surfactant based product with a 1 minute contact time
- Ready to use solution
- Dries quickly, ideal for use between patients

SK- 109184 Wipes Tub 200 SK- 109186 Liquid 1 Ltr SK- 109190 Wipes Refill 200 SK- 109193 Liquid 5 Ltrs

# Neutral Detergent Wipes



Cleans soil and contaminant from surfaces prior to disinfection.



HS-NWIP\_\_ 220 Wipes Tub/Refill

### Mikrozid Sensitive Premium Wipes Features and Benefits:

- Alcohol-free
- Ready-to-use, disinfectant-impregnated wipes
- Very fast-acting and effective
- For alcohol-sensitive surfaces
- Good cleaning performance

### schülke -}



K-70806 50 Pack

# Optim<sup>®</sup> Blue Cleaning & Disinfecting Wipes & Liquid

#### **Features and Benefits:**

- Optim Blue disinfects using Accelerated Hydrogen Peroxide® (AHP®) that has virtually no odour
- Kills broad spectrum of organisms with short contact times (including Polio virus in 1 min)
- Does not contain VOCs (Volatile Organic Compounds) and is non-irritating to eyes and skin
- H<sup>2</sup>O<sup>2</sup> decomposes to water and oxygen
- Wipes remain thoroughly wet and do not dry out quickly so contact times are met
- Wipes are durable and strong with no lint residue

SCI-OPTIMW1 Wipes Tub of 160 SCI-OPTIM1L Liquid 1 Ltr SCI-OPTIM4L Liquid 4 Ltrs

### **Plastex Clinical Detergent**

### **Features and Benefits:**

- Plastex is a clinical detergent specially formulated for cleaning even the most sensitive plastics
- Applications include:
  - Curing Lights Overhead Light Covers
  - Oral Cameras Chairs & Upholstery

MJ-PLS500 Liquid 500mls

### \* MAJAC



### Alcohol-Based Cleaning & Disinfection Solutions & Wipes

### Mikrozid Jumbo Wipes

### **Features and Benefits:**

- Blend of two alcohols for a very high performance against viruses within 1 min contact time
- Cleaning & disinfection, 2 complete actions with 1 product
- Added surfactants for boosted cleaning performance
- Effective against bacteria (inc. TB), fungi and viruses (inc. HBV, HIV, HCV and Rota virus)

SK-109158 Tub of 200 SK-10915 Refill 200 Bag



schülke -}

### FD 333 Wipes

#### Features and Benefits:

- Scent-free
- Fast disinfectant effect
- Fast, residue-free drying
- Practically dispensed from a dispenser tub

DU-CDF333A0150 Wipes Tub DU-CDF333A1250 Wipes Refill



### Caviwipes® Disinfecting Wipes K△VO ■ Kerr : Cavicide® Disinfecting Liquid

#### **Features and Benefits:**

- Does not contain aldehydes, phenols, bleach, or other toxic chemicals
- Large and thick, 7-layer, non-woven, textured cloth
- Low percentage of alcohol (17%) allows for a fungicidal, bactericidal, virucidal and tuberculocidal disinfection in only 3 minutes



KE- 1311011 Regular Wipes 15.2cm x 17.2cm Tub of 160 KE-1311501 X- Large Wipes 30cm x 25cm Tub of 65

### K4VO Kerr

### **Features and Benefits:**

- Broad-spectrum medical device surface disinfectant effective against viruses and bacteria including TB, HBV, HCV, HIV, MRSA, VRE and fungi
- Easy to use disinfectant for non-porous hard surfaces subject to standard infection control procedures

KE- 060210 Liquid 709mls KE-060104 Liquid 3.8Ltrs



### **Environmental Cleaning**

### Terralin® Protect Cleaner and **Disinfectant for Surgery Floors**

### **Features and Benefits:**

- Excellent cleaning performance
- Broad microbiocidal effectiveness
- Pleasant and fresh fragrance
- Excellent material compatibility including sensitive materials such as acrylic glass (incubators)
- Economic in use from 0.25% concentration
- CE marked in line with European directives

SK- 181621 Liquid 2 Ltrs

### schülke ->



### Clini-Sorb® Spill Management

### **Features and Benefits:**

- Super absorbent absorbs 30 times its own weight immediately
- Solidifies fluids into a gel quickly and effectively for easy disposal
- Retains fluid under external pressure, unlike sponges, tissues or powdered minerals
- Safe contains an anti-bacterial agent for safe handling

MJ- CS200 200gms





### **Impression Decontamination**

### Perform® ID Cleaner and schülke -> **Disinfectant for Impressions**



### **Features and Benefits:**


- Modern, Aldehyde free formulation
- Tested material compatibility and dimensional accuracy for polyether, silicone and alginate impressions
- · Comprehensive microbiological effectiveness
- Excellent cleaning ability, dissolves dirt, blood and saliva
- Safe and simple to use

SK-122305 Tub 900gm SK-144341 Bath 3 Ltrs SK-122313 30 x 40gm Sachets

### **Instrument Cassettes**

### **Exam Cassette - 5 Instuments**





### **Dimensions:**

- Exterior Dimensions: 3" x 8" x 1.25" (76mm x 203mm x 32mm)
- Holds 5 instruments

HF-IM6050 Red HF-IM6055 Yellow HF- IM6051 Gray HF-IM6056 Purple HF-IM6052 White HF-IM6058 Blue HF-IM6053 Orange HF-IM6059 Green

### Small Cassette - 8 & 10 Instuments





### **Dimensions:**

• Exterior Dimensions: 5.5" x 8" x 1.25" (140mm x 203mm x 32mm)

| Holds 10          |
|-------------------|
| Instruments       |
| HF-IM5100 Red     |
| HF-IM5101 Gray    |
| HF-IM5102 White   |
| HF- IM5103 Orange |
| HF-IM5105 Yellow  |
| HF-IM5106 Purple  |
| HF-IM5108 Blue    |
| HF-IM5109 Green   |

### Holds 8 instruments with Accessory Area HF-IM5080 Red HF-IM5081 Gray HF-IM5082 White HF-IM5083 Orange HF-IM5085 Yellow

HF-IM5086 Purple HF-IM5088 Blue

HF-IM5089 Green

### Large Cassette - 16 & 20 Instruments with Accessory Area





### **Dimensions:**

- Exterior Dimensions: 8" x 11" x 1.25" (203mm x 280mm x 32mm)
- Holds 16 instruments with an accessory area
- Also available in 8" x 10.5" x 1.25" size (IM416 x Q) to fit unique ultrasonic tanks

# Holds 20 Instruments HF-IM4200 Red HF-IM4201 Gray HF-IM4202 White HF-IM4203 Orange HF-IM4205 Yellow HF-IM4206 Purple HF-IM4208 Blue HF-IM4209 Green

### Holds 16 instruments with Accessory Area HF-IM4160 Red HF-IM4161 Gray HF-IM4162 White HF-IM4163 Orange HF-IM4165 Yellow HF-IM4166 Purple HF-IM4168 Blue HF-IM4169 Green

## Signa-Stat Cassette - 12 Instruments with Accessory Area





### **Dimensions:**

- Exterior Dimensions: 6.5" x 10.5" x 1.25 (165mm x 267mm x 32mm)
- Holds 12 instruments with an accessory area
- Fits in SciCan® STATIM® or the Midmark M3 autoclave

HF-IM8120 Red HF-IM8125 Yellow HF-IM8121 Gray HF-IM8126 Purple HF-IM8122 White HF-IM8128 Blue HF- IM8123 Orange HF-IM8129 Green

### Symmetry IQ® Handpiece Cassette





### **Dimensions:**

- Exterior Dimensions Single: 1.5" x 8" x 1.25" (38mm x 203mm x 32mm)
- Holds 1 Symmetry IQ® Handpiece
- Exterior Dimensions: Double 3" x 8" x 1.25" 76mm x 203mm x 32mm)
- Holds 2 Symmetry IQ<sup>®</sup> Handpieces

HF-IM30IQ8 Single HF-IM60IQ8 Double

### Symmetry IQ® Wrench Cassette





### **Dimensions:**

- Exterior Dimensions: 5.63" x 1.95" x 1.77" (143mm x 50mm x 45mm)
- Accommodates 3 Symmetry Piezo Tips in Guardian Wrenches, Blue Rail
- Wrenches not included

HF-IMIQ38 Blue

### **Ultrasonic Insert / Handpiece Cassette**







### **Dimensions:**

- Exterior Dimensions **Small:** 1.5" x 8" x 1.25" (38mm x 203mm x 32mm)
- Holds 2 Ultrasonic Inserts or 1 Handpiece
- Exterior Dimensions: Large 3" x 8" x 1.25" (76mm x 203mm x 32mm)
- Holds up to 3 high-speed Handpieces or 5 Ultrasonic Inserts

#### **Small**

HF-IM30HU0 Red HF-IM30HU1 Gray HF-IM30HU2 White HF-IM30HU3 Orange HF-IM30HU5 Yellow HF- IM30HU6 Purple HF-IM30HU8 Blue HF-IM30HU9 Green

### Large

HF-IM60H0 Red HF-IM60H1 Gray HF-IM60H2 White HF-IM60H3 Orange HF-IM60H5 Yellow HF-IM60H6 Purple HF- IM60H8 Blue HF-IM60H9 Green

### **Oral Surgery Cassette**





### **Dimensions:**

- Exterior Dimensions: 8" x 14.5" x 1.25" (203mm x 368mm x 32mm)
- Accommodates 12 instruments, 4 elevators, 3 hinged instruments in the lid and a separate accessory area
- Has a built-in needle cap holder and slot for a winged syringe
- Fits in most 10" chambers or larger sterilisers, most instrument washers and larger ultrasonics

HF-IMOS0 Red HF-IMOS5 Yellow HF-IMOS1 Gray HF-IMOS6 Purple HF-IMOS3 Unite HF-IMOS9 Green HF-IMOS9 Green

### Double-Decker® Cassette 8 Instruments





#### **Dimensions:**

- Exterior Dimensions: 3" x 8" x 1.5" (76mm x 203mm x 38mm)
- Holds 8 instruments

HF-IM9080 Red HF-IM9085 Yellow HF-IM9081 Gray HF-IM9086 Purple HF-IM9082 White HF-IM9088 Blue HF-IM9083 Orange HF-IM9089 Green

### Double-Decker® Cassette 14 Instruments





#### **Dimensions:**

- Exterior Dimensions: 4.5" x 8" x 1.5" (114mm x 203mm x 38mm)
- Ortho Holds up to 4 pliers, a Mathieu Needle holder and two instruments

| Holds 14 In | strument |
|-------------|----------|
| HF-IM9140   | Red      |
| HF-IM9141   | Gray     |
| HF-IM9142   | White    |
| HF-IM9143   | Orange   |
| HF-IM9145   | Yellow   |
| HF-IM9146   | Purple   |
| HF-IM9148   | Blue     |
| HF-IM9149   | Green    |

# Ortho HF-IM9140-OR Red HF-IM9141-OR Gray HF-IM9142-OR Whit HF-IM9143-OR Orange HF-IM9145-OR Yellow HF-IM9146-OR Purple HF-IM9148-OR Blue HF-IM9149-OR Green

### Double-Decker® Cassette 20 Instruments with Accessory Area





### **Dimensions:**

- Exterior Dimensions: 6" x 9.7" x 1.5" (152mm x 246mm x 38mm)
- · Holds 20 instruments with an accessory area
- Has a built-in needle cap holder and slot for a winged syringe

 HF-IM9200 Red
 HF-IM9205 Yellow

 HF-IM9201 Gray
 HF-IM9206 Purple

 HF-IM9202 White
 HF-IM9208 Blue

 HF-IM9203 Orange
 HF-IM9209 Green

# **Small Infinity Series Cassette 8 Instruments with Accessory Area**





### **Dimensions:**

- Exterior Dimensions 5.5" x 8" x 1.25" (140mm x 203mm x 3 mm)
- · Cassette holds 8 instruments with an accessory area

| HF-IMN5080  | Red        | HF-IMN5085 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN5081  | Gray       | HF-IMN5086 | Lavender |
| HF-IMN50811 | Ocean Blue | HF-IMN5087 | Black    |
| HF-IMN5082  | White      | HF-IMN5088 | Blue     |
| HF-IMN5083  | Orange     | HF-IMN5089 | Green    |
| HF-IMN5084  | Purple     |            |          |

### Small Infinity Series Cassette 10 Instrument





### **Dimensions:**

- Exterior Dimensions 5.5" x 8" x 1.25" (140mm x 203mm x 32mm)
- Cassette holds 10 instruments

| HF-IMN5100  | Red        | HF-IMN5105 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN5101  | Gray       | HF-IMN5106 | Lavender |
| HF-IMN51011 | Ocean Blue | HF-IMN5107 | Black    |
| HF-IMN5102  | White      | HF-IMN5108 | Blue     |
| HF-IMN5103  | Orange     | HF-IMN5109 | Green    |
| HF-IMN5104  | Purnle     |            |          |

# Large Infinity Series Cassette 12 Instrument with Accessory Area





### **Dimensions:**

- Exterior Dimensions 8" x 11" x 1.25" (203mm x 280mm x 32mm)
- Cassette holds 12 instruments with an accessory area

| HF-IMN4120  | Red        | HF-IMN4125 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN4121  | Gray       | HF-IMN4126 | Lavender |
| HF-IMN41211 | Ocean Blue | HF-IMN4127 | Black    |
| HF-IMN4122  | White      | HF-IMN4128 | Blue     |
| HF-IMN4123  | Orange     | HF-IMN4129 | Green    |
| HF-IMN4124  | Purple     |            |          |

## Large Infinity Series Cassette 16 Instrument with Accessory Area





### **Dimensions:**

- Exterior Dimensions 8" x 11" x 1.25" (203mm x 280mm x 32mm)
- Cassette holds 16 instruments with an accessory area

| HF-IMN4160  | Red        | HF-IMN4165 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN4161  | Gray       | HF-IMN4166 | Lavender |
| HF-IMN41611 | Ocean Blue | HF-IMN4167 | Black    |
| HF-IMN4162  | White      | HF-IMN4168 | Blue     |
| HF-IMN4163  | Orange     | HF-IMN4169 | Green    |
| HF-IMN4164  | Purple     |            |          |

## Large Infinity Series Cassette 20 Instrument





### **Dimensions:**

- Cassette holds 20 instruments
- Exterior Dimensions 8" x 11" x 1.25" (203mm x 280mm x 32mm)

| HF-IMN4200<br>HF-IMN4201 | Red<br>Gray | HF-IMN4205<br>HF-IMN4206 | Yellow<br>Lavender |
|--------------------------|-------------|--------------------------|--------------------|
|                          | Ocean Blue  | HF-IMN4207               | Black              |
| HF-IMN4202               | White       | HF-IMN4208               | Blue               |
| HF-IMN4203               | Orange      | HF-IMN4209               | Green              |
| HF-IMN4204               | Purple      |                          |                    |

### **Infinity Series Signa-Stat Cassette** 12 Instrument with Accessory Area





### **Dimensions:**

- Exterior Dimensions 6.5" x 10.5" x 1.25" (165mm x 267mm x 32mm)
- Cassette holds 12 instruments with an accessory area and fits in the SciCan STATIM or the Midmark M3 autoclave

| HF-IMN8120<br>HF-IMN8121<br>HF-IMN81211 | Red<br>Gray      | HF-IMN8125<br>HF-IMN8126<br>HF-IMN8127 | Yellow<br>Lavender<br>Black |
|---|------------------|--|-----------------------------|
| HF-IMN8122                              | White            | HF-IMN8128                             | Blue                        |
| HF-IMN8123<br>HF-IMN8124                | Orange<br>Purple | HF-IMN8129                             | Green                       |

### Infinity Series Double-Decker Cassette Hu-Friedy **8 Instrument**





### **Dimensions:**

- Exterior Dimensions 3" x 8" x 1.5" (76mm x 203mm x 38mm)
- Cassette holds 8 instruments

| HF-IMN9080  | Red        | HF-IMN9085 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN9081  | Gray       | HF-IMN9086 | Lavender |
| HF-IMN90811 | Ocean Blue | HF-IMN9087 | Black    |
| HF-IMN9082  | White      | HF-IMN9088 | Blue     |
| HF-IMN9083  | Orange     | HF-IMN9089 | Green    |
| HF-IMN9084  | Purple     |            |          |

### **Infinity Series Double-Decker Cassette** 14 Instrument





#### **Dimensions:**

- Exterior Dimensions 4.5" x 8" x 1.5" (114mm x 203mm x 38mm)
- Cassette holds 14 instruments

| HF-IMN9140   | Red        | HF-IMN9144 | Purple   |
|--------------|------------|------------|----------|
| HF-IMN9140-0 | Ocean Blue | HF-IMN9145 | Yellow   |
| HF-IMN9141   | Gray       | HF-IMN9146 | Lavender |
| HF-IMN91411  | Ocean Blue | HF-IMN9147 | Black    |
| HF-IMN9142   | White      | HF-IMN9148 | Blue     |
| HF-IMN9143   | Orange     | HF-IMN9149 | Green    |

### **Infinity Series Double-Decker Cassette 20 Instrument with Accessory Area**





### **Dimensions:**

- Exterior Dimensions 6" x 9.7" x 1.5" (152mm x 246mm x 38mm)
- Cassette holds 20 instruments with an accessory area

| HF-IMN9200  | Red        | HF-IMN9205 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN9201  | Gray       | HF-IMN9206 | Lavender |
| HF-IMN92011 | Ocean Blue | HF-IMN9207 | Black    |
| HF-IMN9202  | White      | HF-IMN9208 | Blue     |
| HF-IMN9203  | Orange     | HF-IMN9209 | Green    |
| HF-IMN9204  | Purple     |            |          |

### Infinity Series Exam Cassette 5 Instrument





### **Dimensions:**

- Exterior Dimensions 3" x 8" x 1.25" (76mm x 203mm x 32mm)
- Cassette holds 5 instruments

| HF-IMN6050  | Red        | HF-IMN6055 | Yellow   |
|-------------|------------|------------|----------|
| HF-IMN6051  | Gray       | HF-IMN6056 | Lavender |
| HF-IMN60511 | Ocean Blue | HF-IMN6057 | Black    |
| HF-IMN6052  | White      | HF-IMN6058 | Blue     |
| HF-IMN6053  | Orange     | HF-IMN6059 | Green    |
| HF-IMN6054  | Purple     |            |          |

## Large Infinity Series Q Cassette 16 Instrument with Accessory Area





- Exterior Dimensions 8" x 10.5" x 1.25" to fit unique ultrasonic tanks
- · Cassette holds 16 instruments with an accessory area

HF-IMN41600 Red HF-IMN41650 Yellow
HF-IMN416110 Ocean Blue HF-IMN41660 Lavender
HF-IMN41610 Gray HF-IMN41670 Black
HF-IMN41620 White HF-IMN41680 Blue
HF-IMN41630 Orange HF-IMN41690 Green
HF-IMN41640 Purple

# Infinity Series Oral Surgery Cassette with Accessory Area





- Exterior Dimensions 8" x 14.5" x 1.25" (203mm x 368mm x 32mm)
- Accommodates 12 instruments, 4 elevators, 3 hinged instruments in the lid and has a separate accessory area

| HF-IMNOS0  | Red        | HF-IMNOS5 | Yellow   |
|------------|------------|-----------|----------|
| HF-IMNOS1  | Gray       | HF-IMNOS6 | Lavender |
| HF-IMNOS11 | Ocean Blue | HF-IMNOS7 | Black    |
| HF-IMNOS2  | White      | HF-IMNOS8 | Blue     |
| HF-IMNOS3  | Orange     | HF-IMNOS9 | Green    |
| HF-IMNOS4  | Purple     |           |          |

### Infinity Series Double-Decker Small Orthodontic Cassette





### **Dimensions:**

- Exterior Dimensions 4.5" x 8" x 1.5" (114mm x 203mm x 38mm)
- Cassette holds up to 4 pliers, 1 Mathieu needle holder, and 2 hand instruments

| HF-IMN91411-0 | Ocean Blue | HF-IMN9145-0 | Yellow   |
|---------------|------------|--------------|----------|
| HF-IMN9141-0  | Gray       | HF-IMN9146-0 | Lavender |
| HF-IMN9142-0  | White      | HF-IMN9147-0 | Black    |
| HF-IMN9143-0  | Orange     | HF-IMN9148-0 | Blue     |
| HF-IMN9144-0  | Purple     | HF-IMN9149-0 | Green    |

Note: Instruments are not included with cassettes

**Dimensions:** 

### Stainless Steel Instrument Cassette - 5 Instrument









### **Dimensions:**

- Exterior Dimensions: 72.5mm x 182.5mm x 32mm includes 3 blue silicone rails
- Instrument cassette is fully autoclavable at 135°C and is washer disinfector safe

HS-9008543

### Stainless Steel Instrument **Cassette - 10 Instrument**





### **Dimensions:**

- Exterior Dimensions: 145mm x 182.5mm x 32mm) includes 3 blue
- Instrument cassette is fully autoclavable at 135°C and is washer disinfector safe

HS-9008544




### **Dimensions:**

- Exterior Dimensions: (279 x 182,5 x 32 mm) includes 6 blue silicone rails
- Instrument cassette is fully autoclavable at 135°C and is washer disinfector safe

HS-9008545

### Coloured Rails for Henry Schein Henry Schein **Stainless Steel Cassettes**





HS-9008550 Green Silicone Rails 3 Pieces For 5 Cassette HS-9008551 Orange Silicone Rails 3 Pieces For 5 Cassette HS-9008552 Red Silicone Rails 3 Pieces For 5 Cassette HS-9008553 Yellow Silicone Rails 3 Pieces For 5 Cassette HS-9008554 Blue Silicone Rails 3 Pieces For 10 & 16 Cassette HS-9008555 Green Silicone Rails 3 Pieces For 10 & 16 Cassette HS-9008556 Orange Silicone Rails 3 Pieces For 10 & 16 Cassette HS-9008557 Red Silicone Rails 3 Pieces For 10 & 16 Cassette HS-9008558 Yellow Silicone Rails 3 Pieces For 10 & 16 Cassette

### **Cleaning Accessories**

### Miltex Instrument Brush





### **Features and Benefits:**

- For cleaning hard to reach debris from instruments prior to ultrasonic cleaning or sterilsation
- Nylon bristles won't scratch or harm delicate instrument surfaces
- Plastic handle measures 7" long
- Autoclavable

HS953-9664 18.4cm, Nylon Bristles, Each

### Miltex Bur Brush





### **Features and Benefits:**

- Ideal for serrations, rasps, files and burs
- Plastic handle measures 7" long
- Autoclavable

HS953-2703 18.4cm, Steel Bristles, Each

### **Hand Scrub Brush**





### Features and Benefits:

- Scrubs instruments of debris before sterilising or ultrasonic cleaning
- Stiff bristles scrub the difficult areas where debris can hide
- Autoclavable, aseptic high impact plastic

HS100- 8055 Nylon, 4" x 1.5" x 1.25", Bristle Length .5"

### **Bottle Cleaning Brush**





### **Features and Benefits:**

- · High grade nylon for performance and durability
- Designed to enter the narrow mouth of tubes and bottles
- Expands to allow the brush to thoroughly clean the bottom and corners of the glassware without scratching
- White nylon with twisted wire handle with hanging loop

HS100-6161 Large, Pkt of 6

### **Aquasorb Lint Free Cloth**





### **Features and Benefits:**

- Lint Free: Knitted nylon/polyester microfibre does not shed lint
- Washable: Reprocess manually, or in a commercial linen service
- Autoclavable: Steam sterilisable to 130°C. Only when used as a tray liner
- Super-absorbent: Fibre will absorb 315% of its own weight

MJ-AQ1 Small, 55 x 22.5cm, Each MJ-AQ2 Medium, 55 x 42.5cm, Each MJ-AQ3 Large 65 x 50cm, Each

### Do you really know the Rotary Endo File Reprocessing Guidelines?

As recommended by NZ Dental Council Infection Prevention & Control Practice Standard 2016\*

Rotary nickel titanium endodontic files can be reprocessed facilitating subsequent reuse if the outlined protocol is followed.

 Insert the files into a scouring sponge soaked with 0.2% chlorhexidine gluconate aqueous solution immediately after use.
 The rubber stoppers should be discarded and care should be taken to avoid the material adherent to the file from drying out.

- Clean the files by using 10 vigorous 'in-and-out' strokes in the sponge
- Soak the files in a suitable enzymatic cleaning solution for 30 minutes
- Ultrasonically clean the files for 15 minutes in the enzymatic cleaning solution
- Drain and rinse in running water for 20 seconds
- The files can now be sterilised.



### **Endodontic File Reprocessing Products**

### **Chlorhexidine**





### **Features and Benefits:**

 Ideal cleaning solution for compliance with the recommended rotary endodontic file reprocessing guidelines as outlined in the current Infection Prevention Code of Practice

CP-478547 Chlorhexidine 0.2%, 200ml

### ENDOFOAM-PVP Endodontic Foam



**ADM** 

### **Features and Benefits:**

- Triple-layered consisting of a thin low density foam layer, an extended 4mm thick scourer layer and a 25mm thick high density foam layer
- ENDOFOAM-PVP is multiple use when used within a sterilising environment for the cleaning of multiple nickel titanium endodontic files as part of the Parashos Validated Protocol
- Can be used on standard hand files

A7-ENDOFOAMPVP Pack of 49

### FileFoam<sup>™</sup> File Cleaner





### **Features and Benefits:**

- Patented, triple-layer foam construction consisting of a thin, low density foam layer, an extended 4mm thick scour layer, and a 25mm thick high density foam layer
- FileFoam™ provides a triple-action solution; cleaning, scouring and disinfecting endodontic instruments as part of the Parashos Validated Protocol
- Can be autoclaved up to 137°C if required, prior to use

HS-7910308 Pack of 49

<sup>\*</sup>Appendix F Reprocessing rotary nickel titanium endodontic files

### **Instrument Cleaning and Disinfecting Solutions**

### **Clinidet® Clinical Detergent**

### **Features and Benefits:**

- Excellent detergency with rapid wetting, detergency is critical to removing soil from the surface being cleaned
- Rapid wetting is important for ultrasonic cleaning
- Very low foaming under dynamic conditions - a major safety feature enabeling the operator to view sharps being cleaned
- Corrosion inhibitor incorporated into the solution. Helps protect delicate and high quality instruments
- Synergistic blend of multi-functional sequesterents that greatly improve the solubilisation of blood, proteins and fats. Helps suspend soil in solution and prevent redeposition of soil onto the cleaned surface
- Excellent performance in hard water
- Safe to use, store and transport
- The Benchmark for clinical detergents

MJ-CL5 5 Litre MJ-CL8P 8ml Pump For 5L Container MJ-CLE500 Easypack 500ml MJ-CLDIS 500ml Office Dispenser







### **Optim 28**

### **Features and Benefits:**

- Economical
- Safe for staff and the environment
- Fast acting 5-minute cleaning time
- Bacteriostatic



**SciCan** 

### Gigasonic<sup>®</sup> Ultrasonic Cleaner and Disinfecting Solution

### **Features and Benefits:**

- Broad spectrum microbiological efficacy
- Aldehyde free
- Pleasant fresh fragrance

SK- 112401 2 Ltrs SK-112402 5 Ltrs



### UC30 Biosonic® General **Purpose Cleaning Solution**

#### **Features and Benefits:**

- UC30 is a non-ionic detergent and is used for general cleaning
- Removes most organic soils like blood, saliva and debris

W- UC30 473ml, Diluted Yield 29 Ltrs

**## COLTENE** 

### **UC31 Biosonic® Cleaner and Rust Inhibitor Solution**

### **Features and Benefits:**

- UC31 combines the proven cleaning action of a BioSonic® solution with corrosion protection in the ultrasonic bath
- Used for general cleaning

W- UC31 20 x 15ml, Diluted Yield 76 Ltrs

### **##** COLTENE

### **UC32 Biosonic® Enzymatic Solution**

### **Features and Benefits:**

- UC32 combines the non-ionic detergency of the UC30 with enzymes that act to break down stubborn dirt particles
- Effectively breaks down and loosens debris hardened onto instruments
- Minimises hand scrubbing

W- UC32 236ml, Diluted Yield 29 Ltrs



**## COLTENE** 

SCI-28004 Litre

### **Ultrasonic Cleaners**

### **BioSonic Ultrasonic Cleaning System**

### **UC150\* Ultrasonic Cleaning System with Basket**

### **Features and Benefits:**

- The BioSonic UC150 can sit directly on the countertop or be recessed
- The USB interface allows for easy downloading or capture of cleaning data
- The BioSonic UC150 offers a 6 ltr stainless steel tank
- A sleek wear resistant housing that will not dent or stain

Tank Capacity 6.6 litres Fill Capacity 5.7 litres Overall Dimensions ( $I \times w \times h$ )  $34.5 \times 27.6 \times 35.5$  cm Tank Dimensions ( $I \times w \times d$ )  $28.7 \times 16.0 \times 17.8$  cm



**#** COLTENE

### 



UC150 Accessories CW-60022040 Basket CW-60022041 600ml Beaker Holder CW-60022042 Accessory Kit, 1 x Basket 2 x Beakers & Holders

W-UC325 Cassette Rack for UC300



W-UC311 Half FingerGuard Basket with 2 Section Inserts for UC300



W-UC310 Large Basket for UC300



Beaker 600ml W-UC153 Positioning Cover

W-UC53



### CW-60021990



### UC50DB\* Ultrasonic Cleaning System with Basket

Tank Capacity 3.2 litres Fill Capacity 2.0 litres Overall Dimensions ( $I \times w \times h$ ) 27.3  $\times$  21.9  $\times$  20.3 cm Tank Dimensions ( $I \times w \times d$ ) 19.7  $\times$  14.6  $\times$  10.5 cm

### W-UC50DB240

\*2 Year Warranty

### UC125\* Ultrasonic Cleaning System With Basket

Tank Capacity 4.8 litres Fill Capacity 3.8 litres Overall Dimensions ( $I \times w \times h$ )  $42 \times 26.7 \times 35$  cm Tank Dimensions ( $I \times w \times d$ )  $28.6 \times 14.3 \times 12.7$  cm

W-UC125B-240



Tank Capacity 14.4 litres Fill Capacity 11 litres Overall Dimensions ( $I \times w \times h$ ) 48.2  $\times$  30.5  $\times$  31.1 cm Tank Dimensions ( $I \times w \times d$ ) 35.6  $\times$  22.9  $\times$  17.1 cm

W-UC300-230

### ULTICLEAN Ultrasonic Bath

#### **Features and Benefits:**

- State of the art 4 litre capacity ultrasonic bath with optional heat
- Wave technology gives a uniform cleaning pattern that eliminates hot spots and weak areas, giving constant power and consistent output
- Smooth, quiet operation and easy to clean housing
- Hinged automatic locking lid prevents opening during cycle
- Easy to use built in drain facility
- Touch pad controls with 60 minute timer



BA-BATH 4L Ultrasonic Bath BA-BATHBASKET Accessory Basket Kit

### **Handpiece Cleaners**

### **QUATTROcare PLUS Handpiece Maintenance Unit**

#### **Features and Benefits:**

- QUATTROcare PLUS offers excellent ease of use and maximum efficiency for optimal care of your handpieces
- Time Saving: 1 minute processing cycle for 4 instruments
- Safety: Program-controlled automatic care cycle
- Efficiency: Cleaning with both integrated detergent and foaming oil (mechanical expansion/low temperature boiling). Effective and optimal lubrication through automated, factory pre-set dosing
- Universal: May be confidently used with all established handpiece brands





KV-1.008.3805\* QUATTROcare KV-1.005.3843 QUATTROcare Spray

\*Includes 1 x KaVo Highspeed and 3 x Lowspeed fittings.

### **Lubrina Handpiece Maintenance Unit**















### **Features and Benefits:**

- Services 4 handpieces in 100 seconds
- Built in chuck maintenance system for Morita handpieces only
- Purges and lubricates most types of handpieces
- Unique dual lubrication system
- Step air blow function
- Accommodates most handpiece attachments and spray can brands
- Atomised spray for maximum penetration

M-181244\* Morita Lubrina M-5010212 Morita AR Sprav M-5071340 Morita LS Spray

\*Includes any 4 x fittings of your choice.



### **Features and Benefits:**

- Rotates 4 handpieces during cleaning cycle. This allows oil to spread to every part of the handpiece
- Filtered air and purging action to flush out all residual solution
- Solution level indicator. With full tank, up to 3,500 handpieces are cleaned and lubricated
- · Chuck cleaning function ensures that the chuck is thoroughly cleaned and lubricated
- Economical, efficient and timesaving

NK-Y10012787 NSK iCare 3\* NK-Z016130 iCare Maintenance Oil, 1 Litre

\*Includes 2 x Lowspeed Fittings, 2 x Highspeed fittings of your choice.

### **Sterilisation Pouches**

### Maxima Self Seal Sterilisation Pouch





#### **Features and Benefits:**

- Superior quality pouch with dual indicators that comply with the recommendations of the CDC
- Preprinted internal and external indicators change colour when processed with steam and EO
- Made with medical-grade paper and an easy-to-open thumb notch
- Lead free and latex free

HS- 9006153 57 x 100mm (2.25 x 4") Box 200 HS- 9006154 70 x 230mm (2.75 x 9") Box 200 HS- 9006155 90 x 140mm (3.5 x 5.5") Box 200 HS- 9006156 90 x 255mm (3.5 x 10") Box 200 HS- 9006157 133 x 255mm (5.2 5x 10") Box 200 HS- 9006158 190 x 330 mm (7.5 x 13") Box 200 HS- 9006159 230 x 380 mm (9 x 15") Box 200

## Bagette® Self-Seal Sterilisation Pouch





#### **Features and Benefits:**

- Constructed of high quality materials and thick seals to help prevent instrument and cassette tear through
- 4 sizes have been specifically designed to optimally fit Hu-Friedy cassettes while 4 additional sizes are available for all other sterilisation needs
- Built in internal and external process indicators provide time and temperature verification
- For single-use only

HF-IMS123

### Self Seal Sterilisation HALAS Pouch



HALAS Sterilisation Pouches Self Sealing with internal and external indicator to ensure temperatures reach the proper level inside the pouch.

HS-0100 57 x 100mm (2.25" x 4") Box 200 HS-0200 57 x 70mm (2.25" x 2.75) Box 200 HS-0500 70 x 255mm (2.75" x 10") Box 200 HS-2000 90 x 255mm (3.5 x 10") Box 200 HS-3000 133 x 255mm (5.25" x 10") Box 200 HS-4000 190 x 330mm (7.5" x 13") Box 200

### PeelVue+ Sterilisation Pouch



#### **Features and Benefits:**

- With Closure Validation (patent pending)
- Easy to close and peel open, even with gloves
- Self-sealing
- Patented internal and external Processing Indicator
- With practical colour coding

VR-316\_\_\_

### Assure Plus™ Self-Seal Sterilisation Pouch





### **Features and Benefits:**

- Unique peel-off indicator
- No need to purchase secondary indicator strips
- Dual external heat and gas sterilisation indicators
- Quality seals for added strength and a clean, fibre-free peel

SU-830\_\_

### **Sure-Check Sterilisation Pouches**







CDC infection control guidelines make it clear: to ensure effective sterilisation, items must be exposed to all 3 sterilisation criteria - time, temperature and pressure (steam).

### **Features and Benefits:**

- Highest sterility assurance Built-in internal/external multi-parameter indicators only change colour when all 3 criteria for sterilisation have been met - unlike traditional sterilisation pouches with a single process indicator that measures only temperature.
- Simplified process Eliminates the extra step of adding separate multi-parameter indicators to each pouch.
- Meets CDC Guidelines for Disinfection and Sterilisation in Healthcare Facilities, 2008.
- Meets ANSI/AAMI ST79 Comprehensive Guide to Steam Sterilisation and Assurance in Health Care

CX-SC

### Duo-Check CROSSTEX® Sterilisation Pouches



### **Features and Benefits:**

- Secure Seal
- Medical Grade Paper
- Blue Tinted Transparent Film
- Wide Self-Seal Adhesive Strip
- Perforated Fold
- Complete Paper-Film
- Process Indicators

CX-SC\_

### **Sterilisation Wrap**

### **IMS Cassette Wrap**

### **Features and Benefits:**

- Engineered to withstand the harsh conditions of sterilisation
- Unique compilation of textiles is designed to be strong and soft for easy wrapping and durability
- · Low-lint, memory free wrap
- · Available in two grades to meet your needs

Autoclave Wrap offers the best combination of softness and added strength.

HF-IMS-1210  $24'' \times 24''$  (610 mm x 610 mm) For Large Cassettes, 500 Sheets

HF-IMS-1212  $\,$  20"  $\,$  x 20" (508 mm x 508 mm ) For Large Cassettes, 500 Sheets

HF-IMS-1211 15" x 15" (380 mm x 380 mm) For Small Cassettes, 1,000 Sheets

HF-IMS-1213 12" x 12" (305 mm x 305 mm) For Small Cassettes, 1,000 Sheets Universal Wrap provides great performance at a lower cost.

HF-IMS-1214  $24" \times 24"$  (610 mm x 610 mm) For Large Cassettes, 500 Sheets

HF-IMS-1216 20" x 20" (508 mm x 508 mm) For Large Cassettes, 500 Sheets

HF-IMS-1215 15" x 15" (380 mm x 380 mm) For Small Cassettes, 1,000 Sheets

HF-IMS-1217 12" x 12" (305 mm x 305 mm) For Small Cassettes, 1,000 Sheets




How the best perform


### **Kimguard Sterile Wrap**

### **Features and Benefits:**

- POWERGUARD\* technology for enhanced microbial barrier protection
- Pre-vacuum steam sterilant penetration
- Bacterial Filtration Efficiency (BFE) from 98.9% to 99.9%
- Low linting
- Made in the USA
- Tear, flame, lint and abrasion resistant

KC-10824 60 x 60cm Box 500









### **Sterilisation Accessories**

### Van R Bur Holder

### **Features and Benefits:**

- Magnetic bur block for accessible storage of all types of rotary cutting instruments
- Made of sturdy polypropylene plastic
- Can be ultrasonically cleaned and are compatible with most chemical disinfectants

VR- 182\_ Small, (8 Burs), Assorted Colours VR- 183\_ Large (24 Burs), Assorted Colours



### Miltex Mesh Diamond Bur Holder



### **Features and Benefits:**

- Securely holds diamond burs and other small items
- For use in autoclaves and ultrasonics
- Stainless steel

HS953-2143

### **Aluminium Bur Block**

### **Features and Benefits:**

- Designed for individual composition of instruments for specific treatment methods and preparation techniques
- Ergonomic swivel lid protects instruments against damage
- Can be disinfected, cleaned and sterilised



### **Stainless Steel Bur Block**

### **Features and Benefits:**

- Instruments in this bur block can be cleaned and disinfected in instrument bath, ultrasonic bath or thermo disinfector
- Blue silicone plugs hold instruments safely in place, even when upside down
- Autoclavable



K-A100S (15 Burs) 10 FG/5RA

K-9949L3 (24 Burs either FG/RA) including Silicone Plugs

### **Performance Testing**

### **Air Removal and Steam Penetration Tests**

An Air Removal Test is required each morning before processing of your first load.

### Valisafe Helix Test Device

Air Removal Device

 Reads on completion of the cycle to ensure that your Type B table top steriliser is performing correctly. Has the capability to provide you with effectively sterilised items in your practice.



### **TST Helix Test Device**

Air Removal Device

• The Browne TST Control Helix is a Hollow Load Process Challenge Device (PCD) and has been developed and validated for testing the air removal (steam penetration) capability of small Type B steam sterilisers.



Browne

ST-7BR00001 Includes 250 Chemical Indicators

### **Bowie Dick Helix Test Set**



Air Removal Device

- Conforms to ISO 11140 1 class 2 EN 867 part 5 class B
- Indicator marked with red printing for identification
- Clear transition of colour to pink
- Easy to interpret easy to document
- Non-Toxic / Lead-fFree.



MJ-S1010020250 Includes 250 Chemical Indicators

### Statim Quad Helix Device\*

Air Removal Device - Cassette

• This process challenge device has been designed as a monitoring device for use within a specific cycle, HOLLOW / WRAPPED (S) 134°C 3.5 min.



SCI-01108341 Quad Helix for Statim SCI-99108332 Sterilisation Emulators 250 Pcs/Box

\*For use only with the Statim Sterilisers. Device included with purchase of machine.

All performance testing devices shown above include Class 2 chemical indicators.

### **Process Indicators**

Visual indication that load has been processed only.

### Two Line Label Gun 💥 MAJAC with Steam Indicator





MJ-S6020020001 Label Gun MJ-S60201190003 Roll of 750 Indicators

### **Comply™ Lead Free Steam Indicator Tape**



An exposure monitor that securely seals packs, and provides a visual assurance that the package has been exposed to the steam sterilisation process.

MD-1322- 18 & 24mm

### **IMS Autoclave Tape**





Strong tape that can be written on and visually differentiates processed from non-processed cassettes.

HF-IMS1250 3/4" (19 mm), 55m per roll

### **Cycle Monitoring**

### Multi Parameter Indicators - Class 4. 5. 6

For machines without recording capabilities.

### Comply™ SteriGage® 31M




Class 5 Steam Chemical Integrators

 Use for pack control monitoring of all 118-138°C (245-280°F) steam sterilisation cycles.



MD-1243B 100 Integrators

### Comply™



 For use in determining whether the sterilisation process conditions were met inside each pack. Can be used in all 121 -134°C (250-273°F) sterilisation cycles.

MD-1250 240 Strips

### **STEAMPlus**<sup>™</sup>



Class 5 Integrators

• Used with steam sterilisation cycles to provide a definitive pass or fail result allowing for the release of all non-implant loads before the results of the spore test are received. Moving-front technology provides an immediate readout for use in gravity, prevacuum, and flash cycles. Enhance sterility assurance by monitoring steriliser performance.



CX-SSI-100 Box of 100

### Comply<sup>™</sup> Steam-Clo **31** Indicator



Class 6 Chemical Indicator

 Monitors multiple parameters of time, temperature and presence of saturated steam to provide better accuracy of results.

MD-0104



### Sterintech MAJAC **Emulator**



Used to show the heat parameters for sterilisation have been achieved.

MJ-S106103500S Box 500

### ConFirm® 24 In-Office **Biological Monitoring System**



- Read out in just 24 hours
- Easy to use, easy to read results
- One of the most cost-effective methods of biological monitoring
- Universal BI ConFirm® 24 biological indicators can be processed in most incubators
- Built-in Ampoule Crusher



CX-CSBI25 25 Box CX-NDB-060E Dry Block Incubator For Use With 24 hr BI 220V

### **Biological Indicators**

To be used as part of the Validation process. (Only required annually if fully validated process is implemented)

### Attest™ Biological Indicators



PCD: Monitoring 250°F/121°C gravity and 270°F/132°C vacuum-assisted steam sterilisation processes.

MD-1261P Brown Pack of 25 MD-1261 Brown Pack of 100





 PCD: Use the Attest 1261 Biological Indicator to monitor: 1. 270°F (132°C) gravity displacement steam sterilisation cycles.

MD-1262P Blue Pack of 25 MD-1262 Blue Pack of 100

### Attest™ Steam Incubator







3M Attest Incubator that provides optimal incubation conditions for the accurate readout of 3M<sup>™</sup> Attest<sup>™</sup> Biological Indicators.

MD-I118 3M Attest MD-1266 Log Book

### **Instrument Sterilisation**

### **Statim G4 Autoclave**

### **Experience the Statim G4 Difference**

Gentleness you know but now with interactivity you've never seen before!

#### **Features and Benefits:**

- Protect option of pin-protected load release safety and traceability
- Collect integrated USB datalogger, cycle history never lost again
- Fast unmatched in cycle speed for fast instrument turnaround meaning less instrument inventory

SCI-G4152314 Statim 2000 G4 SCI-G4222313 Statim 5000 G4





### B Futura Steriliser Vacuum Chamber Autoclave



### **Mocoms finest performance!**



B Futura is an intelligent steriliser, destined to be an indispensable work tool. Only years of research and experience can produce such an exceptional result. A high performance and user-friendly product.

The new B Futura steriliser offers even faster performance, has an easy to use interface and enables full traceability of each sterilisation cycle. There are two chamber sizes available, 17 and 22 litres. All sizes have the same overall machine dimensions and feature a colour touch-screen LCD display.

ST-7MOCKP 17 Litre Steriliser ST-7MOCKG 22 Litre Steriliser

### B Classic Steriliser Vacuum Chamber Autoclave



### An assured success!



The B Classic was designed to make work easier by optimising time and costs. Exceptional quality and practicality ensure exceptional safety and reliability.

Easy to use, the B Classic provides full traceability of every single sterilisation cycle. There are two chamber sizes available, 17 and 22 litres. All sizes have the same overall machine dimensions and feature a monochromatic LCD display and keypad.

ST-7MOCKM 17 Litre Steriliser ST-7MOCKB 22 Litre Steriliser



### **HYDRIM C61w G4 Instrument Washer**

### **SciCan**

### A new wave in instrument washing & disinfection!

### **Features and Benefits:**

- Countertop or under the counter, steri-center installation.
- G4 web portal includes data storage, service logs and remote access capabilities
- Validated cycles to ensure repeatable washing performance
- Improved serviceability, reliability and minimal maintenance requirements
- Vivid colour touch screen
- Robust stainless steel exterior
- USB port for data storage
- Active hot air drying

SCI-C61WDD02



### **Tethys H10 Hybrid Disinfection Device**

### The evolution of disinfection!

The Tethys H10 Hybrid Disinfection Device expands the traditional role of the thermal washer disinfector, adding ultrasonic cleaning as an initial step and delivering instruments dry and ready for packaging and sterilisation at the end of a 35-minute cycle.

### **Features and Benefits:**

- Fast 35 minute cycle
- Thermal disinfection, ultrasound and drying in one device
- User-friendly color touch screen operation
- Internal memory as well as USB and Ethernet connectivity for paperless record keeping
- Delayed cycle start to perform vacuum an helix tests prior to the start of the work day









SCT-7M0CWD



### **Hu-Friedy Signature Series™ Tubs**





### **Features and Benefits:**

- Together, procedure tubs and cassettes create the most efficient and organised system for managing
  instruments and consumable products. The time savings achieved by using IMS (Instrument
  Management System) cassettes and procedure tubs will allow the staff to spend time on more value
  added activities in the dental office
- Signature Series Tubs are designed to match IMS cassette rail colours

### Signature Series™ Complete Tub Includes:

- Signature Series Tub
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover

Signature Series™ Complete Tub

- Divided Slide Tray
- · Locking Tub Cover

### Signature Series™ Tub

| HF-IMS-1400 | Red      | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
|-------------|----------|---|
| HF-IMS-1401 | Gray     | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
| HF-IMS-1402 | White    | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
| HF-IMS-1403 | Orange   | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
| HF-IMS-1405 | Yellow   | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
| HF-IMS-1406 | Lavender | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
| HF-IMS-1408 | Blue     | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |
| HF-IMS-1409 | Green    | 12 <sup>3</sup> /4" x 11 <sup>3</sup> /16" x 4" |

| g                    |   |
|----------------------|---|
| HF-IMS-1410 Red      | 12 <sup>1</sup> / <sub>4</sub> " x 10 <sup>7</sup> / <sub>8</sub> " x 2 <sup>3</sup> / <sub>4</sub> " |
| HF-IMS-1411 Gray     | 12 <sup>1</sup> /4" x 10 <sup>7</sup> /8" x 2 <sup>3</sup> /4"  |
| HF-IMS-1412 White    | 12 <sup>1</sup> / <sub>4</sub> " x 10 <sup>7</sup> / <sub>8</sub> " x 2 <sup>3</sup> / <sub>4</sub> " |
| HF-IMS-1413 Orange   | 12 <sup>1</sup> /4" x 10 <sup>7</sup> /8" x 2 <sup>3</sup> /4"  |
| HF-IMS-1415 Yellow   | 12 <sup>1</sup> /4" x 10 <sup>7</sup> /8" x 2 <sup>3</sup> /4"  |
| HF-IMS-1416 Lavender | 12 <sup>1</sup> /4" x 10 <sup>7</sup> /8" x 2 <sup>3</sup> /4"  |
| HF-IMS-1418 Blue     | 12 <sup>1</sup> / <sub>4</sub> " x 10 <sup>7</sup> / <sub>8</sub> " x 2 <sup>3</sup> / <sub>4</sub> " |
| HF-IMS-1419 Green    | 12 <sup>1</sup> / <sub>4</sub> " x 10 <sup>7</sup> / <sub>8</sub> " x 2 <sup>3</sup> / <sub>4</sub> " |















To order telephone 0800 808 855, 8am to 5.30pm OR freefax 0508 808 555, 24 hours www.henryschein.co.nz