



STERI EQUIPMENT OFFERS



1ST OCTOBER - 30TH NOVEMBER 2024

TETHYS H10+



SAVE
\$1,890



NOW ONLY \$9,300 ~~\$11,190~~

SCT-7MOCWD

Tethys H10 Plus makes the reprocessing process even simpler and more functional – decontamination, washing, thermal disinfection and drying are all performed in a single, rapid, automatic process.

1. DECONTAMINATION

Reduction of instruments' bacterial count through washing at low temperatures.

2. ULTRASOUND WASHING

Elimination of organic and inorganic residues from instruments. The ultrasound's force eliminates even the most resistant bioburden.

3. THERMAL DISINFECTION

Drastic reduction of the number of live microorganisms on instruments with certified thermal disinfection.

4. HOT AIR DRYING

A stream of hot air enables the complete removal of water from instruments of any material.

HYDRIM C61 WD G4



SAVE
\$1,300



NOW ONLY \$11,500 ~~\$12,800~~

SCI-C61WDD02

This bench top washer-disinfector is user-friendly, safe and effective, allowing for thorough cleaning as well as thermal disinfection at 93°C.

FAST AND EFFICIENT

The HYDRIM C61WD G4 has a 6 full size cassette capacity and can turn over 120 instruments in 60 minutes, start-to-dry.

SAFER INSTRUMENT HANDLING

Automated washer-disinfectors reduce physical handling of instruments, decreasing the risk of exposure and injuries. HYDRIM G4 ensures that every set of instruments receive the same validated cleaning process by dosing the correct amount of detergent with each cycle.

STATIM 2000 G4**COLTENE**
SciCan**SAVE**
\$901

The STATIM 2000 G4 autoclave with a 1.8L chamber, has the reprocessing capability of 4 pouches or 1kg load. It sterilises wrapped instrument loads in 35 mins.

NOW ONLY \$8,999 ~~\$9,900~~

SCI-G4152313

STATIM 5000 G4**COLTENE**
SciCan**SAVE**
\$1,861

The STATIM 5000 G4 autoclave with a 5.1L chamber, has the reprocessing capability of 10 pouches or 1.5kg load. It sterilises wrapped instrument loads in 40 mins.

NOW ONLY \$11,999 ~~\$13,860~~

SCI-G4222314

STATIM 6000B G4+**COLTENE**
SciCan**SAVE**
\$1,000

The STATIM B G4+ autoclave with a 6L chamber has a small footprint like the STATIM 5000. The reprocessing of 2 cassettes or up to 12 pouches offers enough capacity to meet the needs of any size practice. It perfectly sterilises and dries wrapped instrument loads in as little as 27 minutes.

NOW ONLY \$9,500 ~~\$10,500~~

SCI-G4700009

B FUTURA 22L



SAVE
\$2,790



THE FIRST STERILISER CAPABLE OF USING TAP WATER

The B Futura steriliser offers even faster performance with full traceability of each sterilisation cycle, also with a modern colour touch-screen perfectly embodying what we call Evolved Simplicity.

- New WiFi option supplied as standard
- Intuitive user interface for easy use
- Unique door locking mechanism
- Filter changes from the front of steriliser
- Programmable vacuum and helix test time
- Large HF Cassette load
- User management system allows users to be tracked to cycle data via personal PIN codes. Complies with AU & NZ standards and means no paper printout is required
- Optional cloud-based software to manage machine performance, cycle data and provide remote access for technical departments

NOW ONLY \$12,900 ~~\$15,690~~

22L

ST-M-BFUT22

3 year Warranty

Available in 17, 22 and 28L

**For more information, or to request a quote,
please contact your Henry Schein Equipment Specialist.**



EDWARD HARRISON

021 796 069

edward.harrison@henryschein.co.nz



DANIEL ELS

021 750 609

daniel.els@henryschein.co.nz



MICHAEL FINDLAY

022 016 5898

michael.findlay@henryschein.co.nz



PAUL MACDOUGALL

022 342 4900

paul.macdougall@henryschein.co.nz



Prices are effective from 1/10/2024 to 30/11/2024 and exclude GST, Freight, (E&OE) Installation and software unless otherwise specified.

Prices or offers in this publication are not applicable to government or corporate accounts.