



HIGHLIGHTS

April - June 2026

NEW ZEALAND



WHITENING

Call Courtney Gerlach: 021 400 323

Take-home whitening

3+1^x OR **10+5^x**
BUY 3 GET 1 FREE BUY 10 GET 5 FREE

POLA DAY

7.5% and 9.5% HYDROGEN PEROXIDE

- Fast acting: 30 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



POLA NIGHT

10%, 16% and 22% CARBAMIDE PEROXIDE

- From 45 minutes to overnight wear
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



POLA DAY CP

35% CARBAMIDE PEROXIDE

- Super fast: Just 15 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



POLA LIGHT

6% HYDROGEN PEROXIDE

- From just 20 minutes a day
- Portable, re-chargeable and easy to wear



POLA FOR ALIGNERS

10% HYDROGEN PEROXIDE and 35% CARBAMIDE PEROXIDE

- Straighten and brighten together
- The greatest patient comfort in a whitening kit



POLA DAY/NIGHT 4 X 1.3g MINI KIT
 POLA DAY/NIGHT 10 X 1.3g SYRINGE KIT
 POLA DAY/NIGHT 50 X 3g DISPENSER PACK
 POLA LIGHT 6% HP
 POLA FOR ALIGNERS

P
O
L
A

In-office whitening

POLA RAPID

38% HYDROGEN PEROXIDE

- Fast and simple to use
- Non-stick gel for faster application and removal
- Enhanced blue gel for superior application visibility



POLA RAPID 1 & 3 PATIENT KITS



POLA RAPID 3 PATIENT KIT
 POLA RAPID 3 PATIENT KIT W/O
 POLA RAPID 1 PATIENT KIT
 POLA RAPID 1 PATIENT KIT W/O

2+1^x
BUY 2 GET 1 FREE

BPA & HEMA FREE



STELA

THE FUTURE OF COMPOSITES

50 years of dental material experience has engineered a two step restorative that is stronger, gap free, easier to place and perfectly cures every time.

- ✓ Unlimited depth of cure (self cure)
- ✓ Gap-free interface
- ✓ Only 2 steps: 15 sec prep
- ✓ High strength
- ✓ Choose: capsule or automix



STELA AUTOMIX 1 SYR REFILL
 STELA AUTOMIX INTRO KIT
 STELA CAPSULE REFILL 50 CAP
 STELA CAPSULES INTRO KIT
 STELA PRIMER REFILL 5ML BT

BUY A STELA REFILL KIT GET A STELA PRIMER

FREE*

FREE



AUTOMIX REFILL - 8640001



CAPSULE REFILL - 8640003



PRIMER REFILL - 8640006

LUNA 2 | LUNA FLOW BUY 2 - GET 1 FREE



BPA
FREE

HEMA
FREE

SUPERIOR STRENGTH & AESTHETICS

LUNA 2

UNIVERSAL COMPOSITE

- Improved non-stick, non-slump handling and reliable mechanical properties
- Logical Shade Matching Technology for easier and accurate shade selection
- BPA free
- Very high radiopacity of 250% Al
- Flexural strength similar to dentine for Class III, IV and V restorations and proximal boxes
- Opalescent and fluorescent like the tooth for lifelike results



Photos courtesy:
Dr Neus Sancho Lecina - SPAIN

LUNA 2 4g SYRINGE
LUNA 2 20 x 0.25g COMPLETS
LUNA 2 INTRO KIT



Scan QR code for more information



HEMA
FREE



PERFECT FLOW PERFECT STAY

LUNA FLOW / LUNA FLOW LV

UNIVERSAL FLOWABLE COMPOSITE

- Exceptionally high compressive strength of 421 MPa and 128 MPa of flexural strength
- Flows on demand, holds its shape and does not string for maximum control over application
- Outstanding radiopacity of 265% Al for faster and more precise diagnosis
- Logical Shade Matching Technology for easier and accurate shade selection
- Two viscosity options for all clinical indications: Luna Flow offers medium/high viscosity and Luna Flow LV offers low viscosity




Photos courtesy:
Dr Carla Cohn - CANADA

LUNA FLOW/LV 2g SYRINGE
LUNA FLOW/LV BULK KIT 5 X 2g SYRINGES
LUNA FLOW/LV INTRO KIT 4 X 2g SYRINGES



Scan QR code for more information



CEMENTS

SET/SET PP

DUAL-CURED, SELF-ETCHING, SELF-ADHERING RESIN CEMENT

- Chemically bonds to tooth structure and substrates
- Available in capsules and a paste/paste automix syringe system (seT PP)
- Fast and easy - etch, bond and fill with one component



SET 50 CAPSULES
SET PP INTRO KIT
SET PP REFILL

BUY 2, GET 1 FREE*

ACCESSORIES

POINTS

DISPOSABLE BRUSH APPLICATORS

- Points disposable brush applicators are used to quickly transfer small amounts of material into areas that are difficult to reach. Points are resistant to common solvents in dental materials.
- Solvent resistant, non-absorbent, non-linting fibers
- Boxes of 400 in medium, fine or superfine



BOX OF 400 POINTS
DISPENSER BOX

**PLUS FREE
DISPENSER**



BUY 3, GET 1 FREE*

SEALANTS

CONSEAL F

LIGHT CURED FLUORIDE RELEASING RESIN PIT AND FISSURE SEALANT

- Extremely low viscosity - quickly and deeply penetrate small pits and tight fissures
- High Fluoride releasing
- BPA & HEMA free
- Accurate and controlled delivery with an ultrafine tip
- Opaque white shade
- 5mL bottles and 1g syringes

CONSEAL F SYRINGE KIT
CONSEAL F BULK SYRINGE KIT
CONSEAL F BOTTLE KIT
CONSEAL F 5ML BOTTLE



BUY 3, GET 1 FREE*

ADHESIVES

ZIPBOND

THE ALL-IN-ONE
ADHESIVE SOLUTION
FOR EVERY SITUATION

IMPROVED BOND STRENGTH AND
DEFENDS AGAINST CLINICAL VARIABILITY

- Consistent bond strengths to dentin and enamel
- Wet and dry dentin bonding
- Predictable bond strengths to indirect substrates
- Effective hybridization
- Exceptional control
- Highly esthetic



DISCOUNTED TO \$100.00

HENRY SCHEIN®
DENTAL

Rely on Us™

free phone 0800 808 855
orders@henryschein.co.nz
www.henryschein.co.nz

*All prices exclude GST. Offers cannot be combined. Valid on invoices dated April 1st - June 30th 2026. All claims for free goods must be made within 30 days from the end of promotion date. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Bonus goods must be of equal or lesser value. Offers are valid in New Zealand only.

Courtney Gerlach
Sales Product Specialist -
New Zealand
SDI Limited
m: +64 21 400 323
e: Courtney.Gerlach@sdi.com.au

SDI

YOUR
SMILE.
OUR
VISION.