


VistaCool™ - Direct-to-Drain Systems for Sterilizer Wastewater

VistaCool

Never empty your condenser bottle again.

VistaCool, the only direct-to-drain system designed from the ground up for your STATIM® 2000 and 5000 cassette autoclaves and other selected sterilizers.



Sophisticated wax motor allows cold water into the system ONLY when the sterilizer is in use.

VistaCool Double System
Code: **V7502**

VistaCool Single System
Code: **V7501**

FAST
Drains sterilizer exhaust safely and automatically.

EFFICIENT
Eliminates wasted staff time emptying condenser bottle.

COST EFFECTIVE
Dedicate more resources to patient care and less on maintenance.

FEATURES

UNIVERSAL SYSTEM

VistaCool will work with STATIM 2000 or 5000 and other select sterilizers.

SMART DESIGN

Triangular water tank is shaped to take up less space under your sink. Double unit services 2 sterilizers. No counter drilling or air gap (in most cases).

ECO DESIGN

Water tank is made of recycled plastic material.

EASY TO USE

No consumable parts nor daily maintenance is required.

ROBUST

System will continuously cool your sterilizers hot exhaust and drain it safely – no matter how many cycles you run.

COMPLIANT

All components and system set-up are compliant with most regional plumbing codes.

VALIDATED

VistaCool is the only direct-to-drain product designed to safely and effectively drain the steam exhaust from STATIM 2000 or 5000 and other select sterilizers.

COMPATIBLE

Can be used with all drainage pipe material types.

VistaCool is a trademark of Vista Research Group, Ashland, OH. Manufactured by Vista Research Group. PATENT PENDING. STATIM is a registered trademark and Your Infection Control Specialist is a trademark of SciCan Ltd. SD-463 CA EN FR R1.